



# ORDERING INFORMATION

## SINGLE FLOWMETERS – DIRECT MOUNTED



PLEASE SELECT OPTIONS FROM EACH SECTION:

1. MAIN FLOWMETER















|                                                                                    |                  |                                                                                      |               |
|------------------------------------------------------------------------------------|------------------|--------------------------------------------------------------------------------------|---------------|
|  |                  |  |               |
| MD 33530                                                                           | OXYGEN FLOWMETER | MD 33535                                                                             | AIR FLOWMETER |

# ORDERING INFORMATION

## SINGLE FLOWMETERS – DIRECT MOUNTED



### 2. INLET PROBE

|                 |          | FOR OXYGEN                                                                          |          | FOR AIR                                                                               |
|-----------------|----------|-------------------------------------------------------------------------------------|----------|---------------------------------------------------------------------------------------|
| British         | MD 33540 |    | MD 33547 |    |
| German          | MD 33541 |    | MD 33548 |    |
| Japanese        | MD 33542 |    | MD 33549 |    |
| French          | MD 33543 |  | MD 33550 |  |
| Ohmeda          | MD 32208 |  | MD 32210 |  |
| Chemetron       | MD 32203 |  | MD 32205 |  |
| Oxequip         | MD 32213 |  | MD 32215 |  |
| Puritan Bennett | MD 32218 |  | MD 32220 |  |

DODHYS MEDICAL LIMITED




P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.

# ORDERING INFORMATION

## SINGLE FLOWMETERS – DIRECT MOUNTED



### 3. CONNECTION AT OUTLET

|                                                                               |          |                                                                                      |
|-------------------------------------------------------------------------------|----------|--------------------------------------------------------------------------------------|
| Humidifier                                                                    | MD 33566 |   |
| Tubing Nipple – Plastic<br>(Add – W/B/G for white,<br>Black and Green colors) | MD 32277 |   |
| Tubing Nipple – Metal                                                         | MD 30024 |  |

### 4. ACCESSORIES

|            |          |
|------------|----------|
| Adult Mask | MD 30000 |
| Child Mask | MD 30001 |




# ORDERING INFORMATION

## PISTON TYPE REGULATORS




PLEASE SELECT OPTIONS FROM EACH SECTION:

### 1. MAIN REGULATOR




MD 33558

PISTON TYPE REGULATOR WITH  
FLOWMETER AND BULLNOSE SIDE  
CONNECTION




MD 33562

PISTON TYPE REGULATOR WITH  
FLOWMETER AND CGA-540  
CONNECTION



MD 33560

PISTON TYPE REGULATOR WITH  
FLOWMETER AND BULLNOSE  
INTERCHANGEABLE CONNECTION



MD 33561

PISTON TYPE REGULATOR WITH  
FLOWMETER AND PIN-INDEX  
CONNECTION

DODHYS MEDICAL LIMITED




P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.

[dodhys.com](http://dodhys.com)

# ORDERING INFORMATION

## PISTON TYPE REGULATORS

### 2. CONNECTION AT OUTLET

|                                                                              |          |                                                                                       |
|------------------------------------------------------------------------------|----------|---------------------------------------------------------------------------------------|
| Humidifier                                                                   | MD 33566 |    |
| Tubing Nipple – Plastic<br>(Add – W/B/G for white,<br>Black and Green colors | MD 32277 |   |
| Tubing Nipple – Metal                                                        | MD 30024 |  |

### 3. ACCESSORIES

|            |          |
|------------|----------|
| Adult Mask | MD 30000 |
| Child Mask | MD 30001 |




# ORDERING INFORMATION

## DIAPHRAGM TYPE REGULATORS



PLEASE SELECT OPTIONS FROM EACH SECTION:

### 1. MAIN REGULATOR



|          |                         |
|----------|-------------------------|
| MD 33570 | OXYGEN REGULATOR        |
| MD 33571 | NITROUS OXIDE REGULATOR |
| MD 33572 | AIR REGULATOR           |

### 2. INLET CONNECTION

|               | BULLNOSE | PIN INDEX | CGA                |
|---------------|----------|-----------|--------------------|
|               |          |           |                    |
| OXYGEN        | MD 33576 | MD 33579  | MD 33584 (CGA 540) |
| NITROUS OXIDE |          | MD 33580  | MD 33585 (CGA 326) |
| AIR           | MD 33576 | MD 33581  | MD 33586 (CGA 346) |

DODHYS MEDICAL LIMITED

P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.




[dodhys.com](http://dodhys.com)



# ORDERING INFORMATION

## DIAPHRAGM TYPE REGULATORS



### 3. CONNECTION AT OUTLET (Up to 3 types of outlet connection can be selected per regulator)

|                | BRITISH                                                                           | GERMAN                                                                             | FRENCH                                                                              |
|----------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| QUICK COUPLING |  |  |  |
| OXYGEN         | MD 31450                                                                          | MD 31495                                                                           | MD 31485                                                                            |
| NITROUS OXIDE  | MD 31451                                                                          | MD 31496                                                                           | MD 31486                                                                            |
| AIR 4 BAR      | MD 31452                                                                          | MD 31497                                                                           | MD 31487                                                                            |
| AIR 7 BAR      | MD 31453                                                                          |                                                                                    |                                                                                     |

|           |                                                                                     |                                                                                       |
|-----------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| FLOWMETER |  |  |
|           | MD 33530                                                                            | MD 33535                                                                              |
|           | OXYGEN FLOWMETER                                                                    | AIR FLOWMETER                                                                         |

DODHYS MEDICAL LIMITED

P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.




  [dodhys.com](http://dodhys.com)

# ORDERING INFORMATION

## DIAPHRAGM TYPE REGULATORS



### ACCESSORIES FOR THE FLOWMETER

|                                                                              |          |                                                                                      |
|------------------------------------------------------------------------------|----------|--------------------------------------------------------------------------------------|
| Humidifier                                                                   | MD 33566 |   |
| Tubing Nipple – Plastic<br>(Add – W/B/G for white,<br>Black and Green colors | MD 32277 |   |
| Tubing Nipple – Metal                                                        | MD 30024 |  |

### 4. ACCESSORIES

|            |          |
|------------|----------|
| Adult Mask | MD 30000 |
| Child Mask | MD 30001 |





## FEATURES AND BENEFITS

- Available for Oxygen, Nitrous Oxide and Air
- Adjustable Output Pressure Reducer
- Built-in Excess Pressure Relief Valve
- Stainless Steel body
- Tool-free installation
- Durable and robust

A highly versatile, adjustable output pressure regulator, available in a wide range of inlet and outlet configurations

Three simultaneous output connections possible

Two simultaneous output connections possible



Please refer the ordering information sheet for details about part numbers

**DODHYS MEDICAL LIMITED**

P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.

  [dodhys.com](https://www.dodhys.com)



**DISCLAIMER** (a) Please note that this brochure was prepared prior to all necessary Government consents being obtained. (b) The photographs may have been digitally enhanced or altered and are for indicative purpose only. (c) Changes may be made during the development and dimensions, fittings and specifications are subject to change without notice. Fittings and finishes are subject to availability and vendor discretion. (d) The information contained herein is believed to be correct but is not guaranteed. (e) Prospective purchasers should make and must rely on their own enquiries. (f) This brochure is a guide only and does not constitute an offer or contract.



Available in different connections

Please refer the ordering information sheet for details about part numbers

**DODHYS MEDICAL LIMITED**

P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.

  [dodhys.com](http://dodhys.com)





**DISCLAIMER** (a) Please note that this brochure was prepared prior to all necessary Government consents being obtained. (b) The photographs may have been digitally enhanced or altered and are for indicative purpose only. (c) Changes may be made during the development and dimensions, fittings and specifications are subject to change without notice. Fittings and finishes are subject to availability and vendor discretion. (d) The information contained herein is believed to be correct but is not guaranteed. (e) Prospective purchasers should make and must rely on their own enquiries. (f) This brochure is a guide only and does not constitute an offer or contract.




## FEATURES AND BENEFITS

- Factory-preset Fixed Output Pressure Reducer
- Built-in Excess Pressure Relief Valve
- Back Pressure Compensated Flowmeter
- Stainless Steel Valve
- Shatter proof tube cover and tube
- Stainless Steel body
- Large black ABS adjusting knob
- Aesthetically appealing colour scheme
- Tool-free installation
- Durable and robust



Factory pre-set, fixed output pressure regulators with flowmeters for attachment to cylinders only



Please refer the ordering information sheet for details about part numbers

**DODHYS MEDICAL LIMITED**

P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.

**f in dodhys.com**




**DISCLAIMER** (a) Please note that this brochure was prepared prior to all necessary Government consents being obtained. (b) The photographs may have been digitally enhanced or altered and are for indicative purpose only. (c) Changes may be made during the development and dimensions, fittings and specifications are subject to change without notice. Fittings and finishes are subject to availability and vendor discretion. (d) The information contained herein is believed to be correct but is not guaranteed. (e) Prospective purchasers should make and must rely on their own enquiries. (f) This brochure is a guide only and does not constitute an offer or contract.



## FEATURES AND BENEFITS

- Back Pressure Compensated
- Stainless Steel Valve
- Shatter proof tube cover and tube
- Stainless Steel body
- Large black ABS adjusting knob
- Aesthetically appealing colour scheme
- Durable and robust

Flowmeters with inlet connections/probes/adapters according to a wide range of international standards



Flowmeter to serve one patient at a time, attached to a rail on the wall or on a service unit

Please refer the ordering information sheet for details about part numbers

**DODHYS MEDICAL LIMITED**

P. O. Box: 7752, SAIF ZONE, SHARJAH, UNITED ARAB EMIRATES.

  [dodhys.com](https://www.dodhys.com)



**DISCLAIMER** (a) Please note that this brochure was prepared prior to all necessary Government consents being obtained. (b) The photographs may have been digitally enhanced or altered and are for indicative purpose only. (c) Changes may be made during the development and dimensions, fittings and specifications are subject to change without notice. Fittings and finishes are subject to availability and vendor discretion. (d) The information contained herein is believed to be correct but is not guaranteed. (e) Prospective purchasers should make and must rely on their own enquiries. (f) This brochure is a guide only and does not constitute an offer or contract.