

Infusion Stand Adjustable Height



IV-1105-2NP IV-1105-4NP IV-1104-2NC



IV-1104-4NC IV-1004-2P IV-1004-4P

Features:

- * Telescopic tubes made from non-magnetic Stainless Steel material.
- * Five legged Nylon base.
- * Mobile on 5 Twin wheel swivel castors Ø 50 mm.
- * Secure locking system for height adjustment.
- * Height adjustment from 120 cm to 207 cm.
- * 2 or 4 safety hooks.
- * Knock down construction.

Cat No:	Description
IV-1105-2NP	Infusion Stand, Nylon Base & Plastic 2 Hooks
IV-1105-4NP	Infusion Stand, Nylon Base & Plastic 4 Hooks
IV-1104-2NC	Infusion Stand, Nylon Base & S/S 2 Hooks
IV-1104-4NC	Infusion Stand, Nylon Base & S/S 4 Hooks
IV-1004-2P	Infusion Stand, Heavy Nylon Base & S/S 2 Hooks.
IV-1004-4P	Infusion Stand, Heavy Nylon Base & S/S 4 Hooks

Infusion Stand Adjustable Height

Features:

- * Telescopic tubes made from non-magnetic Stainless Steel material.
- * Four / Five legged base.
- * Mobile on swivel castors.
- * Secure locking system for height adjustment.
- * Height adjustment from 120 cm to 207 cm.
- * 2 or 4 safety hooks.
- * Knock down construction.



IV - 1114 - 2SS

IV - 1114 - 4SS

IV - 3222 - 4SS



IV - 1564 - HR

IV - 3255 - SSL

IV - 3265 - SS

Cat No:	Description
IV-1114-2SS	Infusion Stand, with S/S Base & S/S 2 Hooks.
IV-1114-4SS	Infusion Stand, with S/S Base & S/S 4 Hooks
IV-3222-SS	Infusion Stand, with Aluminum Base & S/S 4 Hooks (Suitable for Infusion Pumps)
IV-1564-HR	Infusion Stand, with S/S Base & S/S 4 Hooks (Suitable for Infusion Pumps)
IV-3255-SSL	Infusion Stand, with S/S Base & S/S 4 Hooks (Suitable for Infusion Pumps)
V-3265-SS	Infusion Stand, with S/S Base & S/S 4 Hooks (Suitable for Infusion Pumps)

Infusion Stand Adjustable Height



IV - 995 - 5SS

Features:

- * Telescopic tubes made from non-magnetic Stainless Steel material.
- * Four / Five leged base.
- * Mobile on swivel castors.
- * Secure locking system for height adjustment.
- * Height adjustment from 120 cm to 207 cm.
- * 2 or 4 safety hooks.
- * Knock down construction.

Cat No:

IV-995-5SS

Description

Infusion Stand Adjustable Height, With Extra Heavy Duty Base. (Suitable for Infusion Pumps)



SAFETY HOOKS



HEIGHT LOCKING SYSTEM




PLASTIC HANDLE (OPTIONAL)



SEAMLESS HEAVY BASE WITH 75MM ANTI STATIC WHEELS

Infusion Stand Adjustable Height



IV - 3216 - RB

Features:

- * Structure made from non-magnetic stainless steel material.
- * Stainless Steel Ultra Heavy duty Five legs base.
- * Mobile on swivel casters all with breaks.
- * Secure locking system for height adjustment.
- * Height adjustment from 120cm to 207cm.
- * Stainless steel 4 Safety hooks.
- * Complete Knock down construction.

Cat No:

Description

IV-3216-RB

Infusion Stand Adjustable Height, With Extra Heavy Duty Base. (Suitable for Infusion Pumps)



SAFETY HOOKS



HEIGHT
LOCKING
SYSTEM



PLASTIC
HANDLE
(OPTIONAL)



STAINLESS STEEL
HEAVY BASE WITH
75MM WHEELS



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

Examination Lights



HL-1220-PP



HL-1210-FH



HL-1220-WM

Features:

Non-magnetic Stainless Steel central tubes.

- * Reinforced nylon base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Flexible goose neck arm 700 mm.
- * Step down transformer 230V, 50W.
- * Halogen Bulb 12 V, 35 W, 10° angle of emission with UV block
- * Input voltage : 220/230V 50-60Hz.
- * Knock down construction.



EL-1230-PS

Cat No:

Description

HL-1220-PP

Halogen Examination Light,
Adjustable Height, Mobile.
(Halogen Bulb 12 V, 35 W)

HL-1210-FH

Halogen Examination Light, Fixed
Height, Mobile.
(Halogen Bulb 12 V, 35 W)

HL-1220-WM

Halogen Examination Light, Wall
Mounted.
(Halogen Bulb 12 V, 35 W)

EL-1230-PS

Examination Light Fixed Height,
Mobile.
(Normal Bulb 220 V, 60 W)

Features:

- * Non-magnetic Stainless Steel central tubes.
- * Reinforced nylon base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Flexible goose neck arm 500 mm.
- * Normal Bulb 220 V, 60 W.
- * Input voltage : 220/230V 50-60Hz.
- * Knock down construction.

Examination Lights

Features:

Non-magnetic Stainless Steel central tubes.

- * Reinforced nylon base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Flexible goose neck arm 700 mm.
- * Step down transformer 230V, 50W.
- * Halogen Bulb 12 V, 35 W, 10° angle of emission with UV block
- * Input voltage : 220/230V 50-60Hz.
- * Knock down construction.



LED - 0103 - FH

Features:

- * Non-magnetic Stainless Steel central tubes.
- * Reinforced nylon base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Flexible goose neck arm 500 mm.
- * Built in transformer, Suitable for power source 100V-277VAC.
- * LED Bulb with adjustable focus.
- * Knock down construction.



EL - 0030 - LEB

Features:

- * Non-magnetic Stainless Steel central tubes.
- * Stainless Steel Heavy base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Non flashing flourescent light with Magnifying Lens
- * Input voltage : 220/230V 50-60Hz.
- * Knock down construction



HL - 1310 - FH



HL - 1311 - FHD

Cat No:	Description
HL-1310-FH	Halogen Examination Light, Fixed Height, Mobile. Without Dimmer (Halogen Bulb 12 V, 35 W)
HL-1311-FHD	Halogen Examination Light, Fixed Height, Mobile. With Dimmer (Halogen Bulb 12 V, 35 W)
LED-0103-FH	LED Examination Light, Fixed Height, Mobile. With Adjustable Focus
EL-0030-LEB	Examination Light , Fixed Height, Mobile. With Magnifying Lens

Revolving Stool Adjustable Height



PS - 5638 - SS



PS - 5638 - LT



PS - 5639 - SL

Features:

- * Reinforced nylon base with 5 feet.
- * Mobile on 5 twin wheel swivel castors Ø 50 mm.
- * Gas lift for height adjustment from 490 - 680 mm.
- * Stainless Steel or upholstered Seat.
- * Knock down construction.

Cat No:	Description
PS-5638-LT	Revolving Stool Adjustable Height, Upholstered Seat & Nylon Base.
PS-5639-SL	Revolving Stool Adjustable Height, Upholstered Seat & S/S Base.
PS-5638-SS	Revolving Stool Adjustable Height, Stainless Steel Seat & S/S Base.

Revolving Stool Adjustable Height With Backrest



LN - 5770 - LPU



DC - 5670 - LPULT



PS - 5638 - BR



AS - 5578 - B

Features:

- * Reinforced nylon base with 5 feet.
- * Mobile on 5 twin wheel swivel casters Ø 50 mm.
- * Gas lift for height adjustment from 520 - 710 mm.
- * Seat with high back rest made from imitation leather upholstery with 50 mm thick high density foam padding,
- * Knock down construction.

Cat No:	Description
LN5770-LPU	Laboratory Stool Adjustable Height, With Low Back Rest
DC-5670-LPU	Operating Stool Adjustable Height, With High Back Rest.
PS-5638-BR	Patient Stool Adjustable Height, With Back Rest
AS-5578-B	Anesthesia Adjustable Height, With High Back Rest.

Blood Donor Chairs



BC - 9774 - PU

Features:

- * Reinforced nylon base with 5 feet.
- * Mobile on 5 twin wheel swivel casters Ø 50 mm.
- * Gas lift for height adjustment from 520 - 710 mm.
- * Seat and backrest upholstered with imitation leather with high density foam padding.
- * Padded armrest, adjustable in any position for patient conveniences.
- * Knock down construction.



BC - 9784 - PU

Cat No:

Description

BC-9774-PU

Blood Donor Chair With One Padded Armrest.

BC-9784-PU

Blood Donor Chair With Two Padded Armrest.

Blood Donor Chairs



MR - 5063 - GIV

TECHNICAL DETAILS

- * Dimensions with upright position: cm 74x90x127 h.
- * Dimensions with extended position: cm 74x170x73 h.
- * Seat width: cm 51.
- * Backrest inclination: 100° - 162°.
- * Country of Origin : Italy

Cat No:

Description

MR-5063-GIV

**Relax Armchair Special
without Armrest.**

Features:-

- * Structure made of tubular varnished steel.
- * Armrest Chair with fixed height and equipped also with separate movement of the backrest and the leg rest.
- * The separate inclination is made through gas spring moved by laterals levers.
- * The Seat and the backrest are upholstered with simil leather which is fireproof, washable and disinfectable.
- * backrest is equipped with push handle and inner lumbar support.
- * 125mm wheels all with independent brakes.



BC - 9784 - PU

Examination Couch (Quick Assembling)



ECF - 2970 - SS

Features:-

- * Structure made from non-magnetic Stainless Steel Tube.
- * Imitation leather upholstery with 50 mm thick foam padding.
- * Single hand operated adjustable back rest.
- * Provision for paper roll holder at both ends.
- * Holder for Paper Roll.
- * Quick assembly design.



PAPER ROLL
(OPTIONAL)



EXPORT PACKING



ECF - 2970 - SSW



ECF - 2970 - SSR

Cat No:

Description

ECF-2970-SS

Examination Couch With Roll Holder.
(L 186 cm x W 62 cm x H 75 cm.)

ECF-2970-SSW

Examination Couch With Roll Holder & Wheels.
(L 186 cm x W 62 cm x H 75 cm.)

ECF-2970-SSR

Examination Couch With Roll Holder & Side Rails.
(L 186 cm x W 62 cm x H 75 cm.)

ECF-2978-GS

Examination Couch With Gas Cylinder Operated
Head Rest
& Roll Holder (L 186 cm x W 62 cm x H 75 cm.)



ECF - 2978 - GS

Examination Couch Wooden



EC - 2970 - WD

Features:-

- * Structure made from Solid Wood or MDF (Medium Density Fiber Board).
- * Imitation leather upholstery with foam padding.
- * Fully Adjustable Head rest section.
- * Holder for Paper Roll

Cat No:

Description

EC-2970-WD

Examination Couch With Roll Holder.
(L 186 cm x W 62 cm x H 75 cm.)

EC-2970-WD1

Examination Couch With Gas Cylinder Operated Head Rest, Roll Holder
1 Drawer & Cabin with Sliding Doors (L 186 cm x W 62 cm x H 75 cm.)

EC-2970-WD2

Examination Couch With Gas Cylinder Operated Head Rest, Roll Holder
2 Drawers & 2 Cabins (L 186 cm x W 62 cm x H 75 cm.)



EC - 2970 - WD 2



EC - 2970 - WD 1



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

Examination Couch



ECW - 4070 - 8GS

Features:-

Heavy Duty Type (Stainless Steel).

* Strong Steel Tube Construction.

* Fully Adjustable Head rest section,
By Gas Spring.

* Two Section top with upholstery covered with
washable and flame retardant material.

* Complete welded construction.

Dimensions :- L 186 cm x W 62 cm x H 75 cm

Cat No:	Description
ECW-4070-8GS	Examination Couch with Gas Spring Operated Headrest, Suitable for Patients
ECF-2970-3SC	Examination Couch, 3 Section with Roll Holder, Suitable for Patients.
ECF-2973-3SC	Examination Couch, 3 Section with Roll Holder, Suitable for Patients.

Features:-

* Structure made from non-magnetic
Stainless Steel Tube.

* Imitation leather upholstery 3 Section with
50 mm thick foam padding.

* Single hand operated adjustable both sides.

* Provision for paper roll holder at both ends.

* Holder for Paper Roll.

* Quick assembly design.



ECF - 2973 - 3SC



ECF - 2970 - 3SC

Electric Examination Couch With Wheels

Features:-

- * Structure made from non-magnetic stainless steel
 - * Electrically Height Adjustment with Hand Operated Remote.
 - * Adjustable Head rest section, By Double Gas Spring.
 - * Two Section top upholstery covered with washable and flame retardant material.
 - * Base cover design enables free side access.
 - * Swivel wheel castors to move at short distance.
- Dimensions : L 186 cm x W 62 cm x H 50cm-90cm



EEC - 9250 - HW

Cat No:

Description

EEC-9250-
HW

Electric Examination Couch With Roll Holder.

EEC-9250-
3SC

Electric Examination Couch, 3 Section With Roll Holder

OPTIONS
AVAILABLE

Head Section With Nose Cut
Adjustable Arm Supports
Adjustable Leg Rests
Available With or Without Wheels
Available Assorted Colors



DOUBLE
GAS SPRING



REMOTE
CONTROL



WHEELS WITH
CENTER
LOCKING SYSTEM



MAXIMUM HEIGHT
90.0 CM



MINIMUM HEIGHT
50.0 CM



EEC - 9250 - 3SC

Massage Couch With Nose Cut.

Features:-

- * Structure made from non-magnetic Stainless Steel Tube.
- * Imitation leather upholstery with 50 mm thick foam padding.
- * Single hand operated adjustable back rest, with nose cut.
- * Provision for paper roll holder at both ends.
- * Holder for Paper Roll.



MC - 2000 - SS

Cat No:	Description
MC-1905-SS	Massage Couch With Nose Cut. (L 195cm x W 62cm x H 75cm.)
MC-2000-SS	Massage Couch With Nose Cut. (L 200cm x W 80cm x H 75cm.)

Features:-

- * Structure made from non-magnetic Stainless Steel Tube.
- * Imitation leather upholstery with 50 mm thick foam padding.
- * Single hand operated adjustable back rest, with nose cut.
- * Provision for paper roll holder at both ends.
- * Holder for Paper Roll.



MC - 1905 - SS

Gynecology & Delivery Examination Couch

Features:-

- Structure made from Non-magnetic stainless steel tube.
- * Imitation leather upholstery with 50mm thick foam padded.
- * Single hand operating adjustable backrest.
- * Extendable sliding foot section.
- * Lithotomy poles with PU leg support.
- * Stainless steel tray for drainage.
- * Paper roll holder.
- * Foldable Legs.
- * Quick Assembling Design.

Dimension :

Head Section : 77 cm x 61 cm

Seat : 57 cm x 61 cm

Height: 8+6 cm

:



D V - 8 8 0 0 - T B



G C F - 3 1 6 1 - S S

Cat No:	Description
GCF-3161-SS	Gynecology Examination Couch with Lithotomy Poles, Drainage Tray, Paper Roll Holder
DV-8800-TB	Delivery Table With Lithotomy Poles, Drainage Tray, Paper Roll Holder, Grip Handle and Shoulder Support



G R I P H A N D L E



S H O U L D E R S



P A P E R R O L L
H O L D E R



S T A I N L E S S S T E E L
D R A I N A G E
T R A Y



F O A M P A D D E D



E X P O R T P A C K I N G

Foot Stool Steps

Features:-

- * Frame made from non-magnetic Stainless Steel tube.
- * Step made from MDF, covered with anti slip rubber pad.
- * Rubber foot caps for stability.
- * Knock down construction.



FS - 4047 - DR




FS - 4036 - SR

FSD - 4147 - SR



FSR - 4137 - SR



Cat No:

Description

FS-4047-DR

Foot Stool, Double Step. (L 51 cm x W 40 cm x H 35.5 cm)

FS-4036-SR

Foot Stool, Single Step. (L 46 cm x W 28 cm x H 22 cm)

FS-6840-SR

Foot Stool, Single Step. (L 68 cm x W 40 cm x H 22 cm)

FSR-4137-SR

Foot Stool, Single Step With Hand Rail (L 46cm x W 28cm x H 91cm)

FSD-4147-DR

Foot Stool, Double Step With Hand Rail (L 51 cm x W 40 cm x H 91cm)



FS - 6840 - SR

Folding Ward Screen



PF - 1764 - N



PF - 1763 - N



PF - 1762 - N

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Folding panels, each 61cm wide.
- * PVC washable and flame retardant curtains for each section.
- * Mobile on twin wheel swivel castors Ø50mm, for easy movement.
- * Knockdown Construction.

Cat No:	Description
PF-1762-N	Ward Screen Two Folds, With PVC Curtain. (H175cm x Extended Width 122cm)
PF-1763-N	Ward Screen Three Folds, With PVC Curtain. (H175cm x Extended Width 184cm)
PF-1764-N	Ward Screen Four Folds, With PVC Curtain. (H175cm x Extended Width 245cm)

Instrument Cabinets



IC - 7 6 3 8 - 4 D



IC - 5 3 3 6 - 1 D



IC - 7 5 3 8 - 2 D

Features:-

- * Constructed from enameled steel sheet.
- * Lock-able glass/metal doors.
- * Adjustable glass shelves..
- * Elegant appearance.
- * Knock down construction.

Cat No:	Description
IC-5336-1D	Instrument Cabinet, 1 Door. (L 53 cm x W 36 cm x H 144 cm.)
IC-7538-2D	Instrument Cabinet, 2 Doors.(L 75 cm x W 38 cm x H 169 cm.)
IC-7638-4D	Instrument Cabinet, 4 Doors.(L 75 cm x W 38 cm x H 169 cm.)

Instrument / Storage Cabinets



IC - 1709045 - SS



IC - 1909045 - SS

Features:-

- * Constructed from Non-magnetic stainless steel Sheet.
- * Lock-able Stainless Steel Doors.
- * Adjustable/Fixed shelves.
- * Elegant appearance.
- * Complete welded construction



IC - 1209045 - SW



FC - 1208840 - SW



Cat No:

Description

IC-1709045-SS

Instrument/Storage Cabinet With 5 Adjustable Shelves & S/S Framed Glass Doors (W 90 cm x D 45 cm x H 170 cm)

IC-1909045-SS

Instrument/Storage Cabinet With 4 Fixed Shelves & S/S Framed Glass Doors (W 90 cm x D 45 cm x H 190 cm)

IC-1209045-SW

Instrument/Storage Cabinet With 2 Fixed Shelves & S/S Doors & Wheels (W 90 cm x D 45 cm x H 120 cm)

FC-1208840-SW

Filing Cabinet With S/S Doors, Drawer & Push Handle (W 88 cm x D 40 cm x H 120 cm)

Instrument / Storage Cabinet


Features:-

Structure made from non-magnetic Stainless Steel Materials.

- * Lock-able Stainless Steel Doors with Handle.
- * Sound-Absorbing Stainless Steel Drawer.
- * Stainless Steel Handle for Easy to Move.
- * Stainless Steel Pocket on Side of Cabinet. (Optional)
- * Mobile on Rubber wheel swivel castors Ø100mm / 75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Complete Welded Construction.



IC - 8090 - DHP



IC - 6040 - 2D



IC - 5545 - 2DW



IC - 5546 - DDW



DC - 5546 - DDW

Cat No:	Description
IC-8090-DHP	Instrument/Storage Cabinet with Shelves, Doors Push Handle. (L-80cm x W-50cm x H-100cm)
IC-6040-2D	Instrument/Storage Cabinet with Shelves, Doors Push Handle. (L-60cm x W-45cm x H-100cm)
IC-5545-2DW	Instrument/Storage Cabinet with Shelf, Doors Push Handle. (L-55cm x W-46cm x H-85cm)
IC-5546-DDW	Instrument/Storage Cabinet with Drawer, Doors Push Handle. (L-55cm x W-46cm x H-85cm)
DC-5546-DDW	Dressing/Instrument Cabinet with Drawer, Doors Push Handle. (L-55cm x W-46cm x H-85cm)

Instrument / Storage & Filing Trolley

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Lock-able Stainless Steel Doors on Both Sides.
- * Each compartment has 12 recesses (total capacity of the trolley: 24 recesses)
- * Four recesses are enclosed by a sliding door, complete with key-lock.
- * Zinc Plated Wheels 125mm / 200mm, 2 with independent Brakes.
- * Elegant appearance.
- * Complete Welded Construction.



IC - 1008545 - SS

Cat No:

Description

IC-1008545-SS

Instrument Cabinet With 2 Shelves S/S Doors, & Push Handle
(L 85 cm x W 45 cm x H 100 cm)

PRC-9650-24F

Patient Record Filing Trolley, with Doors and Side Rails.
L 96 cm x W 50 cm x H 96 cm.



PRC - 9650 - 24 F

Kick Bucket / X-Ray View Boxes



KB-1336-PU

KB-1830-PU

Features:-

- * 5 Leged non-magnetic Stainless Steel tube base.
- * 5 Vertical supports with rubber sleeves to hold the Bucket / Bowl.
- * Mobile on 5 twin wheel swivel castors Ø 50 mm.
- * One piece Stainless Steel Bucket / Bowl.

Cat No:	Description
M604043	X-Ray Film Viewer Single (One Panel). Dimensions : 40 cm x 43 cm
M608043	X-Ray Film Viewer Double (Two Panel). Dimensions : 80 cm x 43 cm
M612043	X-Ray Film Viewer Triple (Three Panel). Dimensions: 120 cm x 43 cm
M616043	X-Ray Film Viewer Tetra (Four Panel). Dimensions : 160 cm x 43 cm



M 6 1 6 0 4 3

Features:-

- * Steel Housing with scratch and chip resistant thermosetting epoxy-polyster powder finishing.
- * Unbreakable white monochromatic diffuser with supporting device.
- * Viewing Area of each panel 40cm x 42cm.
- * Secure X-Ray film retaining mechanism.
- * Instant start, Fluorescent cold light.
- * Two-way on/off switch.
- * Input voltage : 220/230V 50-60Hz.
- * Wall mounted model.

Cat No:	Description
KB-1336-PU	Kick Bucket (Base-Ø40cm x Height-36cm)
KB-1830-PU	Kick Bowl (Base-Ø40cm x Height-18cm)



M 6 0 4 0 4 3

M 6 0 8 0 4 3

M 6 1 2 0 4 3

LED X-Ray View Boxes

Features:-

- * The Flat Panel X-Ray illuminator design is truly an innovation.
- * It uses LED Light Source with Isometric Dot-Matrix Layout which results in High Luminance and Uniformity.
- * Slim film illuminator only 45mm thick compared to 5"-6" thickness of normal X-Ray Illuminator.
- * High Brightness of over 4000cd/m2 - Adjustable Brightness by Digital Regulator.
- * Blue light is more suited to X-Ray, MRI, CT Scan and other Films.
- * The Illuminator uses the latest LED Light Source which emits an 9,600 Kelvin color temperature.
- * This produces a blue-type light more suitable for reading black and white x-ray films.
- * Greatly reduced eye-fatigue due to High Frequency LED's of over 50kHz.
- * Adapted with advanced clipping device which features AUTOMATIC ON/OFF which turns on the illuminator when the film is inserted.



JD - A11

Cat No:	Description
JD-A11	LED X-Ray Film Viewer Single (One Panel). Dimensions : 35 cm x 43 cm
JD-A111	LED X-Ray Film Viewer Double (Two Panel). Dimensions : 70 cm x 43 cm



JD - A111



CONTROL PANEL

Pedal Bins With Lid



P B - 6 3 3 0 - 3 0

P B - 4 6 3 0 - 2 0

P B - 3 9 2 7 - 1 2

P B - 3 2 2 6 - 0 5

P B - 2 9 2 7 - 0 3

Features:-

- * Constructed from non-magnetic Stainless Steel sheet.
- * Easy to empty using the removable plastic inner bucket with handle.
- * Inner waste receptacle made from durable plastic.
- * The handle on the waste pedal bin make it extremely easy to move and transport.
- * Foot operated plastic lid rests on plastic ring when closed to avoid noisy collision.

Cat No:	Description
PB-6330-30	Pedal Bin With Lid, 30 Liters. Bucket Dimension:- Ø27 cm x H - 60 cm
PB-4630-20	Pedal Bin With Lid, 20 Liters. Bucket Dimension:- Ø27 cm x H - 39 cm
PB-3927-12	Pedal Bin With Lid, 12 Liters. Bucket Dimension:- Ø23 cm x H - 35 cm
PB-3226-05	Pedal Bin With Lid, 05 Liters. Bucket Dimension:- Ø18 cm x H - 22 cm
PB-2927-03	Pedal Bin With Lid, 03 Liters.. Bucket Dimension:- Ø15 cm x H - 20 cm

Stainless Steel Pedal Bins

Features:-

- * Constructed from non-magnetic Stainless Steel materials.
- * Rubber ring to hold the Garbage bag.
- * Foot operated lid with noiseless closing damper.
- * Stainless steel handle for operation.
- * Heavy Duty construction

Cat No:

Description

PBH-8939-70

Pedal Bin With Lid, Heavy Duty 70 Liters. (W 39 cm x D 35 cm x H 89 cm)

PBH-5939-50


Pedal Bin With Lid, Heavy Duty 50 Liters. (W 39 cm x D 35 cm x H 59 cm)

PBH-5739-40

Pedal Bin With Lid, Heavy Duty 40 Liters. (W 39 cm x D 35 cm x H 57 cm)

PBH-2724-20

Pedal Bin With Lid, Heavy Duty 20 Liters. (W 27 cm x D 24 cm x H 32 cm)



P B H - 2 7 2 4 - 2 0



P B H - 5 7 3 9 - 4 0



P B H - 5 9 3 9 - 5 0



P B H - 8 9 3 9 - 7 0

Dressing /Instrument Trolleys

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Flat trays made from non-magnetic Stainless Steel sheet.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.



Cat No:	Description
DT-5542-2PU	Dressing /Instrument Trolley 2 Shelves. (L 55cm x W 42cm x H 80cm)
DT-7050-2PU	Dressing /Instrument Trolley 2 Shelves. (L 70cm x W 50cm x H 80cm)
PBH-5739-40	Dressing /Instrument Trolley 2 Shelves. (L 90cm x W 55cm x H 80cm)
DT-5542-3PU	Dressing /Instrument Trolley 3 Shelves. (L 55cm x W 42cm x H 80cm)
DT-7050-3PU	Dressing /Instrument Trolley 3 Shelves. (L 70cm x W 50cm x H 80cm)
DT-9055-3PU	Dressing /Instrument Trolley 3 Shelves. (L 90cm x W 55cm x H 80cm)

Dressing /Instrument Trolleys With Side Rails

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Flat trays made from non-magnetic Stainless Steel sheet.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.



Cat No:	Description
DTR-5542-2PU	Dressing /Instrument Trolley 2 Shelves With Side Rail (L 55cm x W 42cm x H 80cm)
DTR-7050-2PU	Dressing /Instrument Trolley 2 Shelves With Side Rail (L 70cm x W 50cm x H 80cm)
DTR-9055-2PU	Dressing /Instrument Trolley 2 Shelves With Side Rail (L 90cm x W 55cm x H 80cm)
DTR-5542-3PU	Dressing /Instrument Trolley 3 Shelves With Side Rail (L 55cm x W 42cm x H 80cm)
DTR-7050-3PU	Dressing /Instrument Trolley 3 Shelves With Side Rail (L 70cm x W 50cm x H 80cm)
DTR-9055-3PU	Dressing /Instrument Trolley 3 Shelves With Side Rail (L 90cm x W 55cm x H 80cm)

Instrument Trolley, 2 Flat Shelves



ITF - 4 6 4 5 - 2 P U

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Flat trays made from non-magnetic Stainless Steel sheet.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.

Cat No:	Description
ITF-4645-2PU	Instrument Trolley, 2 Flat Shelves. (L 46 cm x W 46 cm x H 85 cm)
ITF-6045-2PU	Instrument Trolley, 2 Flat Shelves. (L 60 cm x W 46 cm x H 85 cm)
ITF-7645-2PU	Instrument Trolley, 2 Flat Shelves. (L 76 cm x W 46 cm x H 85 cm)



ITF - 6 0 4 5 - 2 P U



ITF - 7 6 4 5 - 2 P U

Instrument Trolley, 2 Flat Shelves With Side Rails

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Flat trays made from non-magnetic Stainless Steel sheet.
- * Stainless Steel Side Rails.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.



ITR - 4645 - 2PU

Cat No:	Description
ITR-4645-2PU	Instrument Trolley, 2 Flat Shelves. (L 46 cm x W 46 cm x H 85 cm)
ITR-6045-2PU	Instrument Trolley, 2 Flat Shelves With Side Rails (L 60 cm x W 46 cm x H 85 cm)
ITR-7645-2PU	Instrument Trolley, 2 Flat Shelves With Side Rails (L 76 cm x W 46 cm x H 85 cm)



ITR - 6045 - 2PU



ITR - 7645 - 2PU



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Instrument Trolley, 2 Shelves With Raised Edges

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Trays with raised edges, made from non-magnetic Stainless Steel sheet.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.



ITE - 4 6 4 5 - 2 P U

Cat No:	Description
ITE-4645-2PU	Instrument Trolley, 2 Flat Shelves With Raised Edges (L 46cm x W 46cm x H 85cm)
ITE-6045-2PU	Instrument Trolley, 2 Flat Shelves With Raised Edges (L 60cm x W 46cm x H 85cm)
ITE-7645-2PU	Instrument Trolley, 2 Flat Shelves With Raised Edges (L 76cm x W 46cm x H 85cm)



ITE - 6 0 4 5 - 2 P U



ITE - 7 6 4 5 - 2 P U

Dressing /Instrument Trolleys With Drawer



Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Trays with raised edges, made from non-magnetic Stainless Steel sheet.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.

Cat No:	Description
DT-5542-2DW	Dressing /Instrument Trolley 2 Shelves With 1 Drawer (L 55cm x W 42cm x H 80cm)
DT-7050-2DW	Dressing /Instrument Trolley 2 Shelves With 1 Drawer (L 70cm x W 50cm x H 80cm)
DTR-5542-2DW	Dressing /Instrument Trolley 2 Shelves With 1 Drawer & Side Rail (L 55cm x W 42cm x H 80cm)
DTR-7050-2DW	Dressing /Instrument Trolley 2 Shelves With 1 Drawer & Side Rail (L 70cm x W 50cm x H 80cm)



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Dressing /Instrument Trolleys With Drawer



Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Flat trays made from non-magnetic Stainless Steel sheet.
- * Stainless Steel Drawer with SS Handle.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.

Cat No:	Description
DT-5542-2DDW	Dressing /Instrument Trolley 2 Shelves With 2 Drawers. (L 55cm x W 42cm x H 80cm)
DTR-5542-2DDW	Dressing /Instrument Trolley 2 Shelves With 2 Drawers. & Side Rail (L 55cm x W 42cm x H 80cm)

Dressing /Instrument Trolleys With Drawer



Features:-

Structure made from non-magnetic Stainless Steel tube.

* Flat trays made from non-magnetic Stainless Steel sheet.

* Stainless Steel Drawer with SS Handle.

* Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.

* Impact absorbing Buffers for wall protection.

* Knockdown Construction.

Cat No:	Description
DT-5542-2DDW	Dressing /Instrument Trolley 2 Shelves With 2 Drawers. (L 55cm x W 42cm x H 80cm)
DTR-5542-2DDW	Dressing /Instrument Trolley 2 Shelves With 2 Drawers. & Side Rail (L 55cm x W 42cm x H 80cm)

Anesthesia Trolley With Drawers



ANT - 0WD6 - SS

Features:-

- Structure made from non-magnetic Stainless Steel tube.
- * Flat trays made from non-magnetic Stainless Steel sheet.
- * Stainless Steel Drawer with SS Handle.
- * Mobile on twin wheel swivel castors Ø75mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.

Cat No:	Description
ANT-0WD6-SS	Anesthesia Trolley, with 6 Drawers. (L - 70cm x W - 50cm x H - 95cm)
ANT-0WD5-SS	Anesthesia Trolley, with 5 Drawers. (L - 70cm x W - 50cm x H - 95cm)
ANT-0WD3-SS	Anesthesia Trolley, with 3 Drawers. (L - 90cm x W - 50cm x H - 95cm)



ANT - 0WD3 - SS



ANT - 0WD5 - SS

Instrument Trolley with Drawers



ITR - 5546 - 4DB

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Four Medium Capacity Drawers.
- * Two Removable S/S Bowls.
- * Rubber wheels 100 mm, all with independent Brakes.
- * Impact absorbing Buffers for wall protection.
- * Complete welded construction.

Dimensions:- L - 55cm x W - 46cm x H - 80cm

Cat No:	Description
ITR-5546-4DW	Instrument Trolley, with Side Rail, 4 Drawers and Double Bowls
ITR-5546-3DW	Instrument Trolley, with 3 Drawers and Side Rails on Bottom Shelf.

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Three Lock-able Medium Capacity Drawers.
- * Rubber wheels 100 mm, all with independent Brakes.
- * Impact absorbing Buffers for wall protection.
- * Complete welded construction.

Dimensions:- L - 55cm x W - 46cm x H - 80cm



ITR - 5546 - 3DB



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Utility Trolley



UT - 9045 - 3WS



UT - 9045 - 2WS

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Four Medium Capacity Drawers.
- * Two Removable S/S Bowls.
- * Rubber wheels 100 mm, all with independent Brakes.
- * Impact absorbing Buffers for wall protection.
- * Complete welded construction.

Dimensions:- L - 55cm x W - 46cm x H - 80cm

Cat No:	Description
UT-9045-3WS	Utility Trolley, with 3 adjustable Shelves. (L - 90cm x W - 45cm x H - 80cm)
UT-9045-2WS	Utility Trolley, with 2 adjustable Shelves. (L - 90cm x W - 45cm x H - 80cm)
UT-7645-3SS	Utility Trolley, with 3 adjustable Shelves. (L - 76cm x W - 45cm x H - 80cm)
UT-7645-2SS	Utility Trolley, with 2 adjustable Shelves. (L - 76cm x W - 45cm x H - 80cm)

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Handle on both Sides S/S.
- * Removable S/S Shelves.
- * Mobile on Rubber wheel swivel castors Ø100mm, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Knockdown Construction.



UT - 7645 - 2SS



UT - 7645 - 3SS

Multi Purpose Trolley



MPT-7050-WR

Features:-

- * Structure made from non-magnetic stainless steel materials Grade 304 AISI.
 - * Both Basket made from non-magnetic stainless steel bars.
 - * Mobile on Rubber wheel swivel castors Ø100mm, all with independent brakes.
 - * Impact absorbing Buffers for wall protection.
 - * Complete Welded fabrications
- Dimensions : L - 70cm x W - 50cm x H - 90cm.

Cat No:	Description
MPT-7050-WR	Multi Purpose Trolley, (Heavy Duty)
MPT-7050-RB	Multi Purpose Trolley, with Two Removable Baskets.

Features:-

- * Structure made from non-magnetic stainless steel materials Grade 304 AISI.
 - * Both Removable Basket made from non-magnetic Stainless steel bars.
 - * Mobile on Zinc Plated wheel swivel castors Ø100mm, all with independent brakes.
 - * Impact absorbing Buffers for wall protection.
 - * Complete Welded fabrications.
- Dimensions : L - 70cm x W - 50cm x H - 90cm.



MPT-7050-RB



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Operating Table Hydraulic

Features:-

- * Structure made from non-magnetic Stainless steel Materials.
- * Stainless Steel Sunken Top.
- * Foot Pedal on Both Sides for Height Adjustment.
- * Easily Control Height of table with Hydraulic Pump.
- * Complete Welded Construction.



OT - 7 5 7 0 - F T H

Cat No:

Description

OT-7570-FTH

Operating Table Adjustable Sunken Top, Hydraulic.
Dimensions: (L - 170cm x W - 75cm x H - 78cm - 110cm)

OT-7575-STH

Operating Table Adjustable Flat Top, Hydraulic.
Dimensions: (L - 170cm x W - 75cm x H - 78cm - 110cm)

Features:-

- * Structure made from non-magnetic stainless steel materials Grade 304 AISI.
 - * Both Removable Basket made from non-magnetic Stainless steel bars.
 - * Mobile on Zinc Plated wheel swivel castors Ø100mm, all with independent brakes.
 - * Impact absorbing Buffers for wall protection.
 - * Complete Welded fabrications.
- Dimensions : L - 70cm x W - 50cm x H - 90cm.



OT - 7 5 7 5 - S T H



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Instrument / Operating Table

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Flat / Sunken Top made from non-magnetic stainless steel Materials..
- * Sound-Absorbing Working Top.
- * Complete welded construction.

Overall Dimensions:

Width - 70cm

Length - 200cm

Height -86cm

Backsplash - 30cm



WT - 20070 - 2PU



WT - 10065 - SNK

Cat No:

Description

WT-20070-
2PU

Working Table For Kitchen, 2 Shelves
With Backsplash.

WT-10065-
SNK

Working Table, with Sunken Top Shelf.
(L-100cm x W-65cm x H-85cm)

OTT-
09065-2SS


Instrument / Operating Table, with 2
Shelves. (L-90cm x W-65cm x H-85cm)

OTT-11065-
2SS

Instrument / Operating Table, with 2
Shelves. (L-110cm x W-65cm x H-85cm)

OTT-
20065-2SS

Instrument / Operating Table, with 2
Shelves. (L-200cm x W-65cm x H-
85cm)



OTT - 09065 - 2SS

OTT - 11065 - 2SS



OTT - 20065 - 2SS



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Instrument / Operating Table With Wheels

Features:-

- * Structure made from non-magnetic Stainless steel Tube.
- * Flat Tray made from non-magnetic Stainless steel Sheet.
- * Mobile on Twin Wheel Swivel Castors Ø100mm, all with Independent brakes.
- * Impact Absorbing Buffer for Wall Protection.
- * Complete Welded Construction.



ITR - 7060 - 2PU



ITE - 7060 - 2PU

Cat No:

Description

ITR-7060-2PU

Instrument Trolley 2 Flat Shelves With Side Rails on Bottom Shelf. (L 70cm x W 60cm x H 85cm)

ITE-7060-2PU

Instrument Trolley, Top Flat Shelf & Bottom shelf with Raised Edges. (L 55cm x W 42cm x H 80cm)

ITT-7085-2PU

Instrument / Operating Table, 2 Flat Shelves (L - 70cm x W - 50cm x H - 85cm)

ITT-10085-2PU

Instrument / Operating Table, 2 Flat Shelves (L - 100cm x W - 55cm x H - 85cm)

ITT-12095-2PU

Instrument / Operating Table, 2 Flat Shelves (L - 120cm x W - 65cm x H - 95cm)



TT - 7085 - 2PU



ITT - 10085 - 2PU



ITT - 12095 - 2PU



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

Scrub Sink

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Sink made from non-magnetic stainless steel sheet.
- * Sound-Absorbing Working Top.
- * Stainless Steel Drainage System.
- * Complete Welded Construction.



SCK - 8045 - SS

Cat No:

Description

SCK-8045-SS

Single Sink Cabinet, with Double Doors and Backsplash.
(L - 80cm x W - 45cm x H - 85cm)

DBS-26072-2SK

Double Sink for Kitchen, with Bottom Shelf & Backsplash.
(L - 260cm x W - 72cm x H - 82cm)

TBS-30570-3SK

Triple Sink for Kitchen, with Bottom Shelf & Backsplash
(L - 305cm x W - 70cm x H - 87cm)



DBS - 26072 - 2SK



TBS - 30570 - 3SK

Patient Chart Holder Trolley

Features:-

- * Structure made from non-magnetic stainless steel Materials.
- * Sink made from non-magnetic stainless steel sheet.
- * Two Doors with Double Panel welded & sound deadened.
- * Stainless Steel Drainage System.
- * Complete Welded Construction.



CHT - 5 6 4 7 - 2 6 H

Cat No:

Description

CHT-5647-26H

Patient Chart Holder Trolley,
Double with 26 Chart Holders.
L 56 cm x W 47 cm x H 90 cm.

CHT-2947-13H

Patient Chart Holder Trolley,
Single with 13 Chart Holders.
L 29 cm x W 47 cm x H 90 cm.



CHT - 2 9 4 7 - 1 3 H



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Oxygen Cylinder Trolley

Features:-

- * Structure made from non-magnetic stainless steel tube.
- * Capacity for 12 / 9 cylinders.
- * Anti-slip rubber pad.
- * Mobile on PU swivel wheels Ø100mm, all with independent brakes.
- * Impact absorbing buffers for wall protection



CT - 12 SEC - SS



CT - 9 SEC - SS

Cat No:

Description

CT-12SEC-SS

Oxygen Cylinder Trolley, with 12 Sections.

CT-9SEC-SS

Oxygen Cylinder Trolley, with 9 Sections.

Features:-

- * Structure made from non-magnetic Stainless steel tube.
- * Easy loading of cylinder from the front.
- * Secure cylinder retaining mechanism.
- * High quality Ø100mm, rubber wheels ensure smooth ride, and are inoperative when trolley is stationary.

Cat No:

Description

OCT-048-CUF

Oxygen Cylinder Trolley,
24 x 2 CUF, Capacity.

OCT-024-CUF

Oxygen Cylinder Trolley,
24 CUF Capacity.



OCT - 024 - CUF

OCT - 048 - CUF



الاتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Mobile Racks, Urinals & Bedpan Holder

Cat No:	Description
RS-0024-SS	Oxygen Cylinder Trolley, with 12 Sections.

Features:-

- * Structure made from non-magnetic stainless Steel 18 / 10 Materiels.
- * Accommodate 24 Pairs of Shoes.
- * Mobile on Rubber wheel swivel castors Ø100mm, Two with brakes and Two Without Brakes.
- * Impact absorbing Buffers for wall protection. (OPTIONAL)
- * Complete Welded Constructio



RS - 0024 - SS



WLC - 9065 - 7SS

Cat No:	Description
WLC-9065-7SS	Wire Linen Cart with Cover, (L 90cm x W 65cm x H 220cm)

Features:-

- * Structure made from Non-magnetic Stainless Steel Tube.
- * Seven Non-Magnetic Stainless Steel Wire Shelf's.
- * Removable Nylon Cover.
- * Impact absorbing Buffers.
- * Mobile on rubber swivel wheels, all with Independent brake.
- * Complete Welded Construction

Cat No:	Description
UH-11634-SS	Urinals and Bedpan Holder, Wall Mounted. (L 100cm x W 30cm x H 80cm)

Features:-

- * Construction from Non-Magnetic tainless steel Materiél 304 AISI.
- * Holds up to Six Bedpans and Twelve Urinals .
- * Wall Mounted.
- * Complete welded construction.



UH - 11634 - SS



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Tissue Dispenser & Gloves Box Holder

Cat No:


Description

EK-62071

**Tissue Dispenser
Stainless Steel Wall
Mounted.**

Features:-

- * Structure made from Non-magnetic Stainless Steel Sheet.
- * Dust-proof design and maintain the cleanliness of toilet paper to improve the image of Health
- * Popular design, apply to hotels, guesthouses, restaurants, home, public places, and so on.
- * Wall Mounte



E K - 6 2 0 7 1

Cat No:

Description

GBH-2539-3B

**Latex Examination Gloves Box Holder
Triple. L 25.5 cm x W 8.5 cm x H 39 cm**

GBH-3923-3B

**Latex Examination Gloves Box Holder
Triple. (L 39 cm x W 8.0 cm x H 23 cm.)**

GBH-3030-1B

**Latex Examination Gloves Box Holder
Single. (L 11 cm x W 7.5 cm x H 20.5 cm.)**

GBH-8625-1B

**Latex Examination Gloves Box Holder
Single. (L 12 cm x W 8.5 cm x H 24 cm.)**

GBH-2586-1B

**Latex Examination Gloves Box Holder
Single. (L 25.5 cm x W 8.5 cm x H 13 cm.)**



G B H - 2 5 3 9 - 3 B



G B H -
2 5 8 6 -
1 B



G B H -
8 6 2 5 - 1 B



G B H -
3 0 3 0 - 1 B



G B H - 2 7 2 6 - 2 B



Features:-

- * Structure made from non-magnetic Stainless Steel Material.
- * Easy to Clean and Disinfect.
- * Wall Mounting Model with Four Screws for Secure Fixing.
- * Suitable for One Two or Three Gloves Box.

Wash Basin Trolleys

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Stainless Steel basin 6 liter capacity.
- * Mobile on wheel swivel castors, all with independent brakes.
- * Impact absorbing Buffers for wall protection.
- * Complete welded / knockdown construction



SBS - 0801 - FX



DBS - 0802 - FX

Cat No:	Description
SBS-0801-FX	Wash Basin Trolley with Single Basin & Handle.
DBS-0802-FX	Wash Basin Trolley with Double Basin & Handle
SBS-4438-N	Wash Basin Trolley with Single Basin. (L 45cm x W 40cm x H 84cm)
DBS-8438-N	Wash Basin Trolley with Double Basin. (L 80cm x W 40cm x H 84cm)
ST-2294-SS	Double Basin Sponge Trolley. (L 90cm x W 45cm x H 85cm)



SBS - 4438 - N



DBS - 8438 - N



ST - 2294 - SS

Mayo Instrument Trolleys



MT - 5 5 4 0 - N

MT - 5 6 4 0 - L

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Removable Stainless Steel Tray 530 x 400 mm.
- * Height Adjustment from 960 to 1500 mm with Secure locking system.
- * Mobile on twin wheel swivel castors Ø50mm.

Dimensions : L 53 cm x W 40 cm x H 96 - 150 cm.

Cat No:	Description
MT-5540-N	Mayo Instrument Trolley, With 4 Wheels, 2 with independent Brakes.
MT-5640-L	Mayo Instrument Trolley, With 2 Wheels.



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Mayo Instrument Trolley Hydraulic

Features:-

- * Structure made from non-magnetic stainless steel tube.
 - * Removable Stainless Steel Tray 53 x 40cm.
 - * 360° Rotating top with locking system.
 - * Height Adjust with foot pump.
 - * Height Adjustment from 96 to 150cm with Secure Locking System Hydraulic.
 - * Mobile on twin wheels swivel castor 75mm.
- Dimensions : L 53 cm x W 40 cm x H 96 - 150 cm



MT - 9000 - RT

Cat No:	Description
MT-9000-RT	Mayo Instrument Trolley with 360° Rotating Top, Hydraulic Classic.
MT-9002-FT	Mayo Instrument Trolley with Fixed Top, Hydraulic Classic.

Features:-

- Structure made from non-magnetic stainless steel tube.
- * Removable Stainless Steel Tray 53 x 40cm.
 - * Height Adjust with foot.
 - * Height Adjustment from 96 to 150cm with Secure Locking System Hydraulic.
 - * Mobile on twin wheels swivel castor 75mm.
- Dimensions: L 53cm x W 40 cm x H 96 - 150cm



MT - 9002 - FT



MINIMUM &
MAXIMUM
HEIGHT



360°
ROTATING
TOP



HYDRAULIC
PEDAL



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

BPM & Leg Dressing Stand Adjustable Height

Features:-

- * Telescopic tube made from non-magnetic Stainless steel material.
- * Stainless Steel basket for Hold Cuffs & Accessories.
- * Five lagged Nylon base.
- * Plastic Handle for Easy to Move.
- * Mobile on twin wheel swivel castors 50mm.
- * Secure locking system for height adjustment.
- * Height adjustment from 75cm to 125cm.
- * Knock down construction.



OMR-STA

OMR-STB

Cat No:	Description
OMR-STA	Blood Pressure Monitor Stand, Adjustable Height with S/S Bucket
OMR-STB	Blood Pressure Monitor Stand, Adjustable Height with S/S Bucket

Cat No:	Description
LDS-90PB-SS	Leg Dressing Stand, Adjustable Height, Stainless Steel Base.
LDS-92PB-NP	Leg Dressing Stand, Adjustable Height, Nylon Base.

Features:-

- * Telescopic tube made from non-magnetic Stainless steel material.
- * PU pad for leg support.
- * Five lagged heavy stainless steel or Nylon base.
- * Mobile on twin wheel swivel castors 50mm.
- * Secure locking system for height adjustment.
- * Height adjustment from 50cm to 70cm.
- * Knock down construction.



LDS-90PB-SS

LDS-92PB-NP



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Soiled Linen & Linen Changing Trolleys

Features:-

- * Structure made from non-magnetic Stainless Steel tube.
- * Stainless Steel Ring Ø 50cm to hold linen bag.
- * Removable Nylon Linen.
- * Mobile on twin wheel swivel castors, all with independent brakes.
- * Impact absorbing Buffers.



LT - 9550 - PU

Cat No:	Description
LT-9550-PU	Soiled Linen Trolley Single, with Removable Nylon Bag.
LT-9580-PU	Soiled Linen Trolley Double, with Removable Nylon Bags.
LC-5590-2RW	Linen Changing Trolley With Two Shelves & One Linen Bag



LT - 9580 - PU



LC - 5590 - 2RW



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Soiled Linen Trolley With & Without Lid

Features:-

- * Structure made from non-magnetic stainless steel tube.
- * Stainless Steel Ring for holding Nylon Bag.
- * Removable nylon linen Bag washable.
- * Rubber wheels Ø100 mm all independent with brakes.
- * Impact absorbing Buffers for wall protection.
- * Complete welded construction.

Dimensions:- Ring Dia 50 cm x H 90 cm



LTC - 9550 - PU

Cat No:	Description
LTC-9550-PU	Soiled linen Trolley Single Curved, without Lid, with Removable Nylon Bag
LT-9550-WLD	Soiled linen Trolley Single Curved with Lid, and Removable Nylon Bag.
LT-9580-WLD	Soiled linen Trolley Double Curved with Lid, and Removable Nylon Bag.



LT - 9550 - WLD



LT - 9580 - WLD

Patient Trolleys

Features:-

- * Structure made from non-magnetic Stainless Steel Tube.
 - * High pressure laminated flat surface.
 - * Mobile on rubber swivel wheels, all with Independent brake.
 - * High density foam mattress with fluid resistant and flame retardant cover.
 - * Provision For I.V. Pole.
 - * Knock down construction.
- Dimensions : L 200cm x W 62cm x H 75cm



PT - 6 0 2 0 - R W



PT F - 6 5 2 5 - R W

Cat No:	Description
PT-6020-RW	Patient Trolley With Mattress
PTF-6525-RW	Patient Trolley Adjustable Backrest With Mattress




COLLAPSIBLE SIDE
RAILS, STAINLESS
STEEL (OPTIONAL)



GAS CYLINDER
HOLDER,
STAINLESS
STEEL
(OPTIONAL)



STORAGE TRAY,
STAINLESS STEEL
(OPTIONAL)




INFUSION STAND
WITH 2 HOOKS,
STAINLESS STEEL
(OPTIONAL)

Stretcher Transport Trolley With Tilt

Features:-

- * Structure made from non-magnetic Stainless Steel Tube.
- * Quick tilt to trendelenberg position by means of Gas Cylinder.
- * Two Section platform with high ressure laminated surface.
- * Backrest adjustment by means of Gas Cylinder.
- * Push-Pull handle at both ends.
- * Mobile on rubber swivel wheels, all with independent brakes.
- * Two section high density foam mattress with fluid resistant and flame retardant cover.
- * Collapsible Side Rails.
- * Stainless Steel Utility Tray.
- * Stainless Steel Gas Cylinder.
- * Stainless Steel IV Pole.
- * Bumpers Protection.

Dimensions : L 200 cm x W 62 cm x H 75 cm



ST - 6 6 2 5 - A D

Cat No:	Description
ST-6625-AD	Stretcher Transport Trolley With Trendelenberg



TERNDELENBERG



STAINLESS STEEL
CYLINDER HOLDER



INFUSION
STAND S/S
2 HOOKS

Patient Trolley with Mattress

Features:-

- * Structure made from non-magnetic Stainless steel tube.
- * Two section platform with high pressure laminated surface.
- * Backrest adjustment by means of Gas Cylinder.
- * Mobile on Ø20cm, rubber swivel wheels, all with Independent brake.
- * High density foam mattress with fluid resistant and flame retardant cover.
- * Provision for I.V.Pole.
- * Collapsible Side Rails, Stainless Steel.
- * Storage Tray, Stainless Steel.
- * Gas Cylinder Holder, Stainless Steel.
- * Infusion Stand With 2 Hooks, Stainless Steel.
- * Complete welded construction.

Dimensions : L 180 cm x W 60 cm x H 75 cm



PTF - 6525 - AD


Cat No:	Description
PTF-6525-AD	Patient Trolley Two Section Complete with Mattress, Side Rails, Storage Tray, Gas Cylinder & I.V Pole.




STORAGE TRAY &
CYLINDER HOLDER
WHEELS WITH CENTER
LOCKING SYSTEM
(OPTIONAL)



I.V POLE WITH 2
HOOKS



COLLAPSIBLE
SIDE RAIL
WITH LOCKING
SYSTEM



STORAGE TRAY &
CYLINDER HOLDER
WHEELS WITH BRAKE
(OPTIONAL)

Patient Transport and Transfer System

Features:-

- * Structure made from non-magnetic stainless steel tube.
- * High pressure laminated flat surface.
- * Wheel base for sterile & non-sterile area.
- * Transferable platform.
- * High density foam mattress with fluid resistance & flame retardent cover.
- * System to use colour code to distinguish between material used in sterile & non sterile area.
- * Safety catch to link two wheelbase together for Patient Transfer.
- * Platform to move on roller from one wheelbase to another.
- * Platform complete with collapsible side rails & IV Pole.
- * Trendelenburg position achievable by advice by gas spring.
- * Wheel base mounted on 4 wheel swivel castors with locking mechanism.
- * Storage Tray for additional storage.
- * Oxygen cylinder holder.



TT - 2002 - RW



I.V POLE



OXYGEN CYLINDER HOLDER & STORAGE TRAY



MOVING LOCKING SYSTEM



SIDE RAIL LOCK



WALL PROTECTION BUMPERS



WHEEL BRAKE PEDALS

Cat No:	Description
TT-2002-RW	Patient Transport and Transfer System

Patient Chart Holders & Bedside Lockers

Features:-

- * Board made form 4 mm thick Plexiglass, provides smooth work surface whilst being light and easy to clean.
- * Stainless Steel hook to hang on the Bed.
- * Nickel plated spring clip with positive action.
- * Useable for A4 size paper.



CH - 3 4 2 3 - P P



CH - 2 3 3 4 - P P

Cat No:	Description
CH-3423-PP	Patient Chart Holder With Hook to Hang On The Bed. (Vertical) Dimensions : W 23 cm x H 34 cm
CH-2334-PP	Patient Chart Holder With Hook to Hang On The Bed. (Horizontal) Dimensions : W 34 cm x H 23 cm



BL - 5 0 4 6 - P U

Features:-

- * Structure made from medium density fiber board with high Pressure laminate surface.
- * Spacious drawer with roller slides for smooth movement.
- * Lower cupboard with single door and one inner shelf.
- * Mobile on twin swivel wheel castors Ø 50mm for easy handling and maneuverability.



BL - 5 0 4 8 - N

Cat No:	Description
CH-3423-PP	Patient Chart Holder With Hook to Hang On The Bed. (Vertical) Dimensions : W 23 cm x H 34 cm
CH-2334-PP	Patient Chart Holder With Hook to Hang On The Bed. (Horizontal) Dimensions : W 34 cm x H 23 cm

Over Bed Table

Features:-

- * Steel structure with scratch and chip resistant thermosetting epoxy-polyster powder finishing or Stainless Steel.
- * Top made from medium density fiber board with high Pressure laminate surface
- * Mobile on twin swivel wheel castors Ø 40mm for easy handling and maneuverability.
- * Height adjustment with automatic locking at desired height between 79cm - 110cm
- * Knock down construction.




OT - 8040 - SS

Cat No:	Description
OT-8040-SS	Over Bed Table Adjustable Height With Flat Top (Stainless Steel)
OT-9040-EC	Over Bed Table Adjustable Height With Flat Top (Powder Coated)
OT-9141-EC	Over Bed Table Adjustable Height With Tilt-able Top (Powder Coated)



OT - 9040 - EC



OT - 9141 - EC

Infant Cot & Child Cot

Features:-

- * Steel structure with scratch and chip resistant thermosetting epoxy-polyster powder finishing or Stainless Steel.
- * Sliding side rails on both sides.
- * High density foam mattress with fluid resistant and flame retardant cover.
- * Mobile on swivel wheel castors Ø 75mm, two with independent brakes.
- * Knock down construction.



IC - 1536 - RW

Cat No:	Description
IC-1536-RW	Child Cot With Sliding Side Rails, Mattress & Tilt. L 135cm x W 75cm x H 144cm
IC-9155-RW	Infant Cot With Sliding Side Rails & Mattress. L 91cm x W 56cm x H 124cm
CC-1480-RW	Child Cot With Sliding Side Rails & Mattress. L 135cm x W 76cm x H 136cm



CC - 1480 - RW



IC - 9155 - RW

Bassinet Trolleys

Features:-

Features:-

- * Structure made form non-magnetic Stainless Steel tube.
- * Strong Poly Carbonate crib with ribs to improve mattress ventilation.
- * Both sides tilt by simple mechanism.
- * Mobile on Rubber Swivel Wheels Ø 100 mm, all with independent brakes.
- * Mattress with water proof & flame retardant cover.
- * Impact absorbing Buffers for wall protection

Dimensions:-

Overall Length : 750 mm

Overall Width : 440 mm

Overall Height : 940 mm



BT - 7 6 9 3 - P U



B T C - 7 6 9 3 - 3 S S



LOWER SHELF FOR
ADDITIONAL STORAGE,
STAINLESS STEEL
(OPTIONAL)



MESH BASKET,
STEEL EPOXY POWDER
COATED (OPTIONAL)

Crash / Emergency Medical Cart

Features:-

- * Complete Metal Body.
- * Six Drawers with Break-a way lock.
- * Stainless Steel IV Pole with 2 SS hook.
- * One defibrillator shelf.
- * Oxygen tank Holder.
- * Sharp container set & CPR cardiac board.
- * 100 mm wheels 2 with brakes and 2 w/o brakes



MX - 37 E M G



NC - 37 0 6

Cat No:

Description

MX-37-EMG

Crash / Emergency Medical Cart (GREY).
(L - 69 cm x W - 51 cm x H - 110.5 cm)

MX-37-EMG
(R)

Crash / Emergency Medical Cart (RED).
(L - 69 cm x W - 51 cm x H - 110.5 cm)

NC-3706

Narrow Cart, (WHITE).
(L - 51 cm x W - 54.2 cm x H - 97.7cm)



DEFIBRILLATOR SHELF



SLIDING TRAY



2 DRAWERS WITH
PARTITION



IV POLE &
HANDLE



OXYGEN
CYLINDER TANK



CPR CARDIAC
BOARD

Anesthesia Cart

Features:-

- * Complete Metal Body.
- * Six Drawers with Break-a way lock.
- * Stainless Steel IV Pole with 2 SS hook.
- * One defibrillator shelf.
- * 12 Tilt Bin Organizers Small.
- * 5 Tilt Bin Organizers Large.
- * Oxygen tank Holder.
- * Sharp container set & CPR cardiac board.
- * 100 mm wheels 2 with brakes and 2 w/o brakes.



MX - 37 - ANE

Cat No:

Description

MX-37-ANE

Anesthesia Cart, Dark Blue. (L - 69cm x W - 51cm x H - 110.5cm)



IV POLE



TILT BIN ORGANIZERS (S)



TILT BIN ORGANIZERS (L)



ACCESSORY BRIDGE



ACCESSORY BRIDGE



SHARPS CONTAINER SET



WASTE CONTAINER W/LID

Instrument / Double Medical Cart

Features:-

- * Complete Metal Body.
- * Five Drawers Side by Side with Break-a way lock.
- * Stainless Steel IV Pole with 2 SS hook.
- * One defibrillator shelf.
- * Oxygen tank Holder.
- * Sharp container set & CPR cardiac board.
- * 100 mm wheels 2 with brakes and 2 w/o brakes.



D X - 3 4

Cat No:	Description
DX-34	Double Medical Cart, White. DX-34 L 112.5 cm x W 54.5 cm x H 101.0 cm
IT3504D	Instrument Trolley, With Four Drawers White. IT3504D L 100 cm x W 50 cm x H 90 cm

Features:-

- * Complete Metal Body.
- * Five Drawers Side by Side with Break-a way lock.
- * Stainless Steel IV Pole with 2 SS hook.
- * One defibrillator shelf.
- * Oxygen tank Holder.
- * Sharp container set & CPR cardiac board.
- * 100 mm wheels 2 with brakes and 2 w/o brakes.



I T 3 5 0 4 D

Folding Mobile Screen

Features:-

- * Screen and aluminum poles are easy to operate and clean.
- * Screen Material PS + HPL (High Pressure Laminate)
- * 4 Screen color option: White, Cream, Pink and grey Blue.
- * The Wall Clamp is available and tool free.
- * The castor can be taken out for clean.
- * The Castor self-adjustment design can keep stability on uneven ground.



FS165-06-WH

Cat No:

Description

DX-34

Double Medical Cart, White.
DX-34 L 112.5 cm x W 54.5 cm x H 101.0 cm



Same open length for FS & FW series	No. of panels	Part no.	Length (When open)		Height (w/ 3" casters)		Folded size	
			cm	inch	cm	inch	Width cm	Depth cm
	4	FS145-04	120	47.2"	145	57"	30	9.2
		FS165-04			165	65"		
	6	FS145-06	180	70.9"	145	57"	30	13.8
		FS165-06			165	65"		
	8	FS145-08	240	94.5"	145	57"	30	18.4
		FS185-08			165	65"		
	2 panels per set	FSA145-02	60	23.6"	145	57"	30	4.6
		FSA145-02C (w / caster)						
	2 panels per set	FSA165-02	60	23.6"	165	65"	30	4.6
		FSA165-02C (w / caster)						

Transportation Cart


Features:-

- * Metal with Powder Coated.
- * Single Smooth Shutter door with Key Lock.
- * ABS Trays, Wire Baskets, Catheter Hooks and Shelves provide you the organizing flexibility.
- * 3" Iron Casters 75mm 2 with brake & 2 Free.



TC - 3700

Cat No:	Description
TC-3700	Transportation Cart, Dark Blue ABS top with side-out shelf. TC-3700 L 70 cm x W 52 cm x H 110 cm





الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

HOSPITAL FURNITURES

ISSUE 23

Wall Cabinet Shutter Door

Features:-

- * Metal with Powder Coated.
- * Single Smooth Shutter door with Key Lock.
- * ABS Trays, Wire Baskets, Catheter Hooks and Shelves provide you the organizing flexibility.
- * 3" Iron Casters 75mm 2 with brake & 2 Free.



WCS6C

Cat No:


Description

WCS6C


Wall Cabinet Shutter Door W / Key Lock White.
WCS6C L 70 cm x W 65.5 cm x H 204 cm

- Total height of tray/basket per wall cabinet: 60" for single / 60" x 2 for double

- Interchangeable ABS trays, wire baskets, catheter hooks and adjustable shelves to adapt to future needs.



DIVIDER SYSTEM PACKAGE



Castors
Standard 5" steel casters
(2 w/brake, 1 w/direction
control, 1 free)



WHL

Push handle
(2pcs per set), can
assemble on left or
right or both sides

OPTIONAL



الإتقان للمعدات الطبية
AL ITQAN MEDICAL EQUIPMENT

CONTACT US:



info@itqanmedical.com



+971 7 221 1606 / +971 56 667 1855



www.alitqanmedical.com