



WATERLOO HEALTHCARE

Better Carts for Better Outcomes



PRODUCT CATALOG

emergency carts | anesthesia carts | isolation carts | bedside | treatment/procedure carts

WHY YOU SHOULD CHOOSE WATERLOO

MEDICAL STORAGE IS ALL WE DO

At Waterloo Healthcare our singular purpose is to address the ever-changing requirements of our long-time customers in healthcare. As requirements for security and the need for rapid access to critical supplies change we are poised and ready to create a purpose built cart for you.

Our sole focus is on manufacturing the best possible aluminum and steel procedure carts for the medical industry. Without the distraction of other unrelated product lines we have been able to accomplish this goal for well over 50 years.

Our in-house staff of engineers pride themselves on their ability to create custom designs for our customers. If we don't have a cart or accessory that serves your unique needs, chances are we can design something for you.

Waterloo offers a bumper to bumper **FIVE YEAR WARRANTY** that stands out in the industry. More than likely your facility already has Waterloo carts that have been serving you for years. We know of Waterloo carts manufactured in the 1970's that are still in service today...a testament to the quality we build into each and every cart we manufacture!

Let's face it, carts aren't rocket science, but getting the right configuration that will facilitate your staff's needs can seem like it. That's why our customer support team, national distributors and field sales reps are available to meet with you to help you match the right Waterloo cart for the job.

FREE DRAWER CUSTOMIZATION

Select ANY combination of drawer sizes and colors at no extra charge. Other cart companies charge for this privilege.

FREE PLASTIC TOP

Other companies charge extra for plastic tops, but with Waterloo carts, plastic tops are always **FREE!**

CARTS ARRIVE FULLY ASSEMBLED

Lift the top off the box and your carts are ready to roll into service! Drawer fronts, drawer handles and casters are all factory installed. The only assembly required are the two push handles, since you may have a preference on which side you'd like them. We even include the tool to mount them.
(Please note: optional accessories may require field installation)

TURNKEY DELIVERY AVAILABLE

For a very small "per cart" fee we can fully assemble ALL required accessories and deliver your carts palletted, shrinkwrapped and on dedicated trucks so they can be immediately rolled into service! For large cart projects (*Crash carts, Anesthesia carts, Bedside carts, etc*) this represents a huge labor cost savings to your facility and allows for immediate cart deployment!

WATERLOO IS THE MANUFACTURER

Other cart companies may appear to manufacture their own carts when in reality they have them manufactured for them. At Waterloo we ARE the manufacturer which allows us a greater level of control and flexibility on your cart orders. Want something that you don't see in our catalog? Give us a call – we love custom work!

LOCAL SALES REPRESENTATIVES

Local reps are available to bring in a cart for your evaluation, just give us a call to schedule an appointment!




TABLE OF CONTENTS

ALUMINUM OR STEEL?4-5

STANDARD FEATURES & SIZE OPTIONS6-7

SECURITY LOCKING SYSTEMS8-9

BUILD YOUR PART NUMBER10

COLOR OPTIONS11

MEDICAL CARTS

EMERGENCY CARTS12-17

ANESTHESIA CARTS.....18-23

ISOLATION CARTS & ORGANIZERS24-29

BEDSIDE/NURSE SERVER CARTS30-35

SPECIALTY/TREATMENT/PHLEBOTOMY CARTS.....36-39

ACCESSORIES

DIVIDERS & TRAYS.....40-41

UNICART ACCESSORIES42-46

CROSS-REFERENCE SPECIFICATION CHART47

ALUMINUM OR STEEL?

WHICH CART IS RIGHT FOR YOU?

ADVANTAGES OF ALUMINUM

Our ELITE Aluminum carts fly from procedure to procedure almost as if they had wings! That's because they're made of the same high grade aluminum as commercial airliners.

This cart was purpose built for maximum speed and efficiency when your work demands it.



UTRLA-333369-RED

ALUMINUM OR STEEL?

WHICH CART IS RIGHT FOR YOU?

ADVANTAGES OF STEEL

Our steel cart has been a healthcare industry standard for close to five decades!

It offers many of the same great features as our aluminum carts at a more economical price point.



UTGKU-333369-RED

STANDARD FEATURES ON ALL WATERLOO CARTS



SELF-CLOSING BALL BEARING SLIDES
High-Quality "auto closing" drawer glides with 48 ball bearings per slide provide easy drawer access.



STABILIZING BASE
Our "one piece" nylon reinforced polycarbonate stabilizer base provides a wider more stable platform that resists scuffs/scrapes while protecting your cart.



INCLUDED REMOVABLE TOP
Waterloo carts come standard with a removable/replaceable plastic top.



UNICART CONSTRUCTION
Fully welded unibody frame provides unsurpassed durability.



DURABLE POWDER-COATED SURFACES
Durable powder painted surfaces look great while protecting your cart. Cart exteriors available in Red, Electric Blue, and Light Gray. Drawer fronts available in 16 different colors!



ERGONOMIC PUSH HANDLES
Our "one-piece" handles fit naturally in the palm of your hand and provide better control than traditional horizontal push handles.



NO EXTRA CHARGE FOR CUSTOM DRAWER CONFIGURATIONS
Want two 3" drawers instead of a single 6" drawer? No problem and more importantly...**NO EXTRA CHARGE!** Create your own custom drawer configuration with all Waterloo carts. Your choice of 3", 6", 9" and 12" drawers.



FULLY ASSEMBLED
Drawer fronts and drawer handles installed at our factory mean your cart arrives ready to roll out of the box! The only thing you need to attach are the two push handles...and we even include the wrench!

SIZE OPTIONS



STANDARD WIDTH

29" Width x 24.5" Depth
(includes stabilizer base)

Drawer Interior Dimensions
22" Wide x 16.5" Deep

JUNIOR WIDTH

23" Width x 24.5" Depth
(includes stabilizer base)

Drawer Interior Dimensions
16" Wide x 16.5" Deep

Don't see the cart dimensions you need shown above?
Contact Customer Service for more options!

SECURITY LOCKING SYSTEMS

SELECTING THE RIGHT ACCESS METHOD

Waterloo has you covered when it comes to security. Whether you need controlled cart access or simple tamper evident locks, we have a locking configuration to meet your needs! All key locks are “Keyed differently” unless customer requested to be “Keyed alike.” Key override locks for electronic or mechanical keyless locks are “Keyed alike.”

KEY LOCKING SYSTEM



Standard key-lock that includes two keys. On multiple cart orders, we will key all carts differently unless you request them to be keyed alike.

LEVER/BREAKAWAY LOCKING SYSTEM



Our centralized lever-lock provides easy access for your tamper evident breakaway seals. It also serves as a “Transport Lock”, keeping drawers in place during rapid transit.

MECHANICAL PUSH BUTTON



Go keyless with our reliable mechanical push-button lock. No batteries required and the combination is user adjustable.

GATE LOCKING SYSTEM



Individual gate locks allow you to independently secure drawers with either breakaway or key locks. You can even group drawers together by department with larger gate locks.

(Available in 3” increments up to a full 30” bar)



UTGQA-333369-FWG

UTRLA-33669-RED
Shown with optional gatelocking system

SECURITY LOCKING SYSTEMS

SELECTING THE RIGHT ACCESS METHOD

Waterloo also offers Keyless Locking options with and without WiFi Capability, allowing you to remotely change passcodes and access usage activity from a central server. We can even offer customized solutions that would allow you to lock certain drawers independently.

ELECTRONIC KEYPAD



- 250 Individual Users
- Tracks last 1,500 unlock events
- Manually reprogram access codes at the cart
- Quiet auto-relock after a preset amount of time

PROXIMITY READER WITH ELECTRONIC KEYPAD



- 250 Individual Users
- Tracks last 1,500 unlock events
- Manually reprogram access codes at the cart
- Quiet auto-relock after a preset amount of time
- Proximity/Employee Badge Access

ELECTRONIC KEYPAD WITH CONTROLLED SUBSTANCE DRAWER (WIFI CAPABLE)



- 3,000 Individual Users
- Tracks last 15,000 unlock events
- Reprogramming access codes via WiFi through Lockmanager
- Quiet auto-relock after a preset amount of time
- Access code opens all drawers except top “controlled substance” drawer

ELECTRONIC KEYPAD (WIFI CAPABLE)



- 3,000 Individual Users
- Tracks last 15,000 unlock events
- Reprogramming access codes via WiFi through Lockmanager
- Quiet auto-relock after a preset amount of time

PROXIMITY READER WITH ELECTRONIC KEYPAD (WIFI CAPABLE)



- 3,000 Individual Users
- Tracks last 15,000 unlock events
- Reprogramming access codes via WiFi through Lockmanager
- Quiet auto-relock after a preset amount of time
- Proximity/Employee Badge Access

PROXIMITY READER WITH ELECTRONIC KEYPAD & CONTROLLED SUBSTANCE DRAWER (WIFI CAPABLE)



- 3,000 Individual Users
- Tracks last 15,000 unlock events
- Reprogramming access codes via WiFi through Lockmanager
- Quiet auto-relock after a preset amount of time
- HID card opens all drawers except top “controlled substance” drawer
- Proximity/Employee Badge Access

EASY AS 1, 2, 3


BUILD YOUR PART NUMBER

SELECT CART OPTIONS

- 1. Choose your cart width and height:
Width: STANDARD, JUNIOR OR WIDE
Height: SHORT, MEDIUM, TALL OR X-TALL
- 2. Select any combination of 3", 6", 9" OR 12" drawers that fit within your selected carts configurable drawer space
- 3. Select drawer front colors

Each Alpha or Numeric Digit has a meaning

EXAMPLE: **UTGKA-33669-RED**



SECTION 1

1st Letter (Cart Width)	2nd Letter (Cart Height)	3rd Letter (Exterior Color)	4th Letter (Lock System)	5th Letter (Cart Material)
U = Standard J = Junior W = Wide	S = Short M = Medium T = Tall X = X-Tall	B = Elec. Blue G = Gray R = Red	K = Key L = Lever S = Simplex E = Electronic P = Proximity	A = Aluminum U = Steel

SECTION 2

1st Number (Height of 1st Drawer)	2nd Number (Height of 2nd Drawer)	3rd Number (Height of 3rd Drawer)	4th Number, etc. (Height of 4th Drawer)
--------------------------------------	--------------------------------------	--------------------------------------	--

The total for Section 2 must equal the total Configurable Drawer Space for your Cart Size:

- Short Cart = 21" Total
- Medium Cart = 24" Total
- Tall Cart = 27" Total
- X-Tall Cart = 30" Total

SECTION 3

3 letters signify color choice for drawer fronts (Refer to p.11 for a color swatch)

WHT = White	PNK = Pink	SPG = Spring Green
LTG = Light Gray	VIL = Violet	FWG = Fairway Green
NTL = Natural	DKB = Dark Blue	GRA = Charcoal
YEL = Yellow	ELB = Electric Blue	SLT = Salt Ash
ORG = Orange	LTB = Light Blue	CHY = Cherry
RED = Red	TLG = Teal	WNT = Walnut
WNE = Wine		

COLOR OPTIONS

DURABLE FINISHES

All Waterloo carts are coated with a highly durable powder paint that is designed to withstand the demands of the medical environment. With 16 colors to choose from and the ability to mix and match drawer colors on a single cart you can create a unique cart to fit all your specific needs. You can also upgrade to one of our 3 wood grain options.



White WHT	Light Gray LTG	Natural NTL	Yellow YEL	Orange ORG	Red RED	Wine WNE	Pink PNK	Violet VIL	Dark Blue DKB
						OPTIONAL WOOD GRAINS (UPGRADE)			
Electric Blue ELB	Light Blue LTB	Teal TLG	Spring Green SPG	Fairway Green FWG	Charcoal GRA				

ELITE ALUMINUM
EMERGENCY CARTS

Rapid and organized response to a cardiac emergency is essential. The Waterloo aluminum cart is 30% lighter than a standard steel cart, offering faster response times. The included steering caster allows the cart to take tight corners further improving maneuverability.

Work surface and drawer organization are vital during a Code. Our dual slide out shelves virtually double your cart top work surface, allowing you to move equipment and supplies within easy reach.



UTRLA-333369-RED



FEATURE	ELITE ALUMINUM	COMPARISON	TRADITIONAL STEEL
Unicart Construction	☒		☒ Unicart Construction
Self-closing Ball Bearing Slides	☒		☒ Self-closing Ball Bearing Slides
Easy-roll Casters	☒		☒ Easy-roll Casters
Dual Pullout Work Shelves	☒		☒ Single Pullout Work Shelf
Universal Accessory Mounting System	☒		☐ Universal Accessory Mounting System
Includes Drawer Divider Set	☒		☐ Includes Drawer Divider Set
Lightweight Construction	☒		☐ Lightweight Construction

TRADITIONAL STEEL
EMERGENCY CARTS

Our economical steel crash carts offer a robust work platform that can be configured to address all of your code cart needs. Included side shelf provides additional work surface while optional trays and dividers can further organize the contents of your cart.

If you are looking at adding the optional plastic trays/dividers, you may wish to consider our aluminum carts. The included dividers and dual pull out side shelves that come standard on our aluminum carts may offset the additional cost to upgrade.



UTRLU-333369-RED

EMERGENCY CARTS

Crash carts come in all shapes and sizes, and Waterloo manufactures a wide variety of cart sizes designed to fit your specific needs.

If you'd like to color code drawers to make them easy to identify quickly (example: White airway drawer), we are happy to accommodate that request at no extra charge! Just tell us what drawer front and what color and we can personalize your carts for you!

We can also outfit your crash carts with individual breakaway gatelocks over various drawers. These not only save re-stocking time, but add extra security to drawers that need it (See page 8 for lock options).



EMERGENCY CARTS



ER-2000
The Waterloo ER-2000 takes an innovative approach to crash carts and is designed to offer full access to all drawers in a code with the push of a button.

One press of the centrally located latch opens the dual front doors exposing all drawers for rapid access.
Drawer trays can be removed to be positioned where needed during a code.

RECOMMENDED ACCESSORIES

EMERGENCY CARTS (FOR MORE ACCESSORIES SEE PAGE 40)



IV POLE
Chrome IV pole and mounting bracket. Mounts to side of all Standard Carts, Junior and Medtrix. [IV-2](#)



WRAP-AROUND ALUMINUM RAILING
[R-3 \(Steel Standard Carts\) \(pictured\)](#)
[R-4 \(Aluminum Standard Carts\)](#)
[R-JR \(Steel Junior Carts\)](#)
[R-JRA \(Aluminum Junior Carts\)](#)



METAL OXYGEN TANK BRACKET
Can be mounted to all Standard and Junior Carts. [OTB-5](#)



SUCTION SHELF
For Standard or Junior Carts. Fits most brands of portable suction pumps. 18"(w) x 8"(d) x 3"(h) [SS-M](#)



3" STANDARD CART DIVIDER TRAY SYSTEM
Includes: White tray with adjustable clear sub-dividers. Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers. 21.5"(w) x 15.25"(d) x 2.3"(h) [DIVTRAY-STD](#)



5" STANDARD DIVIDER TRAY SYSTEM
White tray with adjustable clear sub-dividers. Fits in Steel and Aluminum full-size Standard only. 21.5"(w) x 15.25"(d) x 4.5"(h) [DIVTRAY-STD5](#)



3" GRAY PLASTIC TRAY WITH ADJUSTABLE DIVIDERS
Includes: 4 short dividers and 4 long dividers. Fits 2 per drawer in Steel and Aluminum Standard Carts. Fits 1 per drawer in a Steel or Aluminum Junior Cart. 9.25"(w) x 15"(d) x 2.5"(h) [WDC-2025](#)
Sold by the each



5" GRAY PLASTIC TRAY WITH ADJUSTABLE DIVIDERS
Includes: 4 short dividers and 4 long dividers. Perfect for 6", 9", 12" drawers. Fits 2 per drawer in Steel and Aluminum Standard Carts. Fits 1 per drawer in a Steel or Aluminum Junior Cart. Sold by the each. 9.25"(w) x 15"(d) x 4.5"(h) [WDC-2050](#)
Sold by the each



SHARPS CONTAINER MOUNTING BRACKET
Mounts to side of all Standard Carts, Junior and Medtrix. Designed for 5 quart container. **NON-LOCKING** [SB-2](#)
LOCKING (pictured) [SB-2L](#)



THIN CARDIAC BOARD BACKSIDE MOUNT
Includes cardiac board and brackets for backside of Steel or Aluminum Standard Carts. 21"(w) x 0.5(d) x 15.25"(h) [CBB-9](#)



SWIVELING MONITOR SHELF
Includes an adjustable strap to secure items in place. Maximum weight capacity 25 lbs. Can be locked in place. [ER-SHLF \(Steel Carts\)](#)
[ER-SHLFA2 \(Aluminum Carts\) \(pictured\)](#)
[ER-SHLF2000 \(ER2000\)](#)



INDIVIDUAL GATELOCKS WITH BREAKAWAY SEAL
Individual gate locks allow you to independently secure drawers with either breakaway or key locks. [GLB-3](#)
See page 8 for more details



CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Score the plastic strips and snap them to create your own divider system. *Includes:* (6) 22" long strips and 8 clips. Additional clips available in a set of 6. (To fit into Steel Junior carts cut strips down to 16.5" long) [DIV-CC2 \(Dividers & Clips\)](#)
[DIV-CLIP \(Clips\)](#)



LARGE CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Includes: 3 front to back panels, 2 side to side panels and 12 clips. Additional panels are available to create more compartments if needed. (To fit into a Steel Junior cart cut strips down to 16.5" long) [DIV-UC6 \(6" Drawers\)](#)
[DIV-UC9 \(9" & 12" Drawers\) \(pictured\)](#)



DIVIDER KIT FOR 3" ALUMINUM STANDARD DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider. [DIV-AU3](#)

Front to back divider [Y700286A1](#)
Side to side divider [Y700287A1](#)



DIVIDERS FOR 6", 9" & 12" ALUMINUM STANDARD DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider. [DIV-AU6](#)

Front to back divider [Y700286A2](#)
Side to side divider [Y700287A2](#)

DIVIDERS & TRAYS

EMERGENCY CARTS

ELITE ALUMINUM
ANESTHESIA CARTS

Organization and work surface are important features on an anesthesia cart and we offer both. With included metal drawer dividers in every drawer, our aluminum cart is organized right out of the box. Pair that with a variety of above cart tilt bins and shelving options and you have all the flexibility needed to organize equipment and supplies for your operating room.



HOW
“CONTROLLED”
ARE YOUR
MEDICATIONS?

- Dual slide out work surfaces virtually double your work surface!
- Plastic top with cutaway front designed to be wrist-friendly
- Large stylish and contemporary drawer handles offer a more updated look and enhanced functionality
- Adjustable metal drawer dividers are included in all aluminum cart drawers (One long and One short)
- Scratch-resistant, nylon reinforced polycarbonate bumper with a wider wheel base
- 5" Premium "Easy-roll" sealed bearing casters with precision bearing swivels. Thread guards minimize build-up around the axles, keeping carts rolling smoothly.




UTGQA-333369-DKB

FEATURE	ELITE ALUMINUM	COMPARISON	TRADITIONAL STEEL
Unicart Construction	●	●	Unicart Construction
Self-closing Ball Bearing Slides	●	●	Self-closing Ball Bearing Slides
Easy-roll Casters	●	●	Easy-roll Casters
Dual Pullout Work Shelves	●	●	Single Pullout Work Shelf
Universal Accessory Mounting System	●	○	Universal Accessory Mounting System
Includes Drawer Divider Set	●	○	Includes Drawer Divider Set
Lightweight Construction	●	○	Lightweight Construction

TRADITIONAL STEEL
ANESTHESIA CARTS

Economical steel anesthesia carts offer a robust work platform that can be configured with accessories to facilitate your work in the OR. If you are looking at adding the optional plastic trays/dividers, you may wish to consider our aluminum carts. The included dividers and dual pull out side shelves that come standard on our aluminum carts may offset the additional cost to upgrade.



- Single pull-out shelf provides additional work surface
- Smooth, easy-to-clean, removable work surface
- Pre-drilled holes to mount accessories
- Scratch-resistant, nylon reinforced polycarbonate bumper with a wider wheel base
- 5" Premium "Easy-roll" sealed bearing casters with precision bearing swivels. Thread guards minimize build-up around the axles, keeping carts rolling smoothly.

UTGSU-333369-DKB

ANESTHESIA CARTS

Anesthesia carts come in all shapes and sizes and various security requirements. Waterloo offers a wide variety of locking options — From a simple key lock to a WiFi capable proximity card reader, we can provide the right level of security for your operating rooms.

Further enhance cart security with an independently locking “controlled substance” drawer to keep drugs secure. One swipe of an approved employee badge allows access to all drawers, with the exception of the top drug drawer. A secondary key-code allows access to the controlled substance drug drawer keeping important drugs secure.



Aluminum carts only



UTGQA-333369-FWG

UMGSA-333369-VIL

USGKA-3369-LTG

ANESTHESIA CARTS



UTGKU-33669-DKB

UXGSU-333669-ELB

WXGKU-333669-FWG

Overall Dims: 36”(w) x 24.5”(d) x 45”(h)
Drawer Dims: 29”(w) x 17”(d)

Available ONLY in X-Tall, Steel with
FWG or LTG drawer front colors

RECOMMENDED ACCESSORIES

ANESTHESIA CARTS

(FOR MORE ACCESSORIES SEE PAGE 40)



TILT BINS
Clear plastic tilt bins with white exterior. Mounts to base rail system (BRS-2, BRS-3 OR BRS-JR). Available in 3, 4, 5, or 6 bin units.
TB-3 (pictured)
TB-4
TB-5
TB-6 (pictured)



HORIZONTAL RAIL SYSTEM
Includes two vertical and three horizontal rails.
BRS-2 (Standard Carts) (pictured)
BRS-JR (Junior Carts)
BRS-2W (Wide Carts)

(Shown on aluminum cart)



CLEAR PLASTIC ORGANIZER
SMALL contains 3 compartments. Mounts to accessory railing. 9"(w) x 5"(d) x 3"(h)
ORG-15



CATHETER HOLDER WITH CLEAR ACRYLIC FRONT
Catheter holder for Standard or Junior carts. Clear acrylic front slides in and out and provides viewing of contents.
4"(w) x 4"(d) x 19"(h)
CH-2



3" ANESTHESIA DIVIDER TRAY SYSTEM
Includes: White plastic tray with dividers. Fits in Steel Standard Carts ONLY. Perfect for 3" drawers.
21.5"(w) x 15.25"(d) x 2.3"(h)
DIVTRAY-ANES



3" AMPULE HOLDER DIVIDER TRAY SYSTEM
Clear dividers are adjustable. Fits into Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
21.5"(w) x 15.25"(d) x 2.3"(h)
DIVTRAY-AMP



3" WHITE PLASTIC TRAY WITH 13 BINS
Can be easily reconfigured. Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
22"(w) x 15"(d) x 2.5"(h)
TT-1A



3" WHITE PLASTIC TRAY WITH 10 BINS
Can be easily reconfigured. Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
22"(w) x 15"(d) x 2.5"(h)
TT-1B



TAPE DISPENSER
Mounts to accessory railing or can attach to the top of a cart with included suction cups.
TD-1



WASTE BASKET
Mounts to accessory rail.
SMALL (6 quart)
9"(w) x 6.25"(d) x 10"(h)
BSK-1

LARGE (20 quart)
13"(w) x 8.5"(d) x 15"(h)
BSK-2 (pictured)



WIRE NOTEBOOK MOUNT
Mounts to accessory railing.
12.5"(w) x 2.5"(d) x 10"(h)
WNH-1



CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Score the plastic strips and snap them to create your own divider system.
Includes: (6) 22" long strips and 8 clips
Additional clips available in a set of 6. (To fit into Steel Junior carts cut strips down to 16.5" long)
DIV-CC2 (Dividers & Clips)
DIV-CLIP (Clips)



LARGE CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Includes: 3 front to back panels, 2 side to side panels and 12 clips.
Additional panels are available to create more compartments if needed. (To fit into a Steel Junior cart cut strips down to 16.5" long)
DIV-UC6 (6" Drawers)
DIV-UC9 (9" & 12" Drawers) (pictured)



DIVIDER KIT FOR 3" ALUMINUM STANDARD DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
DIV-AU3

Front to back divider
Y700286A1
Side to side divider
Y700287A1



DIVIDERS FOR 6", 9" & 12" ALUMINUM STANDARD DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
DIV-AU6

Front to back divider
Y700286A2
Side to side divider
Y700287A2

DIVIDERS & TRAYS

ANESTHESIA CARTS

ELITE ALUMINUM
ISOLATION CARTS

Access to isolation supplies has never been easier or more efficient! Our isolation carts offer storage for caps, gowns and masks. Further accessorize your cart with accessories such as glove box holders, waste baskets, and coat hooks. They also include dual slide out side shelves, that virtually double your work surface, and metal drawer dividers in every drawer.

HOW ARE YOU
PROTECTING
YOU AND YOUR
PATIENTS FROM
INFECTION?

TRADITIONAL STEEL
ISOLATION CARTS


Our economical steel isolation carts offer a robust work platform that can be configured with accessories to accommodate your isolation supplies. If you are looking at adding the optional plastic trays/dividers, you may wish to consider our aluminum carts. The included dividers and dual pull out side shelves that come standard on our aluminum carts may offset the additional cost to upgrade.



USGKA-399-YEL



FEATURE	ELITE ALUMINUM	COMPARISON	TRADITIONAL STEEL
Unicart Construction	●	●	Unicart Construction
Self-closing Ball Bearing Slides	●	●	Self-closing Ball Bearing Slides
Easy-roll Casters	●	●	Easy-roll Casters
Dual Pullout Work Shelves	●	●	Single Pullout Work Shelf
Universal Accessory Mounting System	●	○	Universal Accessory Mounting System
Includes Drawer Divider Set	●	○	Includes Drawer Divider Set
Lightweight Construction	●	○	Lightweight Construction



USGKU-399-YEL

ISOLATION CARTS


Waterloo carts come in both standard and narrow widths along with four different cart heights so you can create the perfect balance between cart size and storage space!



ISOLATION STATIONS

Our Isolation Stations are perfect anywhere a wall/door mounted storage solution for personal protective supplies is needed.

- Pre-drilled keyholes for easy wall mounting
- Door hanger cut out on back of organizer (ISOHANG-2C)
- Made of lightweight yet durable aluminum
- Lifetime Warranty



RECOMMENDED ACCESSORIES

ISOLATION CARTS (FOR MORE ACCESSORIES SEE PAGE 40)



GLOVE & MASK BOX HOLDER

Adjusts to many sizes of boxes. Mounts to accessory rail.

GMH-1



GLOVE BOX HOLDER

Available in 3 sizes to fit any facilities needs.

TWO Box Size - 12"(w) x 4"(d) x 10"(h)

BH-2

THREE Box Size - 18"(w) x 4"(d) x 10"(h)

BH-3 (pictured)



CHART HOLDER

Mounts to accessory rail.

SMALL (6 quart)

9"(w) x 6.25"(d) x 10"(h)

CHT-1



WASTE BASKET

Mounts to accessory rail.

SMALL (6 quart)

9"(w) x 6.25"(d) x 10"(h)

BSK-1

LARGE (20 quart)

13"(w) x 8.5"(d) x 15"(h)

BSK-2 (pictured)



SHARPS CONTAINER MOUNTING BRACKET

Mounts to side of all Standard Carts, Junior and Medtrix. Designed for 5 quart container.

NON-LOCKING

SB-2



FLUID/EQUIPMENT STORAGE TRAYS

For Steel, Junior or Medtrix/Halftrix carts. Includes two trays.

Each tray measures:
17.52"(w) x 4.04"(d) x 4.5"(h)

FT-1 (Steel Carts)

FTA-1 (Aluminum Carts)



SELF-ADHESIVE CHROME HOOKS

Set of two self-adhesive chrome hooks that mount easily on any cart.

ABH-2



SELF-ADHESIVE LABEL HOLDERS

1"(h) x 3"(w), package of 100.

Holders are loaded from the side.

LH-1



3" STANDARD CART DIVIDER TRAY SYSTEM

Includes: White tray with adjustable clear sub-dividers.

Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.

21.5"(w) x 15.25"(d) x 2.3"(h)

DIVTRAY-STD



5" STANDARD DIVIDER TRAY SYSTEM

White tray with adjustable clear sub-dividers.

Fits in Steel and Aluminum full-size Standard only. Perfect for 6" drawers or larger.

21.5"(w) x 15.25"(d) x 4.5"(h)

DIVTRAY-STD5



3" GRAY PLASTIC TRAY WITH ADJUSTABLE DIVIDERS

Includes: 4 short dividers and 4 long dividers.

Fits 2 per drawer in Steel and Aluminum Standard Carts. Fits 1 per drawer in a Steel or Aluminum Junior Cart.

9.25"(w) x 15"(d) x 2.5"(h)

WDC-2025

Sold by the each



5" GRAY PLASTIC TRAY WITH ADJUSTABLE DIVIDERS

Includes: 4 short dividers and 4 long dividers. Perfect for 6", 9", 12" drawers.

Fits 2 per drawer in Steel and Aluminum Standard Carts. Fits 1 per drawer in a Steel or Aluminum Junior Cart. Sold by the each.

9.25"(w) x 15"(d) x 4.5"(h)

WDC-2050

Sold by the each



CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS

Score the plastic strips and snap them to create your own divider system.

Includes: (6) 22" long strips and 8 clips
Additional clips available in a set of 6. (To fit into Steel Junior carts cut strips down to 16.5" long)

DIV-CC2 (Dividers & Clips)

DIV-CLIP (Clips)



LARGE CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS

Includes: 3 front to back panels, 2 side to side panels and 12 clips.

Additional panels are available to create more compartments if needed. (To fit into a Steel Junior cart cut strips down to 16.5" long)

DIV-UC6 (6" Drawers)

DIV-UC9 (9" & 12" Drawers) (pictured)



DIVIDER KIT FOR 3" ALUMINUM STANDARD DRAWERS

Includes: 1 front-to-back and 1 side-to-side divider.

DIV-AU3

Front to back divider

Y700286A1

Side to side divider

Y700287A1



DIVIDERS FOR 6", 9" & 12" ALUMINUM STANDARD DRAWERS

Includes: 1 front-to-back and 1 side-to-side divider.

DIV-AU6

Front to back divider

Y700286A2

Side to side divider

Y700287A2

DIVIDERS & TRAYS

ISOLATION CARTS

ELITE ALUMINUM
BEDSIDE/NURSE SERVER CARTS

Improve your nursing staff's time spent with patient by keeping commonly used patient care/nursing supplies at the bedside. Waterloo aluminum bedside carts can fit into tight patient room areas and even under existing counters. Whatever your storage needs are, Waterloo has a bedside cart designed to fit your available space.

ARE YOU TREATING
YOUR PATIENT
OR HUNTING FOR
SUPPLIES?

TRADITIONAL STEEL
BEDSIDE/NURSE SERVER CARTS

Economical steel bedside carts offer a robust work platform that can be configured with accessories to accommodate your bedside cart supplies. If you are looking at adding the optional plastic trays/dividers, you may wish to consider our aluminum carts. The included dividers and dual pull out side shelves that come standard on our aluminum carts may offset the additional cost to upgrade.



FEATURE	ELITE ALUMINUM	COMPARISON	TRADITIONAL STEEL
Unicart Construction	●	●	Unicart Construction
Self-closing Ball Bearing Slides	●	●	Self-closing Ball Bearing Slides
Easy-roll Casters	●	●	Easy-roll Casters
Dual Pullout Work Shelves	●	●	Single Pullout Work Shelf
Universal Accessory Mounting System	●	○	Universal Accessory Mounting System
Includes Drawer Divider Set	●	○	Includes Drawer Divider Set
Lightweight Construction	●	○	Lightweight Construction



BEDSIDE/NURSE SERVER CARTS

Waterloo bedside carts are available in a wide variety of sizes and configurations. If your bedside storage requirements are more than our Junior carts accommodate, consider one of our Standard size carts which provide a larger storage capacity.

Keyless locking options for our bedside carts include: mechanical push button, electronic keypad and even proximity card readers for swipe in/swipe out convenience!



BEDSIDE/NURSE SERVER CARTS

All carts with a "Low profile" package include 3" casters mounted directly under the cart (no bumpers) and are designed to tuck neatly under an existing counter or desk. Providing the perfect storage solution for an unused undercounter alcove!

Optional heavy duty thermofoil cart tops are also available
(Call for details).



RECOMMENDED ACCESSORIES

BEDSIDE CARTS (FOR MORE ACCESSORIES SEE PAGE 40)



ACCESSORY RAIL FOR STEEL CARTS
For Steel Standard or Junior Carts ONLY.
[SR-1 \(Side Mount\)](#) (pictured)
[CB-1 \(Back Mount\)](#)



GLOVE & MASK BOX HOLDER
Adjusts to many sizes of boxes. Mounts to accessory rail.
[GMH-1](#)



GLOVE BOX HOLDER
Available in 3 sizes to fit any facilities needs.
[TWO](#) Box Size - 12"(w) x 4"(d) x 10"(h)
[BH-2](#)
[THREE](#) Box Size - 18"(w) x 4"(d) x 10"(h)
[BH-3](#) (pictured)



CHART HOLDER
Mounts to accessory rail.
14.25"(w) x 3.75"(d) x 10"(h)
[CHT-1](#)



JUNIOR CART AMPULE HOLDER DIVIDER TRAY SYSTEM
Includes: White plastic tray, 12 clear ampule dividers and 4 clear sub-dividers.
For Steel and Aluminum Junior Carts ONLY. Fits 1 per drawer. Perfect for 3" drawers.
15.25"(w) x 15.25"(d) x 2.3"(h)
[WMT-4](#)



JUNIOR CART DIVIDER TRAY SYSTEM
Includes: White plastic tray and 4 adjustable clear sub-dividers.
For Steel and Aluminum Junior Carts ONLY. Fits 1 per drawer. Perfect for 3" drawers.
15.25"(w) x 15.25"(d) x 2.3"(h)
[WMT-5](#)



CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Score the plastic strips and snap them to create your own divider system.
Includes: (6) 22" long strips and 8 clips
Additional clips available in a set of 6. (To fit into Steel Junior carts cut strips down to 16.5" long)
[DIV-CC2 \(Dividers & Clips\)](#)
[DIV-CLIP \(Clips\)](#)



HALF DRAWER NARCOTICS BOX
Keyless mechanical push-button narcotics box can be used to safely secure patient's personal items.
[SD-3S](#)



FLUID/EQUIPMENT STORAGE TRAYS
For Steel, Junior or Medtrix/Halftr carts. Includes two trays.
Each tray measures:
17.52"(w) x 4.04"(d) x 4.5"(h)
[FT-1 \(Steel Carts\)](#)
[FTA-1 \(Aluminum Carts\)](#)



UTILITY RECEPTACLES
Mounts to accessory rail.
5"(w) x 5"(d) x 7"(h)
[UR-2 \(1.2 Liters\)](#)
7.25"(w) x 7.25"(d) x 9"(h)
[UR-3 \(3.5 Liters\)](#)



SELF-ADHESIVE CHROME HOOKS
Set of two self-adhesive chrome hooks that mount easily on any cart.
[ABH-2](#)



LARGE CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Includes: 3 front to back panels, 2 side to side panels and 12 clips.
Additional panels are available to create more compartments if needed. (To fit into a Steel Junior cart cut strips down to 16.5" long)
[DIV-UC6](#) (6" Drawers)
[DIV-UC9](#) (9" & 12" Drawers) (pictured)



DIVIDER KIT FOR 3" ALUMINUM JUNIOR DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
[DIV-AJ3](#)
Front to back divider
[Y700286A1](#)
Side to side divider
[Y700287A5](#)



DIVIDERS FOR 6", 9" & 12" ALUMINUM JUNIOR DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
[DIV-AJ6](#)
Front to back divider
[Y700286A2](#)
Side to side divider
[Y700287A6](#)

DIVIDERS & TRAYS

BEDSIDE CARTS

SPECIALTY CARTS

CUSTOM CONFIGURED CARTS

By combining our wide array of accessories you can custom configure your own specialty cart. Looking for a "Difficult airway cart"? Simply add a scope holder and a waste basket. Would you like an "IV Cart"? No problem, just add a glove box holder, waste basket and a sharps container bracket.

Give us a call and our customer service team can assist you!

HERE ARE A FEW SPECIALTY CART EXAMPLES



CAST CART
UTGKA-33669-GRA
CAST-1 Casting Shelf

DIFFICULT AIRWAY CART
UTGKA-333369-LTB
SH-2 Scope Holder
BSK-2 Waste Basket

SUTURE CART
UMGKA-33369-LTG
SUT-1 Suture Shelf
AMP-2 Back Accessory Panel
BSK-2 Waste Basket
SB-1 Univ Sharps Holder



IV CART
JSGKA-3369-WHT
SB-1 Univ Sharps Holder
BSK-2 Waste Basket
BH-2 Glove box holder
AMP-1 Side panel
IV-2 IV Pole

MALIGNANT HYPERTHERMIA CART
UXGKA-36H-DKB (Cart Only)
79-2105 (Fridge Only)
ACCUCOLD 1 cu. ft. interior refrigerator with an external microprocessor temperature controller. Temperature alarm sounds if the unit's interior goes out of range. Key locking front door and hospital grade plug come standard.

NOTE: Must order both part numbers to build a hyperthermia cart

WXGKU-333669-FWG
Overall Dims: 36"(w) x 24.5"(d) x 45"(h)
Drawer Dims: 29"(w) x 17"(d)
Available ONLY in X-Tall, Steel with
FWG or LTG drawer front colors

SPECIALTY CARTS

CUSTOM CONFIGURED CARTS

SPECIALTY CARTS

PHLEBOTOMY CARTS

This very popular cart offers a variety of storage compartments and tilt bins for all your phlebotomy supplies. It was purpose built to allow phlebotomists the opportunity to spend more time on patient floors and less time heading back to the lab for fresh supplies.



SHORT PHLEBOTOMY CART
MSWA-3469-WHT
Internal sharps container, internal trash container, tilt bins, drawers with removable divider trays
Accommodates a 12 QT Kendall 85221R sharps container (not included)



TALL PHLEBOTOMY CART
MTWA-34696-WHT
Internal sharps container, internal trash container, tilt bins, drawers with removable divider trays
Accommodates a 12 QT Kendall 85221R sharps container (not included)

ACCESSORIES

DIVIDERS & TRAYS



3" GRAY PLASTIC TRAY WITH ADJUSTABLE DIVIDERS
Includes: 4 short dividers and 4 long dividers.
Fits 2 per drawer in Steel and Aluminum Standard Carts. Fits 1 per drawer in a Steel or Aluminum Junior Cart.
9.25"(w) x 15"(d) x 2.5"(h)
[WDC-2025](#)
Sold by the each



5" GRAY PLASTIC TRAY WITH ADJUSTABLE DIVIDERS
Includes: 4 short dividers and 4 long dividers. Perfect for 6", 9", 12" drawers.
Fits 2 per drawer in Steel and Aluminum Standard Carts. Fits 1 per drawer in a Steel or Aluminum Junior Cart. Sold by the each.
9.25"(w) x 15"(d) x 4.5"(h)
[WDC-2050](#)
Sold by the each



3" WHITE PLASTIC TRAY WITH 13 BINS
Can be easily reconfigured.
Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
22"(w) x 15"(d) x 2.5"(h)
[TT-1A](#)



3" WHITE PLASTIC TRAY WITH 10 BINS
Can be easily reconfigured.
Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
22"(w) x 15"(d) x 2.5"(h)
[TT-1B](#)



3" STANDARD CART DIVIDER TRAY SYSTEM
Includes: White tray with adjustable clear sub-dividers.
Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
21.5"(w) x 15.25"(d) x 2.3"(h)
[DIVTRAY-STD](#) (pictured)
Also, available in 5" high version.
21.5"(w) x 15.25"(d) x 4.5"(h)
[DIVTRAY-STD5](#)



3" COMPLETE AMPULE HOLDER DIVIDER TRAY SYSTEM
Includes: 3 rows of 17 adjustable clear ampule dividers.
Fits in Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
[WMT-3](#)



JUNIOR CART AMPULE HOLDER DIVIDER TRAY SYSTEM
Includes: White plastic tray, 12 clear ampule dividers and 4 clear sub-dividers.
For Steel and Aluminum Junior Carts ONLY. Fits 1 per drawer. Perfect for 3" drawers.
15.25"(w) x 15.25"(d) x 2.3"(h)
[WMT-4](#)



JUNIOR CART DIVIDER TRAY SYSTEM
Includes: White plastic tray and 4 adjustable clear sub-dividers.
For Steel and Aluminum Junior Carts ONLY. Fits 1 per drawer. Perfect for 3" drawers.
15.25"(w) x 15.25"(d) x 2.3"(h)
[WMT-5](#)



3" AMPULE HOLDER DIVIDER TRAY SYSTEM
Clear dividers are adjustable.
Fits into Steel and Aluminum Standard Carts ONLY. Perfect for 3" drawers.
21.5"(w) x 15.25"(d) x 2.3"(h)
[DIVTRAY-AMP](#)



3" ANESTHESIA DIVIDER TRAY SYSTEM
Includes: White plastic tray with dividers.
Fits in Steel Standard Carts ONLY. Perfect for 3" drawers.
21.5"(w) x 15.25"(d) x 2.3"(h)
[DIVTRAY-ANES](#)



WHITE DIVIDER BOX WITH CLEAR LID
Includes: White plastic tray with clear lid. Can be sealed with clear plastic seal. Steel and Aluminum Standard Carts ONLY.
21.75"(w) x 15.75"(d) x 2.5"(h)
[EMER-BOX](#)

NOTE: EMER-DV is available as a divider set for the EMER-BOX



SEALABLE, PERFORATED BAGS
Fits over divider trays.
Package of 20.

23"(w) x 17"(d) x 5.5"(h)
[BAG-3](#)
23"(w) x 17"(d) x 7"(h)
[BAG-6](#) (pictured)



5" ADJUSTABLE, CLEAR DIVIDER SYSTEM WITH SUB-DIVIDERS
For full-size Steel Standard Carts ONLY.
21.5"(w) x 15.25"(d) x 5"(h)
[DIV-DRW5](#)

NOTE: To purchase this divider for the ER2000 Instant Access Cart, please order [DIV-DRW5ER](#)




3" ADJUSTABLE CLEAR DIVIDER SYSTEM WITH ROW OF AMPULE HOLDERS
For Steel Standard Carts ONLY.
21.5"(w) x 15.25"(d) x 3"(h)
[DIV-DRW12](#)

NOTE: To purchase this divider for the ER2000 Instant Access Cart, please order [DIV-DRW12ER](#)



3" ADJUSTABLE, CLEAR DIVIDER SYSTEM WITH SUB-DIVIDERS
For Steel Standard Carts ONLY.
21.5"(w) x 15.25"(d) x 3"(h)
[DIV-DRW34](#)

NOTE: To purchase this divider for the ER2000 Instant Access Cart, please order [DIV-DRW34ER](#)



LONG OR SHORT, CLEAR PLASTIC DIVIDERS
For use in a Waterloo white plastic divider tray.
[DIV-10L](#) (Small Sub Divider)
[DIV-10S](#) (Small Ampule Divider)




CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS (3"H DRAWERS)
Score the plastic strips and snap them to create your own divider system.
Includes: (6) 22" long strips and 8 clips
Additional clips available in a set of 6. (To fit into Steel Junior carts cut strips down to 16.5" long)
[DIV-CC2](#) (Dividers & Clips)
[DIV-CLIP](#) (Clips)




LARGE CUT & CLIP, FRONT TO BACK DIVIDER SYSTEM FOR STEEL CARTS
Includes: 3 front to back panels, 2 side to side panels and 12 clips
Additional panels are available to create more compartments if needed. (To fit into a Steel Junior cart cut strips down to 16.5" long)
[DIV-UC6](#) (6" Drawers)
[DIV-UC9](#) (9" & 12" Drawers) (pictured)



CLEAR DIVIDABLE TRAYS
Trays fit in Standard drawers and Instant Access Series ER2000.
[ER-TRAY3](#)
[ER-TRAY6](#)
For dividing systems, purchase [DIV-DRW12ER](#), [DIV-DRW34ER](#) or [DIV-DRW5ER](#)



SMALL SUB DIVIDER
Clear Sub Divider with Label Holder for 3" & 4" drawers. For use with metal dividers in our aluminum carts.
[M14089](#)



LARGE SUB DIVIDER
Clear Sub Divider with Label Holder for 6" , 9" , 11" & 12" drawers. For use with metal dividers in our aluminum carts.
[M14089A1](#)



DIVIDER SET FOR 3" ALUMINUM STANDARD DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
[DIV-AU3](#)

Front to back divider
[Y700286A1](#)
Side to side divider
[Y700287A1](#)



DIVIDER SET FOR 6", 9" & 12" ALUMINUM STANDARD DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
[DIV-AU6](#)

Front to back divider
[Y700286A2](#)
Side to side divider
[Y700287A2](#)



DIVIDER SET FOR 3" ALUMINUM JUNIOR DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
[DIV-AJ3](#)

Front to back divider
[Y700286A1](#)
Side to side divider
[Y700287A5](#)



DIVIDER SET FOR 6", 9" & 12" ALUMINUM JUNIOR DRAWERS
Includes: 1 front-to-back and 1 side-to-side divider.
[DIV-AJ6](#)

Front to back divider
[Y700286A2](#)
Side to side divider
[Y700287A6](#)

ACCESSORIES

UNICART ACCESSORIES



**LIGHTWEIGHT BACKSIDE
CARDIAC BOARD MOUNT**
Includes plastic cardiac board and
metal mounting brackets for backside
of Standard Carts.
24"(w) x 1"(d) x 18"(h)
[CBB-6 \(Steel Carts\)](#)
[CBB-7 \(Aluminum Carts\)](#)



**LIGHTWEIGHT CARDIAC
BOARD MOUNT**
Includes plastic cardiac board and side
mounting brackets for backside of
Aluminum Standard and Junior Carts.
24"(w) x 1"(d) x 18"(h)
[CBB-8](#)



**THIN CARDIAC BOARD
BACKSIDE MOUNT**
Includes cardiac board and brackets
for backside of Steel or Aluminum
Standard Carts.
21"(w) x 0.5(d) x 15.25"(h)
[CBB-9](#)



**ER2000 BACKSIDE
CARDIAC BOARD MOUNT**
Includes plastic cardiac board and
mounting bracket for backside of
ER2000 cart.
[CBB-ER](#)



OXYGEN TANK MOUNT
Plastic brackets with adjustable
Velcro straps. Can be mounted to all
Standard and Junior Carts.
[OTB-4 \(pictured\)](#)
[O-2 \(ER2000\)](#)



**METAL OXYGEN
TANK BRACKET**
Can be mounted to all Standard and
Junior Carts.
[OTB-5](#)



**OXYGEN TANK
CENTER MOUNTING RAIL**
Includes 2 mounting rails. Allows
OTB-5 to be center mounted on
Steel carts ONLY.
[AMP-1S](#)
(OTB-5 sold separately)



IV POLE
Chrome IV pole and mounting bracket.
Mounts to side of all Standard Carts,
Junior and Medtrix.
[IV-2](#)



**SHARPS CONTAINER
MOUNTING BRACKET**
With adjustable velcro strap.
Mounts to side of all Unicarts, Med
Jr, ER-2000 and Medtrix.
[SB-1](#)



**SHARPS CONTAINER
MOUNTING BRACKET**
Mounts to side of all Standard Carts,
Junior and Medtrix.
Accommodates a 5QT Kendall 85131
sharps container (not included)

NON-LOCKING
[SB-2](#)
LOCKING [\(pictured\)](#)
[SB-2L](#)



**SELF-ADHESIVE
CHROME HOOKS**
Set of two self-adhesive chrome
hooks that mount easily on any cart.
[ABH-2](#)



WASTE BASKET
Mounts to accessory rail.
SMALL (6 quart)
9"(w) x 6.25"(d) x 10"(h)
[BSK-1](#)

LARGE (20 quart)
13"(w) x 8.5"(d) x 15"(h)
[BSK-2 \(pictured\)](#)



GLOVE BOX HOLDER
Available in 2 sizes to fit any facilities
needs.

TWO Box Size - 12"(w) x 4"(d) x 10"(h)
[BH-2](#)
THREE Box Size - 18"(w) x 4"(d) x 10"(h)
[BH-3 \(pictured\)](#)



**HORIZONTAL GLOVE
BOX HOLDER**
Holds 3 standard sized glove boxes.
[BH-3H](#)



**VERTICAL GLOVE
BOX HOLDER**
Holds 3 standard sized glove boxes.
[BH-3V](#)



**GLOVE & MASK
BOX HOLDER**
Adjusts to many sizes of boxes.
Mounts to accessory rail.
[GMH-1](#)



TILT BINS
Clear plastic tilt bins with white
exterior. Mounts to base rail
system (BRS-2, BRS-3 OR BRS-JR).
Available in 3, 4, 5, or 6 bin units.

Shown from top to bottom:
[TB-6](#)
Overall dims: 23.75"W x 4.5"H x 4.5"D
Inside Dims: 3.5"W x 4"H x 3"D
[TB-3](#)
Overall dims: 24.63"W x 9.5"H x 8.5"D
Inside Dims: 7"W x 8.5"H x 5"D
[TB-5](#)
Overall dims: 23.75"W x 6.5"H x 4.5"D
Inside Dims: 6"W x 4.5"H x 4"D
[TB-4](#)
Overall dims: 23.6"W x 8.25"H x 7.25"D
Inside Dims: 5"W x 5"H x 5"D



UTILITY RECEPTACLES
Mounts to accessory rail.

5"(w) x 5"(d) x 7"(h)
[UR-2 \(1.2 Liters\)](#)

7.25"(w) x 7.25"(d) x 9"(h)
[UR-3 \(3.5 Liters\)](#)



CHART HOLDER
Mounts to accessory rail.
14.25"(w) x 3.75"(d) x 10"(h)
[CHT-1](#)



NOTEBOOK MOUNT
Mounts to accessory railing.
12.5"(w) x 2.5"(d) x 10"(h)
[WNH-1](#)



**EXTRA-LARGE
WIRE BASKET**
Mounts to accessory railing.
15.5"(w) x 9.5"(d) x 5.5"(h)
[WSB-1](#)



**SIDE MOUNT
BUCKET HOLDER**
Holds (2) 165oz cast buckets.
17.5"(w) x 6.5"(d) x 8.8"(h)
[CAST-1](#)

ACCESSORIES

UNICART ACCESSORIES



ACCESSORY SIDE MOUNTING PANEL
For Standard Carts ONLY.
Mounts to SIDE of cart.
[AMP-1](#)



ACCESSORY BACK MOUNTING PANEL
For Aluminum Standard Carts ONLY.
Mounts to BACKSIDE of cart.
[AMP-2](#)



ACCESSORY RAIL FOR STEEL CARTS
17.75" Long

[SR-1 \(Side Mount\)](#) *(pictured)*
Fits Steel Standard or Junior Carts ONLY.
[CB-1 \(Back Mount\)](#)
Fits Steel Standard Carts ONLY.



HORIZONTAL RAIL SYSTEM
Includes two vertical and three horizontal rails.
[BRS-2 \(Standard Carts\)](#) *(pictured)*
[BRS-JR \(Junior Carts\)](#)
[BRS-2W \(Wide Carts\)](#)

(Shown on aluminum cart)



FLAT SHELF
Mounts to horizontal rail system (BRS-2).
[SHLF-2 \(Steel Standard Carts\)](#)
24"(w) x 10"(d) x .75"(h)
[SHLF-2W \(Steel Wide Carts\)](#)
35.5"(w) x 10"(d) x .75"(h)
[SHLF-3 \(Aluminum Standard Carts\)](#)
24"(w) x 10"(d) x .75"(h) *(pictured)*



WRAP-AROUND ALUMINUM RAILING
[R-3 \(Steel Standard Carts\)](#) *(pictured)*
[R-4 \(Aluminum Standard Carts\)](#)
[R-JR \(Steel Junior Carts\)](#)
[R-JRA \(Aluminum Junior Carts\)](#)



SWIVELING MONITOR SHELF
Includes an adjustable strap to secure items in place. Maximum weight capacity 25 lbs. Can be locked in place.
[ER-SHLF \(Steel Carts\)](#)
[ER-SHLFA2 \(Aluminum Carts\)](#) *(pictured)*
[ER-SHLF2000 \(ER2000\)](#)



SUTURE SHELF
Includes (2) vertical rails, (1) top horizontal rail and (2) suture shelves. Aluminum carts require AMP-2. Each shelf measures:
21.5"(w) x 5.5"(d) x 3.25"(h)
[SUT-1](#)



SUCTION SHELF
For Standard or Junior Carts. Fits most brands of portable suction pumps.
[SS-M](#) *(pictured)*
18"(w) x 8"(d) x 3"(h)

[SS-M-ER2000 \(ER2000\)](#)
11"(w) x 7.5"(d) x 2.25"(h)



FLUID/EQUIPMENT STORAGE TRAYS
For Steel, Junior or Medtrix/Halftrix carts. Includes two trays. Each tray measures:
17.52"(w) x 4.04"(d) x 4.5"(h)
[FT-1 \(Steel Carts\)](#)
[FTA-1 \(Aluminum Carts\)](#)



STEEL SCOPE HOLDER
Easily secure to the wall or our Unicarts. Locking Plexiglass doors also included for added security.
Available in two sizes:
[SH-2 \(TWO Scope Holder\)](#)
[SH-4 \(FOUR Scope Holder\)](#)



CATHETER HOLDER
For Standard or Junior Carts.
7.5"(w) x 2"(d) x 23"(h)
[CH-3](#)



CATHETER HOLDER
Mounts to Aluminum Standard Carts ONLY. Mounting to backside of cart requires AMP-2 *(sold separately)*.
4"(w) x 2"(d) x 19"(h)
[CH-2](#)



WIRE BASKET
Mounts to accessory railing.
SMALL
8.5"(w) x 5.5"(d) x 3.5"(h)
[WB-1S](#)

LARGE
15.5"(w) x 5.5"(d) x 3.5"(h)
[WB-1L](#)



CLEAR PLASTIC ORGANIZER
SMALL contains 3 compartments. Mounts to accessory railing.
9"(w) x 5"(d) x 3"(h)
[ORG-1S](#)



TAPE DISPENSER
Mounts to accessory railing or can attach to the top of a cart with included suction cups.
[TD-1](#)



CORD WRAP
Commonly used with outlet strip *(sold separately)*.
[CW-1](#)



OUTLET POWER STRIP
Hospital grade outlet strip.
• 15 foot long cord
• 6 covered outlets
• 15A circuit breaker
• UL 1363A rated
[OS-1](#)



PLASTIC CART COVER
For Standard Carts. Includes 5 bags.
Dimensions: 33"(w) x 23"(d) x 42"(h)
[BAG-COVER5](#)

ACCESSORIES

UNICART ACCESSORIES

ACCESSORIES

UNICART ACCESSORIES



**PLASTIC
BREAKAWAY SEALS**
Can be used with lever locking
systems, breakaway gate locks
and EMER-BOX divider trays.
[PS-100](#)



PROXIMITY CARDS

HID PROXIMITY CARD
[PROX-1](#) *(pictured)*

PROXIMITY TAG
1" diameter
[PROX-2](#) *(pictured)*

ICLASS PROXIMITY CAR
[PROX-3](#)



**SELF-ADHESIVE
LABEL HOLDERS**
1"(h) x 3"(w), package of 100.
Holders are loaded from the side.
[LH-1](#)



LOCKING SECURITY BOX
For use of transportation of narcotics
or permanently mounted inside
drawer. Includes two keys.

15.6"(w) x 11.4"(d) x 2.5"(h)
[SD-3](#)



**3" PUSHBUTTON LOCKING
SECURITY BOX**
For use in transportation of narcotics.
Permanently mounted.
Includes override lock with 2 keys.
10.6"(w) x 15.75"(d) x 3"(h)
[SD-3S](#)

CROSS REFERENCE SPECIFICATION CHART

CART FRAME SPECIFICATIONS					DRAWER SPECIFICATIONS				
FRAME STYLE	OVERALL DIMENSIONS	WEIGHT RANGE (lbs)	EXTERIOR COLORS	PULL-OUT SHELF	DRAWER INSIDE DIMENSIONS	DRAWER SIZES (inches high)	VERTICAL DRAWER SPACE	BUILT-IN DIVIDERS	DRAWER COLORS
ALUMINUM STANDARD									
Short	36"h x 29"w x 24.5"d	65-90	LTG, RED, ELB	Dual	22"w x 16.5"d	3, 6, 9, 12	21"	Yes	16 Color Options
Medium	39"h x 29"w x 24.5"d	70-95	LTG, RED, ELB	Dual	22"w x 16.5"d	3, 6, 9, 12	24"	Yes	16 Color Options
Tall	42"h x 29"w x 24.5"d	75-95	LTG, RED, ELB	Dual	22"w x 16.5"d	3, 6, 9, 12	27"	Yes	16 Color Options
X-Tall	45"h x 29"w x 24.5"d	80-100	LTG, RED, ELB	Dual	22"w x 16.5"d	3, 6, 9, 12	30"	Yes	16 Color Options
STEEL STANDARD									
Short	36"h x 29"w x 24.5"d	85-100	LTG, RED, ELB	Single	22"w x 16.5"d	3, 6, 9, 12	21"	No	16 Color Options
Medium	39"h x 29"w x 24.5"d	85-100	LTG, RED, ELB	Single	22"w x 16.5"d	3, 6, 9, 12	24"	No	16 Color Options
Tall	42"h x 29"w x 24.5"d	90-120	LTG, RED, ELB	Single	22"w x 16.5"d	3, 6, 9, 12	27"	No	16 Color Options
X-Tall	45"h x 29"w x 24.5"d	100-140	LTG, RED, ELB	Single	22"w x 16.5"d	3, 6, 9, 12	30"	No	16 Color Options
ALUMINUM JUNIOR									
Short	36"h x 23"w x 24.5"d	60-80	LTG, 3 wood Options	Dual	16.5"w x 16.5"d	3, 6, 9, 12	21"	Yes	16 Color Options 3 Wood Options
Medium	39"h x 23"w x 24.5"d	60-85	LTG, 3 wood Options	Dual	16.5"w x 16.5"d	3, 6, 9, 12	24"	Yes	16 Color Options 3 Wood Options
Tall	42"h x 23"w x 24.5"d	65-85	LTG, 3 wood Options	Dual	16.5"w x 16.5"d	3, 6, 9, 12	27"	Yes	16 Color Options 3 Wood Options
X-Tall	45"h x 23"w x 24.5"d	70-90	LTG, 3 wood Options	Dual	16.5"w x 16.5"d	3, 6, 9, 12	30"	Yes	16 Color Options 3 Wood Options
STEEL JUNIOR									
Short	36"h x 23"w x 24.5"d	80-95	LTG, 3 wood Options	Single	16.5"w x 16.5"d	3, 6, 9, 12	21"	No	16 Color Options 3 Wood Options
Medium	39"h x 23"w x 24.5"d	80-95	LTG, 3 wood Options	Single	16.5"w x 16.5"d	3, 6, 9, 12	24"	No	16 Color Options 3 Wood Options
Tall	42"h x 23"w x 24.5"d	85-115	LTG, 3 wood Options	Single	16.5"w x 16.5"d	3, 6, 9, 12	27"	No	16 Color Options 3 Wood Options
X-Tall	45"h x 23"w x 24.5"d	90-130	LTG, 3 wood Options	Single	16.5"w x 16.5"d	3, 6, 9, 12	30"	No	16 Color Options 3 Wood Options



Contract Holder
36F79719D0279
Effective thru September 30, 2024

800.833.4419 direct **602.414.3691**

WaterlooHealthcare.com



WATERLOO HEALTHCARE
3730 East Southern Avenue
Phoenix, Arizona 85040-3958

REV0222