



ROMET®

ELECTRONIC LARYNX

SINCE 1979



PRODUCT CATALOG

2021



Digital Electronic Larynx Units



ROMET® CLARITY DIGITAL ELECTRONIC LARYNX

REF	DESCRIPTION
-----	-------------

R700	Programmable Digital Electronic Larynx
------	--

- Polished 6061 Aluminum Barrel
- Electroplated gunmetal finish
- Easy to operate
- Volume and Pitch Controls
- Li-ion Permanent Battery
- 30 hours speaking time
- USB rechargeable adapter & cable incl.
- Safety Lanyard
- Weight: 5.0 oz with battery incl.



ROMET® GENESIS DIGITAL ELECTRONIC LARYNX

REF	DESCRIPTION
-----	-------------

R600	Digital Electronic Larynx
------	---------------------------

- Polished 6061 Aluminum Barrel
- Sandblasted black finish
- Economical great as a backup or starter
- Great sound no pitch or volume controls
- Li-ion Permanent Battery
- 30 hours speaking time
- USB rechargeable adapter & cable incl.
- Safety Lanyard
- Weight: 5.0 oz with battery incl.

Analog Electronic Larynx Units



ROMET® R120 ELECTRONIC LARYNX

- Delrin Barrel
- On/Off Switch (No Volume/Pitch Controls)
- 2 Li-ion Batteries
- Smart Charger
- Safety Lanyard
- Weight: 4.4 oz with battery

REF	DESCRIPTION
R120	Electronic Larynx



ROMET® R210 ANALOG ELECTRONIC LARYNX

- Sandblasted Aluminum Barrel
- Volume & Pitch Controls
- 2 Li-ion Batteries
- Smart Charger
- Safety Lanyard
- Weight: 4.7 oz with battery

REF	DESCRIPTION
R210	Analog Electronic Larynx



ROMET® R310 ANALOG ELECTRONIC LARYNX

- Delrin Barrel with Soft Touch Grip
- Volume & Pitch Controls
- 2 Li-ion Batteries
- Smart Charger
- Safety Lanyard
- Weight: 4.5 oz with battery

REF	DESCRIPTION
R310	Analog Electronic Larynx