

HEALTHCARE FURNITURE MANUFACTURER









INDEX

Upright Cabinets		00
Base & wall Cabinets		00
Narcotic Cabinets		00
Work Tables		00
Sink Tables		00
Trolleys		00
Dust Bins		00
Overbed Table		00
Carts		00
Laundry		00
Shelving		00
Beds		00
Infusion Stands		00
Examination Lights		00
Stools		00
Blood Donor Chair		00
Illuminator		00
Ward		00
Lab		00





UPRIGHT CABINET

UPRIGHT CABINET



UC-4HD
Upright Cabinet 4 Hinged Doors



UC-PG20
Upright Cabinet 20 Pegion Holes



UC-3HD
Upright Cabinet 3 Hinged Doors



UC-2HD
Upright Cabinet 2 Hinged Doors

Upright Cabinets

- Completely made of stainless steel 304 grade #4 finish sheets and tubes.
- With shelves.
- Sliding / hinged doors with lock and keys.
- With Center partetion
- On Height adjustable foot studs.
- Welded fabrication."

Options

- Glass Door
- Wheels
- Powder coated cabinets
- Drawers



UC-3GHD

Upright Cabinet 3 Hinged Glass Doors



UC-2SD

Upright Cabinet 2 Solid Sliding Doors.



UC-2GHD

Upright Cabinet 2 Hinged Glass Doors



UC- 4G3D

Upright Cabinet 4 Hinged Glass Doors with 3 Sliding Drawers"

Options

- * Completely made of non magnetic stainless steel sheets, tubes and rods.
- * Doors with lock and keys
- * Height adjustable shelves
- * On height adjustable foot studs.
- * Welded fabrication

- * Additional Shelves
- * Drawers
- * Without doors
- * Heavy duty swivel castors
- * Ventilation as per requirement
- * Skirting

UPRIGHT CABINET

UPRIGHT CABINET



UC-2GHD

UC-2GHD-4D

UC-2HD

UC-2GHD-R

UC-2HD-4D



UC-2GHD-4SD

Upright Cabinet 2 Hinged Glass Doors with 4 Drawers

- Upright Cabinet 2 Hinged Glass Door with 4 drawers.**
- * Completely made of stainless steel 304 grade
 - * #4 finish sheets and rods-tubes
 - * Two nos stainless steel framed glass hinged doors on top with lock and keys
 - * Four nos of side by side sliding drawers at bottom
 - * Three nos of height adjustable shelves.
 - * Welded fabrication



UC-2GSD-W

Upright Cabinet 2 Glass door with on casters"

Options

- * Additional Shelves
- * Drawers
- * Without doors
- * Heavy duty swivel castors
- * Ventilation as per requirement
- * Skirting



UCC-6R-P

Upright Catheter Cabinet with 6 Rails

Catheter Cabinet

- * Hinged glass doors.
- * One storage shelf at top
- * Bottom cabinet with Six nos of sliding rails, each with Five hooks.
- * Completely made of electro galvanized sheets with powder coated oven baked paint finish .
- * On height adjustable foot studs..

UC-1 GHD-P

Upright Cabinet 1 Hinged door with Powder coated"

Options

- Heavy duty swivel castors
- Custom made design available
- Skirting

LOCKERS



STL-1 DP

Locker with Single Door

GI Powder Coated locker.

- * Door with lock and keys.
- * With door handle and smart label holder.
- * With Hanger.
- * With Ventilation louvers.
- * Completely made of Electro galvanized sheets with powder coated oven baked paint finish.



STL-2DP

Locker with Two Door



STL-3DP

Locker with Three Door



STL-2ZP

Z Type Locker

GI Powder Coated Z'-shape locker.

- * Double door with lock and keys.
- * With door handle and smart label holder.
- * With Hanger.
- * With Ventilation louvers.
- * Completely made of Electro galvanized sheets with powder coated oven baked paint finish.

Options

- Stainless steel lockers
- Different colours
- Custom made lockers



- * High pressure compact laminate (HPL/HPCL) Lockers
- * Different type of HPL lockers
- * Integrated lockers
- * Lockers with integrated bench.

Options

- * Different colours
- * Custom made lockers

Advantage of HPL/HACL

- * High Pressure Compact Laminates are manufactured by fusing layers of kraft papers which are impregnated with Phenolic Resins at a high pressure and heat to form a perfect solid product.
- * Resistant to water & moisture
- * Resistant to abrasion & scratches
- * Easy to clean
- * Resistant to fire, heat and cigarette burns
- * Eco friendly, non toxic & anti-static
- * Extremely durable and does not deform

* Single tier locker	HPL-1
* Two tier locker	HPL-2
* Three tier locker	HPL-3
* Four tier locker	HPL-4
* Z Type locker	HPL-Z
* Locker with integrated Bench	HPL-B
* Lockers with steel base	BSB

HPL LOCKERS

BASE & WALL CABINET



BC-HD



BC-SD

Base Cabinet with Sliding Doors



BC-8D

Base Cabinet with 8 Drawers

Stainless steel Base Cabinet.

- * 40mm thick reinforced top with water proof laminated boards, covered with stainless steel sheets.
- * One under shelf and one mid shelf
- * Hinged doors.
- * On height adjustable foot studs.
- * Completely made of stainless steel 304 grade
- * 4 finish sheets and tubes.

Options

- * Corian Or Marble Tops
- * Stainless Steel Framed Glass Doors
- * Lock & Keys
- * Sliding Drawers
- * Skirting
- * Heavy Duty Swivel Castors
- * Back splash



WC- HGD

Wall Cabinet with Hinged Glass Door



WC- SD

Wall Cabinet with Sliding Door



WC- HD

Wall Cabinet with Hinged Doors

Stainless steel Wall Cabinet

- * With hinged doors
- * One adjustable shelf
- * With wall mounting clamps.
- * Completely made of stainless steel 304 grade
- * # 4 finish sheets and tubes.

Options

- * Sliding door
- * Glass Door
- * Drawers
- * Powder coated



FB

First Aid Box

First Aid Box

- * Completely Made Of Stainless Steel 304 Grade
- * # 4 Finish Sheets & Tubes
- * Door With Locks & Keys
- * Knob On Door
- * Inner Shelf
- * Wall Mounting Type
- * with First Aid Sticker



NC- ELD

Narcotic Cabinet with Number Lock



NC- 315D

Narcotic Cabinet with Double Door

Narcotic Cabinet

- * Wall Mounted type
- * Completely made of Stainless steel 304 grade #4 finish sheets and tubes.
- * Single hinged door with two nos. of lock and keys.
- * Two inner shelves.
- * Inner safety cabinet with lock and key.
- * 220V Buzzer, light and door switch.

Options

- Electronic Number Lock System
- Inner Shelves
- Additional Inner Door
- Custom made design



NC- 315S

Narcotic Cabinet with Single Door



NC- UELD

Upright Narcotic Cabinet with Double Doors and Electronic Number Lock

Upright narcotic cabinet.

- * Free standing type
- * Double hinged doors and electronic number lock.
- * Two nos adjustable inner shelves.
- * Inner safety cabinet with lock and key.
- * 220V Buzzer, light and door switch.
- * Completely made of non magnetic stainless steel sheets and tubes.
- * On height adjustable foot studs with/without bottom skirting



NC- USD

Upright Narcotic Cabinet with Single Door

Options

- * Electronic Number Lock System
- * Inner Shelves
- * Additional Inner Door
- * Free Standing or with Wheels
- * Custom made design

NARCOTIC CABINETS



NC- SDD

Narcotic Cabinet with Full Inner Door



NC- UD-2LP

Upright Narcotic Cabinet with Single Door & Locks (Powder Coated)



NC- DDP

Narcotic Cabinet with Full Inner Door- Powder Coated

Narcotic cabinet.-Powder Coated

- * Free standing type
- * Double/Single hinged doors and electronic number lock
- * Adjustable inner shelves.
- * Inner safety cabinet with lock and key.
- * 220V Buzzer, light and door switch.
- * Completely made of GI powder coated sheets and tubes
- * On height adjustable foot studs with/without bottom skirting



WORK TABLES



WT- POP

Prosthetic Work Table



PT- 120 WA


Packing Table

Packing Table-Double Side.

- * 40mm thick reinforced top with water proof laminated boards, covered with stainless steel sheets.
- * One under shelf.
- * Two tire over shelf with roof shelf and pouch reel holder.
- * One Sliding drawer
- * Provision for holding Wire baskets in one side and stainless steel tray at another side
- * Completely made of stainless steel 304 grade
- * # 4 finish sheets and tubes.
- * On PVC bushes.

Options

- *Stainless steel wire baskets
- * Solid tray
- *Additional drawers
- *Wheels Heavy Duty



WT- US- 150

Work Table with Undershelf

Work Table

- * With/without under shelf and back splash.
- * 40mm thick reinforced top with water proof laminated boards, covered with stainless steel sheets.
- * With/Without drawers
- * Completely made of stainless steel 304 grade # finish sheets and tubes.
- * On Height adjustable foot studs/Wheels



WT- US-100

Work Table with Undershelf



WT- SD- 3

Work Table with 3 Drawer



Inspection Table with Lights.

- * Completely made of stainless steel 304 grade sheets and tubes.
- * Milky white acrylic sheet on the top of the table.
- * Fluorescent light and control switch.
- * Heavy duty swivel castors, two with brakes.
- * Welded fabrication.

INT - 150W

Inspection Table with Light

WORK TABLES

SINK TABLES



Triple bowl sink cabinet

- * With Hinged doors
- * Completely made of 304 grade #4 finish stainless steel sheets.
- * Triple bowl sink
- * With Back Splash
- * On height adjustable foot studs



WC-SGD

Wall Cabinet with Sliding Glass Doors



BC- SSD

Sink Cabinet with Sliding Door

BC- DN - W

Base cabinet with drawer
nest & work space

BC- SD

Base Cabinet with Sliding Door



BC- 2B- SD

Double Bowl Sink Cabinet with Sliding Doors



BC-1B-HD

Sink Cabinet with Hinged Doors



WT- 2B-U

Double Bowl Sink Table



WT- 1B-U

Single Bowl Sink with Undershelf



HWS

Hand Wash Sink- Knee Operated

Sink Cabinet/Table

- * With Hinged doors/under shelves
- * Completely made of 304 grade #4 finish stainless steel sheets.
- * With Sink
- * With Back Splash
- * On height adjustable foot studs /Wheels
- * Without Mixer tap

Options

- * Mixer/gooseneck type water faucet
- * High pressure mixer rinse
- * Extra bowls
- * Drawers

Hand Wash Sink

- * Completely made of stainless steel sheets
- * 300mm high Back Splash
- * Knee operated mixer tap
- * Size: Dia 30X15 cm



DSS- CWS

Double Station Scrub Sink with Digitally Controlled Water & Soap Dispensing System.



DSSU

Double station scrub sink-
Wall Mounted



SSSU

Single Station Scrub Sink

Scrub Sink-Single/Double

- * Free standing/wall mounted unit
- * Perforated drain board
- * High back splash
- * Sensor type soap dispenser
- * Completely made of SS 304 grade
- * # 4 finish sheets and tubes



DSSU

Double Station Scrub Sink

Options

- * Sensor Mixer water faucet
- * Sensor water faucet
- * Sensor type soap dispenser
- * Knee operated water flow system
- * Elbow operated water faucet

SCRUB SINK

DIRTY UTILITY



SH- WS
Sluice Sink



CAS- S
Cast Sink



SH-NS-120
Slop Hopper



JAS
Janitor Sink

- * Completely made of stainless steel 304 grade # 4 finish sheets & rods
- * Single bowl sink with perforated bin
- * Wall mounting/Free standing type
- * Back Splash

Options

- * Sensor Mixer water faucet
- * Health Shower



ESC- 1DR

Endoscopic Cabinet with Rotating Holder



ESC- 2GD

Endoscopic Cabinet with Glass Door



ESC- 2HD

Endoscopic Cabinet with Solid Door

Endoscopic Cabinet

- * Completely made of stainless steel 304 grade
- * # 4 finish sheets and tubes.
- * Rotation/straight holders
- * Removable tray at bottom
- * Ventilation Louvers
- * Solid/Glass doors with lock and keys.
- * On Height adjustable foot studs/wheels.
- * Welded fabrication.

ENDOSCOPY

TROLLEYS



MP- 70

Multipurpose Trolley

- * Made of stainless steel 304 grade rods and pipes.
- * Two tier wire baskets.
- * Pull handle.
- * Front side 4" swivel castors, with brakes.
- * Back side 10" rubber wheels.
- * Welded fabrication.



POP- 70

POP Trolley

- * Completely made of 304 grade #4 finish stainless steel sheets, tubes & rods.
- * Sliding drawers (different size)
- * Top with three sides raised edges.
- * One removable tray at the top.
- * Pull handle.
- * Heavy duty 4" swivel casters, two with brakes.
- * Welded fabrication



WAD- 50

Ward Trolley

- * Made with 304 grade #4 finished stainless steel rods and pipes.
- * Two shelves with three side lipping
- * Three side rails
- * Pull handle.
- * On heavy duty swivel castors, two with brakes..
- * Welded fabrication.



CT- 100
Curved Trolley



ITT- 50
Instrument Transport Trolley



WBT- 4
Wirebasket Trolley



PFT- 90
Platform Trolley



TROLLEYS



MTH- 3
Hydraulic Mayo Trolley



MT- 2S
Mayo Trolley 2 wheel



MT-4S
Mayo Trolley 4 wheel

Mayo Trolley

- * Completely made of stainless steel 304 grade
- * # 4 finish sheets and tubes
- * Stainless steel tray
- * Secure locking system for height adjustment
- * Height adjustment from 80 to 130cm
- * 50 mm-On wheels



SBT- 50
Swab Trolley

Swab Trolley

- * Completely made of stainless steel 304 grade
- * # 4 finish sheets and pipes.
- * Five stainless steel swab rails, each rails to hold 10 swabs.
- * Five number of stainless steel swab holder with 10 spikes
- * Removable stainless steel drip tray on bottom.
- * Welded fabrication.

Options:

- * Two/Four wheel type,
- * Removable tray,
- * Hydraulic height adjustable system



DT- 501-31
Dressing Trolley with 3 Drawers



DT- 501-11
Dressing Trolley with 1 Drawer



DT-501-21
Deressing Trolley 2 Drawers



DT-501-10
Dressing Trolley

Dressing Trolley

- * Completely made of 304 Grade #4 fish stainless steel sheets & tubes
- * Two Shelves with three side raised edges
- * Four nos of heavy duty swivel castors, two with brakes
- * Three side rails on the top and bottom shelf
- * Welded fabrication

Options

- * Pull handle
- * Impact absorbing bumpers
- * Drawers

TROLLEYS



IT- 503-12
Instrument Trolley



IT- 502-21
Dressing Trolley with Two drawers



IT- 502- 10
Dressing Trolley



IT- 503-10
Instrument Trolley

Instrument Trolley

- * Completely made of 304 Grade #4 fish stainless steel sheets & tubes
- * Two shelves with three side raised edges
- * Four nos of heavy duty swivel castors, two with brakes
- * Three side bended downwards
- * Curves corners
- * Welded fabrication

Options

- * Pull handle
- * Impact absorbing bumpers
- * Trash bin

IT- 505.11

Theater Trolley



IT-504.12

Instrument Trolley



- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Two shelves with three sides bend upward and one side bend downward,removable type.
- * Two sets of side pipes for reinforcement
- * Four nos of heavy duty swivel castor, two with brakes.
- * Knock down-Welded.



IT- 507.10

Service Trolley

IT- 507.11

Service Trolley

Service Trolley

- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Two shelves with four sides bend downward.
- * Four nos of heavy duty swivel castors, two with brakes.
- * Welded fabrication.

IT- 505.10

Theater Trolley

Theater Trolley

- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Two shelves with four sides raised edges
- * One side pull handle.
- * Four nos of heavy duty swivel castors, two with brakes.
- * Welded fabrication.

TROLLEYS



IT- 508-11R
Instrument Trolley with Rails



IT- 508. 11
Instrument Trolley



IT- 508.11D
Instrument Trolley with Drawer



IT- 508-12R
Instrument Trolley with 3 Shelf

Instrument Trolley

- * Made of stainless steel 304 grade #4 finish sheets and tubes
- * Two/Three flat Shelves with four sides bend downward.
- * Two side pull handle
- * Four nos of light duty swivel castors, two with brakes.
- * Welded fabrication

Options

- * Side rails
- * Impact absorbing bumpers



IT-509.10
Instrument Trolley



IT-509.11
Instrument Trolley

Instrument Trolley

- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Two/ flat shelves with four sides bend downward.
- * Four nos of light duty swivel castors, two with brakes.
- * Welded fabrication

Options

- * Side rails
- * Impact absorbing bumpers



MT-4W
Monitor Trolley

TROLLEYS



IT- 510.10
Instrument Trolley



IT- 510.21
Instrument Trolley 2 Drawer



IT- 510-5D
Instrument Trolley 5 Drawer



IT- 510.11
Instrument Trolley 1 Drawer

Instrument Trolley

- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Two/flat shelves with four sides bend downward.
- * Rails on each shelf
- * Four nos of light duty swivel castors, two with brakes.
- * Welded fabrication

Options

- * Side rails
- * Drawers
- * Impact absorbing bumpers



CC-2S-BB
Cast Cutting Trolley



WB-SB-6
Single Wash Basin Trolley



NT-506-2D
Nursing Trolley

- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Two/flat shelves with four sides bend downward.
- * Rails on each shelf
- * Four nos of light duty swivel castors, two with brakes.
- * Welded fabrication

Options

- * Side rails
- * Drawers
- * Welded fabrication
- * Impact absorbing bumpers



WB-DB-6
Double Wash Basin Trolley

Basin Trolley

- * Completely made of stainless steel 304 grade #4 finish sheets and tubes
- * Removable stainless steel basin 6 liter capacity,
- * Mobile in twin wheels with brakes
- * Welded fabrication

TROLLEYS

PT- 3D
Procedure Trolley 3 Drawer



CCC- B6
Critical Care Cart

Critical Care Cart

- * Completely made of stainless steel 304 grade #4 finish sheets and tubes
- * Top plain work surface.
- * One sliding drawer.
- * Six nos of removable wire baskets.
- * Heavy duty 4" swivel castors, two with brakes.
- * Welded fabrication



PT- 5D
Procedure Trolley 5 Drawer

Procedure Trolley

- * Completely made of 304 Grade #4 finish stainless steel sheets & pipes.
- * Five nos. of sliding drawers (different size)
- * Three/Five numbers of sliding drawer
- * Top and bottom drawer without divider and middle three drawers with dividers
- * Top with three sides raised edges
- * Pull handle
- * Heavy duty 4" swivel casters, two with brakes
- * Welded fabrication



CTT 24
Cylinder Transport Trolley



CTT- 12
Cylinder Transport Trolley



CTT- 24D
Cylinder Transport Trolley- Double

- Cylinder Transport Trolley.**
- * Frame made of Aluminium/ stainless steel tubes & plates
 - * 12/24/48, cylinders
 - * Height adjustable pull handle
 - * Two fixed castors
 - * Nylon tape to hold the cylinder
 - * Welded fabrication



CTT- D48
Cylinder Transport Trolley

EMC. 1
Emergency Cart



EC- 4D-3T
Emergency Cart with 4 Drawers & 3 Tilt Bin

Carts

- * Made with 304 grade #4 finished stainless steel sheets, rods and pipes.
- * 5 nos of sliding drawers, different sizes with powder coated
- * Three side rails on top.
- * Two side pull handles.
- * Heavy duty 4" swivel casters, two with brakes.
- * Sliding writing board.
- * Four side protection bumpers.
- * Centralized locking system.



ANC.1
Anesthesia Cart

Options

- * Corner defibrillator shelf.
- * Oxygen cylinder holder.
- * Height adjustable IV pole.
- * Waste bin
- * Tilt bin
- * Sharp container holder
- * Gloves box dispenser.
- * Extension cord with 4 plug points at the back.



MC-15-2D
Medication Cart



IC- 4DP
Isolation Cart

Medication Cart

- * Made with non magnetic stainless steel sheets.
- * Twenty five nos of crates on both sides (Total 50)
- * One writing board
- * Three side rails on top.
- * Two side pull handles.
- * Two sliding drawer
- * Heavy duty swivel castors, two with brakes.
- * With corner bumpers.
- * Central locking system

Option:

- * Digital Lock

CART



CSD- 120

CSSD Sterile Transport Trolley



CSD- 90

CSSD Sterile Transport Trolley



CSD- 90P

Powder Coated Transport Trolley

CSSD Transport Trolley

- * Completely made of stainless steel 304 grade # 4 finish sheets and tubes.
- * Two nos of inner solid shelves.
- * One side pull handle.
- * Two hinged door with lock and keys.
- * Heavy duty swivel castors, two with brakes.
- * Welded fabrication.

Options:

- * Powder coated paint finish
- * Wire shelves

Spica Table.

- * With acrylic support on head portion and removable Stainless steel support and teflon rod for leg portion.
- * Head portion is height adjustable and base portion is length adjustable.
- * Completely made of stainless steel 304 grade sheets, flats, rods and pipes.



SPCT

Spica Cast Table



IBDS

Infant Blood Drawing Station

Infant Flood Drawing Station

- * Welded stainless steel frame
- * Top has 2 rims with comfortable padding and straps sitting above the HPL Sheet
- * Removable middle padded section
- * Stainless steel tray on side
- * Height adjustable stand from 90 cm to 110 cm
- * Free standing type
- * Welded fabrication



WB

Wax Bath

- * Double walled body, well insulated.
- * Wax chamber made of 304 grade #4 finish Stainless steel sheet.
- * Wax Chamber: Size: 50cm x 35cm x 25cm.
- * Heaters: Two numbers of 1kw each, special wax Heaters.
- * Heat Control: Thermostat to control temperature from 00 to 85°C
- * Indicator Lamps: Two numbers for Mains & Thermostat."
- * Mounted on heavy duty swivel castors
- * Power: 220-240 V AC. (Load, 2kw).
- * Finish: Externally Powder Coated finish.

CART



FTT- 50

File Transport Trolley

Stainless Steel File Transporting Trolley

- * With two side pull handles.
- * Well type file storage tray 35 cm depth with file holding rods inside and two hinged lids.
- * With one bottom shelf.
- * On heavy duty wheels, two with brakes.
- * Completely made of stainless steel 304 grade # 4 finish sheets and tubes.



SUT- TB

Suturing Trolley

Suturing Trolley

- * Completely made of stainless steel 304 grade sheets and pipes.
- * With two compartment
- * One side with Four sliding drawers
- * Other side one removable wash basin and suture waste bin
- * Top Shelf with four sides inclined upward.
- * Four nos of swivel castors, two with brakes.
- * Welded fabrication.



CLC- 3S

Clean Linen Cart with Sliding Door



CLC- 3H

Clean Linen Cart with Hinged Door

Clean Linen Cart.

- * Completely made of stainless steel 304 grade # 4 finish sheets and tubes.
- * Hinged/sliding doors with lock and keys.
- * Two nos of inner shelves.
- * Pull handle.
- * Four nos of heavy duty swivel castor, two with brakes.
- * Welded fabrication.



PD- 120T

Paper Dispensing Cart on wheel

Paper Dispensing Cart

- * Completely made of stainless steel 304 grade tubes.
- * Four nos of tubes with catchers to hold the papers.
- * One side pull handle.
- * One wire shelf at the bottom.
- * Four nos of heavy duty castors, two with brakes.
- * Welded fabrication.



LT- SB

Loading Trolley with Solid Bottom



LT- WB

Loading Trolley with Wire Bottom



MVC- SS

Mobile Video Cart

Loading trolley.

- * Completely made of stainless steel 304 grade sheets, rods and tubes.
- * Foldable wire shelf/solid shelf
- * On heavy duty swivel castors two with brakes.
- * Solid/wire base.
- * Side protection bumpers.

Mobile Video Cart

- * Completely made of stainless steel 304 grade # 4 finish sheets and tubes.
- * Five shelves and one sliding drawer at the bottom.
- * One control switch in the front.
- * Four nos of three pin plug points with individual switch.
- * Heavy duty swivel castors, two with brakes.
- * Welded fabrication.



DBS- 60

Stainless Steel Dustbin 60 ltr



PC- DB- 60

Powder Coated Dustbin 60 ltr



DBS- 20

Stainless Steel Dustbin 20 ltr



DBRD- 10

Stainless Steel Dustbin Round



PC-DB- 20

Powder Coated Dustbin 20 ltr

Stainless Steel Dustbin

- Completely made of non magnetic stainless steel sheets
- 20,30,40 and 60 litres capacity
- Pedal operated Slow closing lid
- Light weight removable tank

Options

- Round dust bins
- With wheels

Powder Coated Dustbin

- Completely made of powder coated CR sheets
- 20,30,40 and 60 litres capacity
- Pedal operated Slow closing lid
- Light weight removable tank

Options

- Powder coated dust bins
- With wheels
- Different colours

DUST BIN

BBC
Baby Bassinet Cabinet



Baby Bassinet

- * Electro galvanized oven baked paint finish tubes and sheets
- * Both side Dia 25 mm pipe with holder
- * bottom storage cabinet with hinged doors & magnetic catcher
- * Completely made of G.I Powder coated finish sheets and tubes.
- * Four nos of heavy duty swivel castors, two with brakes.
- * Welded fabrication
- * With acrylic Baby Basket

BBT
Baby Bassinet Trolley



OV- BT
Over Bed Table



SC- T
Sharp Container Trolley



DL- SBD2L

Dirty Linen Trolley Double bag with Lid



LS- WL

Linen Skip with Lid

Dirty Linen Trolley - Square with Lid
 * Structure made of stainless steel pipes.
 * Pedal operated slow closing lid.
 * Double/Single removable nylon bag.
 * Mobile on swivel wheels with brakes.
 * Welded fabrication.



DL- SB1L

Dirty Linen Trolley Single Bag Square with Lid

LAUNDRY



DL- SB
Dirty Linen Trolley Single Bag square



DL-SB2L
Dirty Linen Trolley Double Bag square



DL- RS
Dirty Linen Trolley Single Bag Round



DL- RD
Dirty Linen Trolley Double Bag Round

Dirty Linen Trolley - Round
* Structure made of stainless steel pipes.
* Double/Single removable nylon bag.
* Mobile on swivel wheels with brakes.
* Welded fabrication.



CT- W100

Wire Mesh Trolley

Stainless Steel Wire Mesh Trolley

- * Completely made of stainless steel 304 grade rods and tubes.
- * Frame made of 25mm square tubes.
- * All sides with 6mm rod mesh.
- * Top open.
- * On heavy duty swivel castors, two with brakes.



CT- 100

Laundry Trolley



BSLT- 120

Box Sacks Linen Cart

- Bulk linen transport cart
- Completely made of stainless steel 304 grade #4 finish sheets and tubes
- Top side half door
- Front side half folding door
- Welded fabrication
- Four side corner bumpers
- Heavy duty swivel castors

LAUNDRY



GWSU-SS-150

Pathology Grossing Work Station

Pathology Grossing Work Station -PGW-150

Size: 1500x750x1900 mm

Pathology Grossing Work Station –PGW-150

- * Completely made of stainless steel 304 grade #4 finish sheets and tubes.
- * Backdraft ventilation system along the entire work station
- * Fume exhaust pipe with Hepa filter and removable perforated stainless steel cover
- * Tissue holder
- * UV Disinfecting system
- * Polythene dissection board with ruler
- * Single bowl sink unit with fittings and removable perforated working top
- * Goose neck type Mixer water faucet.
- * Sprinkler system for sink type stain station
- * Fixed glass (tempered) covering - half section
- * Spot lights & fluorescent tube lights
- * Bottom left side with removable fresh and used formalin containers holder on wheels
- * Bottom Right side with Storage compartment
- * Free standing unit with height adjustable foot studs

Optional Accessories:

- * Electrical sockets and connections.
- * Document holder
- * Magnetic Tool bar holder
- * Formalin dispensing system
- * Sliding glass (tempered glass) covering for work station
- * Magnifying glass
- * LED Monitor Mount
- * Camera stand



AP-OWS-100

Lead Apron Hanger Two Tier on wheel



AP-OW-100

Lead Apron Hanger on Wheel



AP-WM-110

Lead Apron Hanger wall Mounted

Lead Apron Hanger

- * Completely made of stainless steel 304 grade rods and sheets.
- * Free Standing/Wall mounted
- * Space to hang 5-10 nos of lead aprons.
- * Welded fabrication.

LEAD APRON HANGER

PHARMACY SHELVING



ST- WM.D
Pharmacy Shelving



2B- WT.LBD
Pharmacy Shelving with Lin Bins



PC-WT
Pharmacy cabinet with tray



2B-180
Pharmacy Shelving with Sliding Drawer



HD-RS

Heavy Duty Rack Shelving



PR-100

Pot Rack

Stainless steel Pot Rack. Single tier.
Completely made of stainless steel 304 grade pipes.



SW-HK-120

Single Tier Wall Shelf with Hook

Stainless steel wall shelf Single tier.
* Completely made of stainless steel 304 grade
4 finish sheets
* With wall mounting brackets
* With Hooks



SW-120

Wall Shelf

Stainless steel wall shelf Single tier.
* Completely made of stainless steel 304 grade # 4 finish sheets.
* With wall mounting brackets

STORAGE SHELVING

STORAGE SHELVING



WS-120-W

Wire Shelving on wheel



WS-120

Wire Shelving



WS-120H

Wire Shelving



WST-5

Wire Shelving

Stainless Steel Three tier Height adjustable Wire Shelf Unit on Wheels

- * Shelf made of mat wires Ø3mm rod and shelf frame made of Ø6mm rods act as natural frame to help prevent items from falling down
- * Shelf clip - two piece plastic split sleeve
- * Shelf collar Tapered collar snugly locks on tapered plastic split sleeves.
- * Completely made of Stainless steel rods and tubes.
- * On heavy duty Swivel castors/height adjustable foot studs



LBW Lin Bin Holder Panel- Wall Mounted

- * Lin bin holder. Completely made of GI powder coated sheets.
- * Wall mounted /on heavy duty wheels two with brakes.



SSU - 100
Solid shelving unit



LBH-WL
Lin Bin Holder on Wheel

Solid shelving unit

- Completely made of 304 grade #4 finish stainless steel sheets and tubes.
- Fixed shelves.
- Height adjustable foot studs.
- Welded fabrication.

Options

- On heavy duty swivel castors
- Perforated shelves

STORAGE SHELVING

CCSS

Compact Cabinet Storage System

- * Orkitch Brand Powder Coated Compact Storage Shelving
- * System for Pharmacy- With 11 Movable Units and one fixed unit.



BFSS

Boltfree Shelving System

- * Height Adjustable Bolt free storage shelving unit.
- * Posts with slotted angle.
- * Completely made of GI Powder coated sheets.
- * Optional Colours



EC-BC-601

Examination Couch with Cabinet



EC-601

Examination Couch



Examination Couch

- * Structure made stainless steel 304 grade # 4 finish tubes.
- * Black imitation leather upholstery with 50 mm thick foam padding.
- * Single hand operated adjustable back rest.
- * Paper roll holder
- * Welded fabrication.



EC-604

Massage Couch

Massage Couch

- * Heavy duty structure made from stainless steel 304 grade tube.
- * Black imitation leather upholstery with 50mm thick foam padding.
- * Single hand operated adjustable back rest, with nose cut.
- * Paper roll holder.
- * Welded fabrication.



EC-602

Examination Couch with Side Rails

Examination Couch with Side Rails

- * Structure made of stainless steel 304 grade tubes
- * Black imitation leather upholstery with 50mm thick foam padding.
- * Single hand operated adjustable back rest.
- * Paper roll holder
- * Welded fabrication.
- * PVC foldable side rails



PB-2

Ward Bed

- * Structure made of mild steel perforated sheet and pipes.
- * Epoxy powder coated.
- * Four legs with rubber bush
- * Knock down fabrication.
- * Without mattress



PB-1C

Patient Bed with Single Crank



PB-2C

Patient Bed with 2 Crank



GEC-E
Gynecology Examination
Couch Electrical

GYNAECOLOGY AND OBSTETRIC EXAMINATION CHAIR - 3 MOTORS

- * Designed with gynaecological backrest and trendelenburg
- * Height adjustment with electronic motor and hand control unit
- * Battery
- * Long lived, easily-cleaned, water resistant and antibacterial special furnishing with polyurethane foam
- * Formed by 2 main parts (back & seat).
- * Easily removed legs and support
- * 360° movable castors with linear brake



EC- 605

Gynecology Examination Couch

Gynecology Couch

- * Heavy duty structure made of stainless steel tube.
- * Three sections, black imitation leather upholstery with 50mm thick foam padding
- * Single hand operated adjustable head and foot sections.
- * Paper roll holder
- * Adjustable lithotomy poles with PU leg supports.
- * Removable stainless steel tray



EEC- 610

Electrical 2 Section Couch



ECH-61H

Hydraulic Examination Couch -3 section

Hydraulic Pedal operated height adjustable Examination Couch.

- * Structure made of stainless steel /powder coated tubes.
- * Blue imitation leather upholstery with foam padding.
- * Single hand operated adjustable back rest.
- * Welded fabrication.



ECH-613

**Hydraulic Examination Couch with
Side Rails**



**Hydraulic Examination Couch with Padded
Side Rails**

ECH-614

Hydraulic Examination Couch.

- * Structure made of stainless steel/powder coated tubes.
- * Blue imitation leather upholstery with foam padding.
- * Single hand operated adjustable back rest.
- * Welded fabrication.
- * PVC /Padded side rails.



EEC-615

Electrical Examination Couch- 4 section



EEC-614

Electrical Examination Couch- 3 section



BA- CT
Baby Coat

PEDIATRIC BED(BABY COT)

- * Sliding side rails on both sides.
- * High density foam mattress with fluid resistance cover
- * On Swivel castors, two with brakes.
- * Welded construction
- * Completely made of GI powder coated sheets and tubes



BSL- 50
Bed Side Cabinet

PEDIATRIC BED

Infusion Stand - IV - PB4

- * Telescopic stainless steel pipes.
- * Secure locking system for height adjustment.
- * Height adjustable from 110-195 cm
- * Four nos of stainless steel safety hooks.
- * Five legs nylon base.
- * Five twin wheel swivel castors
- * Knock down fabrication.



IV- PB4

Infusion Stand with Nylon Base



IV-SS4

Infusion Stand with Stainless Steel Base



IV- SPB-4

Heavy Duty Spider Leg Base Infusion Stand



IV- HD

Heavy Duty Infusion Stand

Infusion Stand

- * Telescopic stainless steel pipes.
- * Secure locking system for height adjustment.
- * Height adjustable from 110-195 cm
- * Four nos of stainless steel safety hooks.
- * Five legs nylon base.
- * Five twin wheel swivel castors
- * Knock down fabrication.



IV- SP 4

Syringe Pump Stand

Syringe Pump Stand

- * Telescopic stainless steel
- * 304 Grade pipes & plates.
- * Secure locking system for height adjustment
- * Height adjustable from 120-210cm
- * Four nos of stainless steel safety hooks.
- * Four socket, 3 pin plug point
- * Clinch to clutch the wire
- * Pull handle.
- * Stainless steel heavy duty spider leg base.
- * Five heavy duty swivel castors
- * Welded fabrication



IV-MST

Mobile Suction Trolley



MS-WB

Monitor Stand with Wire Basket
MRI Compactable



IV- MRI

Infusion Stand -MRI Compact

INFUSION STAND



PPEH-P

PPE Holder Powder Coated



GD-PC

Powder coated Gloves Dispenser



TD-P

Tissue Dispenser



GD-SS

Stainless Steel Gloves Dispenser



PPE-AL

Acrylic PPE Holder

PPE Holder-Stainless Steel

* Cabinet with Three gloves dispenser, Two face mask and one gown Holders

* Completely Made of 304 grade #4 finished stainless steel sheets.

Options:

Acrylic
 Powder coated
 Custom design



ML-OW

Magnifying Lamp on Wheel



EL- WM

Examination Light Wall Mounted

LED Examination Light Adjustable Height

- * stainless steel central tubes.
- * Mobile on twin wheel swivel castors ϕ 50mm.
- * Flexible goose neck arm 700mm.
- * LED Bulb
- * Input voltage: 220/230V 50-60Hz
- * Knock down fabrication.



EL- LED

Examination Light LED

Option:

- *Free standing
- *Wall mounting

EXAMINATION LIGHT



FS-SR
Single Step Footstool



FS- DR
Double Step Footstool



FS- SS
Single Step Foot Stool- Stainless Steel



FS- DS
Double Step Foot Stool- Stainless Steel

Single/Double Step Foot Stool

- * Frame made of stainless steel pipes.
- * Step made of MDF covered with antislip rubber pad/SS chequered plate
- * Rubber foot caps for stability.
- * Welded fabrication.



FS- 2R
Double Step Foot Stool with Two Side Rail



RS- SS
Revolving Stool SS Seat



FS- TR
Triple Step Foot Stool



FS- SR1
Single Step Foot Stool with One Side Hand Rail

- Single/Double Step Foot Stool**
- * Frame made of stainless steel pipes.
 - * Step made of MDF covered with antislip rubber pad/SS checkered plate
 - * Rubber foot caps for stability.
 - * With/without side rails
 - * Welded fabrication.



RS- LB

Revolving Stool - Low Backrest



RS- PBF

Revolving Stool



RS- HB

Revolving Stool with High Backrest



RS-SDR

Revolving Stool with Side Rest

Revolving Stool

- * Without back rest /Low or Hight Back Rest
- * Gas lift for height adjustment from 520-710mm.
- * Seat with low back rest made from black imitation leather upholstery with 50mm thick high density foam padding
- * Chrome plated foot support ring.
- * Five leg nylon/ Chrome plated base
- * Mobile on five twin wheel swivel castors
- * Knock down fabrication.



BD-M

Blood Donor Chair



Blood Donor Chair

- * Manually adjustable back rest, can be adjusted up to 180°
- * Manually adjustable leg rest and convert to a bed.
- * Padded armrest, adjustment in any position for patient conveniences.
- * Completely made of stainless steel 304 grade tubes.
- * Black imitation leather upholstery with high density foam padding

BLOOD DONOR CHAIR



XV-HA 1

X-ray film viewer Double Standard



XV- HA 2

X-ray film viewer Single Standard



XV- TB 1

X-ray film viewer Single Slim



XV- TB 2

X-ray film viewer Double Slim

X ray Film Viewer Single slim.

- * Case made of steel powder coated finish.
- * High intensity screen
- * Variable adjustment of screen luminance.
- * CCFL is used for luminance.
- * Wall mounting model.
- * Screen white polymethyl methacrylate plate.



UBH
Urine Bag Holder



WS-3 F
3 Fold Ward Screen



BD-80
Brush Dispenser



WS-4F
4 Fold Ward Screen

Three/Four Folds Ward Screen

- * Frame completely made of stainless steel tubes
- * Folding panels, each 610mm wide
- * PVC washable curtains for three/four sections
- * Easy movement on twin wheel nylon swivel castor

WARD



KB- 6
Kick Bowl



BP- HW
Bed pan/Urinal Holder



KBK- 20
Kick Bucket

Kick Bucket

- * Five leg stainless steel tube base
- * Five vertical supports with rubber sleeves to hold the bucket
- * Mobile on 5 twin wheel swivel castors
- * Base diameter is 400mm
- * Height is 360mm
- * Bucket with lid



SR- FS
Shoe Rack Standard



SR- OW
Shoe Rack on Wheel



WBKC
Wire Basket



WBK
Wire Basket



WB-WMB
Wire Basket with Wall Mounting Bracket



PSB
Perforated Sheet Basket

WARD



SD- EL

Soap Dispenser Elbow Operated



SCD- N

Soap Dispenser- Normal



CHP- A4

Patient Chart Holder -Plastic



SD-TM

Soap Dispenser Table Mounted

- * Made of non magnetic stainless steel sheets & rods
- * Wall mounting type, Elbow operated pressing system, Table mounting type
- * Provision to hang plastic soap dispenser
- * Welded fabrication



CHS-A4

Patient Chart Holder- Stainless Steel

- * Board made from 1mm thick non-magnetic stainless steel sheets, smooth work surface whilst being light & easy to clean
 - * Stainless steel hook to hang on the bed
 - * Nickel plated spring clip
- Code: CHS-A4 Board Size: Length: 230mm Width: 340mm
 Code: CHS-A3 Board Size: Length: 460mm Width: 340mm




CU-A

Adult Crutches- Arm pit




CU-C

Child Crutches-Arm pit



QWS-S

Quad Walking Stick Straight Handle



QWS-C

Quad Walking Stick Curved Handle



WFC-C

Walking Forearm Crutches-
Closed Type



WFC-O

Walking Crutches-
Open Type



WS-F

Walking Stick Foldable

CRUTCHES & WALKING SUPPORTS



WFF-W
Walking Frame Foldable with Wheel



WFF-S
Walking Frame Foldable Standard



WF-S-2W
Walking Frame with Seat on Two Wheel



WF-S-4W
Walking Frame with Seat on Four Wheel



WC- P-NF

Wheel Chair with
Padded Seat Non Foldable.



WC- F

Wheel Chair Foldable



BC-HD

Heavy Duty Bariatric Chair


WHEEL CHAIR

CATHETER BASKET



CB- L

Catheter Basket -Long



CB- W

Catheter Basket -Wiremesh



CB- DV

Catheter Basket- Divided

- Completely made non magnetic stainless steel sheets and rods
- Mounting bracket
- Frame made of 8mm rods
- Basket made of 1.5 mm sheets with perforation
- Welded fabrication

CW.T

Cadaver washing Trolley



Cadaver Washing Trolley

- * Completely made of stainless steel 304 grade sheets.
- * On heavy duty wheels two with breaks
- * With stainless steel bucket
- * Drain system



CCT

Cadaver closed Transport Trolley

Cadaver Trolley

- * Completely made of stainless steel 304 grade sheets.
- * On heavy duty wheels two with breaks
- * With curved stainless steel cover



CWT-S

Autopsy Table

CADAVER TROLLEY



PA-CH-S

Pass through Chamber Sliding Door

Stainless steel Pass through Hatch

- * Mechanical type.
- * Chamber consists of mechanically operated hatches, stainless steel framed clear acrylic/tempered glass sliding doors upward.
- * Only one door can be opened at a time.
- * Completely made of 304 grade # 4 finish sheets and tubes.
- * Wall mounting type.



PA-CH-H

Pass through Chamber Hinged Door

Pass through: hinged door

- * Custom Fabricated stainless steel Pass through Hatch - mechanical type.
- * Chamber consists of mechanically operated hatches, stainless steel framed glass hinged doors.
- * Completely made of 304 grade # 4 finish sheets & rods.
- * Wall mounting type.



PAW-H

Pass through Window Hinged Door - Touchless Sensor



CLADDING

LAB





Electrically controlled dental chair-EC-DNTL

- * LED Operation Light
- * Three way syringe-2 nos
- * X-ray film viewer
- * Central suction system
- * Automatic water heating system
- * Turntable assistant tray
- * Voltage: 220V/50-60Hz



DNTL-BC

Dental Cabinet:

- * Completely made of GI powder coated sheets and tubes.
- * With Single bowl sink
- * Top: with white Korean top
- * Two mobile drawer nest.
- * Two side storage cabinet with hinged doors
- * Four no's of heavy duty swivel castors
- * Welded fabrication.

Option:

- * Water faucet
- * Color
- * Additional Drawers and shelves
- * Customized design

P No.	Items	Code	Measurement (CM)	Discription
14	Upright Cabinet 4 Hinged Doors	UC-4HD	120X45X180	* Four nos of hinged solid doors with eight shelves. * Custom made design available
	Upright Cabinet 3 Hinged Doors	UC-3HD	120x45x180	* Three nos of hinged solid doors with Four shelves. * Custom made design available
	Upright Cabinet 2 Hinged Doors	UC-2HD	90X45X180	* Two nos of hinged solid doors with Four shelves. * Custom made design available
	Upright Cabinet 20 Pegion Holes	UC-PG20	120x45x180	* Two nos stainless steel hinged doors at bottom * 20 nos of pegion holes (front and back open) on top * Custom made design available
15	Upright Cabinet 3 Glass Doors	UC-3GHD	120x45x180	* Three nos of hinged glass doors with Four shelves. * Custom made design available
	Upright Cabinet 2 Glass Doors	UC-2GHD	90X45X180	* Two nos of hinged glass doors with Four shelves. * Custom made design available
	Upright Cabinet 2 Sliding Doors	UC-2SD	120X45X180	* Two nos of sliding solid doors with Four shelves. * Custom made design available
	Upright Cabinet with 3 Drawers	UC- 4G3D	120X45X180	* Four nos stainless steel framed glass hinged doors & three nos of sliding drawers at centre * Custom made design available
16	Upright Cabinet 2 Glass Doors	UC-2GHD	90x45x180	* Two nos stainless steel framed glass hinged doors & four nos of height adjustable shelves. * Custom made design available
	Upright Cabinet Hinged Glass Doors with 4 Drawers	UC-2GHD-4D	90x45x180	* Two nos stainless steel framed glass hinged doors & four nos of sliding drawers at bottom * Custom made design available
	Upright Cabinet	UC-2HD	90x45x180	* Two nos stainless steel hinged doors with four nos of shelves. * Custom made design available
	Upright Cabinet 2 Hinged Glass door with Rails for Basket	UC-2GHD-R	90x45x180	* Two nos stainless steel hinged doors with rails for GI/Plastic containers. * Custom made design available
	Upright Cabinet 2 Hinged Glass Doors with 4 Drawers	UC-2HD-4SD	90x45x180	* Two nos stainless steel framed glass hinged doors on top with four nos of side by side sliding drawers at bottom * Custom made design available
	Upright Cabinet 2 Hinged Glass Doors with 4 Drawers	UC-2GHD-4D	90x45x180	* Two nos stainless steel framed glass hinged doors on top with four nos of full length sliding drawers at bottom * Custom made design available
	Upright Cabinet 2 Glass door with on casters	UC-2GSD-W	100x45x180	* Two nos stainless steel framed glass hinged doors on top with our nos of height adjustable shelves. * On Heavy duty swel castors two with breaks * Custom made design available
17	Upright Catheter Cabinet with 6 Rails for Baskets	UCC-6R-P	120x45x200	* Six nos of sliding rails, each with Five hooks * Custom made design available
	Upright Cabinet 1 Hinged door with Powder coated	UC-1 GHD-P	50x45x180	* Single GI framed glass hinged doors with four nos shelves. * Custom made design available
18	Locker with Single Door	STL-1 DP	40x47.5x180	* Made with oven bake paint finish Mild Steel sheet. * Custom made design available
	Locker with Two Door	STL-2DP	40x47.5x180	* Made with oven bake paint finish Mild Steel sheet. * Custom made design available
	Locker with Three Door	STL-3DP	40x47.5x180	* Made with oven bake paint finish Mild Steel sheet. * Custom made design available
	Z Type Locker	STL- 2ZP	40x45x190	* GI Powder Coated Z'-shape locker. * Double door with lock and keys.
19	HPL Locker			* HPL Lockers * customised design.

P No.	Items	Code	Measurement (CM)	Discription
20	Base Cabinet with Sliding Doors	BC-SD	150x70x90	* 40mm thick reinforced top with water proof laminated boards, covered with stainless steel sheets. * With one mid shelf and hinged doors.
	Base Cabinet with 8 Drawers	BC-8D	100X70X85	* 40mm thick reinforced top with water proof laminated boards, covered with stainless steel sheets. * With drawer nest (8 drawers) and hinged doors.
	Base Cabinet with Hinged door	BC-HD	150x70x90	* 40mm thick reinforced top with water proof laminated boards, covered with stainless steel sheets. * With one mid shelf and hinged doors.
21	Wall Cabinet with Hinged Glass Door	WC- HGD	150X30X90	* With hinged glass doors and one adjustable shelf .
	Wall Cabinet with Sliding Door	WC- SD	150X30X60	*With sliding doors and one adjustable shelf.
	Wall Cabinet with Hinged Doors	WC- HD	150X30X60	*With hinged doors and one adjustable shelf.
	First Aid Box	FB	40x17x60	* Made Of Powder coated /Stainless Steel Sheets & Tubes * Hinged door with singlev shelf and First Aid Sticker
22	Narcotic Cabinet with Number Lock	NC- ELD	60X30X60	* Single hinged door with electronic number lock * Two inner shelves and Inner safety cabinet
	Narcotic Cabinet with Double Door	NC- 315D	80X30X60	* Single hinged door with two nos of locks and keys * Two inner shelves and Inner safety cabinet
23	Upright Narcotic Cabinet	NC- UELD	110X40X180	* Double hinged doors and electronic number lock. * Two nos adjustable inner shelves and Inner safety cabinet
	Upright Narcotic Cabinet with Single Door	NC- USD	60X40X160	* Single hinged doors with locks and keys * Two nos adjustable inner shelves and Inner safety cabinet
	Narcotic Cabinet with Single Door	NC- 315S	60X30X60	* Single hinged doors with locks and keys * Two nos adjustable inner shelves and Inner safety cabinet
24	Narcotic Cabinet with Full Inner Door	NC- SDD	60X30X60	* Powder Coated * Full inner hinged door and adjustable inner shelves and inner safety cabine
	Narcotic Cabinet with Full Inner Door- Powder Coated	NC- DDP	60X30X60	* Single hinged doors and electronic number lock * Adjustable inner shelves and inner safety cabinet
	Upright Narcotic Cabinet with Single Door (Powder Coated)	NC- UD- 2LP	60x45x90	* Single hinged doors and electronic number lock * Adjustable inner shelves and inner safety cabinet
	Upright Narcotic Cabinet with Single Door& Locks (SS)	NC-UD- 2LS	60x45x90	* Single hinged doors and adjustable inner shelves and inner safety cabinet
	Narcotic Cabinet with Inner &Outer Door (Powder Coated)	NC-2D-PL	60x45x90	* Inner and outer doors * Adjustable inner shelves and inner safety cabinet
	Narcotic Cabinet with Inner &Outer Door(SS)	NC-2D-SS	60X30X60	* Inner and outer doors * Adjustable inner shelves and inner safety cabinet
25	Prosthetic Work Table	WT- POP	120x65x90+50	* Table with heigh back splash * With tools hanging provision
	Packing Table	PT- 120 WA	150x70+70x90+70	* Double Side Stainless Steel Packing Table. * Custom made design available
27	Work Table with Undershelf	WT- US- 150	150x70x90+10	* With under shelf and back splash.and 40mm thick reinforced top *Custom made design available
	Work Table with Undershelf	WT- US-100	100x70x90+10	* With under shelf and back splash.and 40mm thick reinforced top * Custom made design available

P No.	Items	Code	Measurement (CM)	Discription
27	Work Table with 3 Drawer	WT- SD- 3	130x70x90	* With under shelf and back splash.and 40mm thick reinforced top and drawers * Custom made design available
	Inspection Table with Light	INT - 150W	150x70x90	* Milky white acrylic sheet on the top of the table. * Fluorescent light and control switch. * Heavy duty swivel castors, two with brakes.
28	Tripple Bowl Sink Cabinet	BC- 3 B	210x70x90+10	* Triple bowl sink cabinet-Stainless steel * With Hinged doors
29	Wall Cabinet with Sliding Glass Door	WC-SGD	150x30x60	* Wall cabinet-Stainless steel * With sliding glass doors
	Sink Cabinet with Sliding Door	BC- SSD	120x70x90+10	* Single Bowl sink cabinet-Stainless steel * With sliding solid doors
	Base Cabinet with Sliding Door	BC- SD	140x70x90	* Base cabinet-Stainless steel * With sliding solid doors
	Base Cabinet with Double Side Drawer Nest and Work Space.	BC- 2DN-WS		* Base cabinet-Stainless steel * With two side drawer nest and central work space
30	Double Bowl Sink Cabinet with Sliding Door	BC- 2B- SD	120x70x90+10	* Single Bowl sink cabinet-Stainless steel * With sliding solid doors
	Sink Cabinet with Hinged Door	BC- 1B- HD	120x70x90+10	* Single Bowl sink cabinet-Stainless steel * With hinged solid doors
	Double Bowl Sink Table	WT- 2B- U	150x70x90+10	* Triple Bowl sink table-Stainless steel * With/without under shelf and with back splash
	Single Bowl Sink with Undershelf	WT- 1B-U	150x70x90+10	* Triple Bowl sink table-Stainless steel * With under shelf and with back splash
	Hand Wash Sink- Knee Operated	HWS	40x40x20+30	* Hand Wash Sink-Stainless steel * 300mm high Back Splash and Knee operated mixer tap * Round Bowl Size: Dia 30X15 cm
31	Double Station Scrub Sink with Digitally Controlled Water & Soap Dispensing System.	DSS- CWS		* Digitally Controlled Water & Soap Dispensing System.
	Double Station Scrub Sink	DSSU	160x65x125	* Wall mounted unit * Perforated drain board and High back splash * Knee operated /sensor type water flow system
	Double Station Scrub Sink	DSSU	160x65x147.5	* Free standing/wall mounted unit * Perforated drain board and High back splash * Knee operated /sensor type water flow system
	Single Station Scrub Sink	SSSU	90x65x147.5	* Free standing/wall mounted unit * Perforated drain board and High back splash * Knee operated /sensor type water flow system
32	Sluice Sink	SH- WS	100x65x90	* One hand wash sink and one conical bowl grating on top. * With Lid and Continuous water flushing system
	Cast Sink	CAS- S	120x60x90+10	* Single bowl sink and SS bin * Wall mounting/Free standing type
	Janitor Sink	JAS	60x60x70	* Small sink bowl on top with water faucet * Bottom large sink bowl with perforated gratings
	Slop Hopper	SH-NS-120	120x65x90	* Heigh Back Splash and Continious flushing rim * Conical bowl with Hinged gratings

P No.	Items	Code	Measurement (CM)	Discription
33	Endoscopic Cabinet with Solid Door	ESC- 2HD	100x45x235	* Straight endoscope holders and removable tray at bottom * Ventilation Luovers
	Endoscopic Cabinet with Glass Door	ESC- 2GD	100x45x235	* Straight endoscope holders and removable tray at bottom * Ventilation Luovers
	Endoscopic Cabinet with Rotating Holder	ESC- 1DR	110x60x235	* Rotating endoscope holders and removable tray at bottom * Ventilation Luovers
34	POP Trolley	POP- 70	60x50x95+20	* One side IV stand. * Two sides removable buckets and Five nos. of sliding drawers * Heavy duty 4" swivel casters, two with brakes.
	Ward Trolley	WAD- 50	60x45x90	*Two shelves with three side lippings and Three side rails *On heavy duty 75 mm swivel castors, two with brakes..
	Multipurpose Trolley	MP- 70	70x50x90	* Two tier wire baskets. * Front side 4" swivel castors, and Back side 10" rubber wheels.
35	Curved Trolley	CT- 100	130x70x90	* Top shelf with lipping on back side * On heavy duty swivel castors, two with brakes.
	Instrument Transport Trolley	ITT- 50	50x50x70	* Six nos of wire basketsand Sliding rails. * Heavy duty castors, two with brakes.
	Wirebasket Trolley	WBT- 4	60x60x130	* With two nos of removable wire basket on both sides * On heavy duty swivel castors, two with brakes.
	Platform Trolley	PFT- 90	96x72x20+74	* One long pull handle, foldable type. * Heavy duty castors two with brakes.
36	Hydraulic Mayo Trolley	MTH- 3	56x45x80+130	* Stainless steel tray 66x40cm * Adjustable height by hydraulic system.
	Mayo Trolley 2 wheel	MT- 2S	56x45x80+130	* Stainless steel tray 56x43cm * Secure locking system for height adjustment
	Mayo Trolley 4 wheel	MT-4S	56x45x80+130	* Stainless steel tray 56x43cm * Secure locking system for height adjustment
	Swab Trolley	SBT- 50	90x40x140	* Five stainless steel swab rails, each rails to hold 10 swabs. * Removable stainless steel drip tray on bottom.
37	Dressing Trolley	DT-501-10	51x50x90	* Two shelves with three side raised edges * Three side rails on the top and bottom shelf.
	Dressing Trolley with 1 Drawer	DT- 501-11	51x50x90	* Two shelves with three side raised edges and three side rails on the top and bottom shelf. * One sliding Drawer
	Deressing Trolley 2 Drawers	DT-501-21	51x50x90	* Two shelves with three side raised edges and three side rails on the top and bottom shelf. * Two sliding Drawer
	Dressing Trolley with 3 Drawers	DT- 501-31	51x50x90	* Two shelves with three side raised edges and three side rails on the top and bottom shelf. * Three sliding Drawer
38	Dressing Trolley	IT- 502- 10	51x50x90	* Two shelves with three side raised edges and
	Dressing Trolley with Two drawers	IT- 502-21	51x50x90	* Two shelves with three side raised edges *Two sliding Drawer
	Instrument Trolley	IT- 503-10	60x50x90	* Two shelves with four sides bended downward. * Curved corners.

P No.	Items	Code	Measurement (CM)	Discription
38	Instrument Trolley	IT- 503-12	100x50x90	* Two shelves with four sides bended downward. * Curved corners.
39	Instrument Trolley	IT-504.12	60x46x90	* Two shelves with three sides bended upward and one side bended downward, * Removable type. * Knock down fabrication
	Theater Trolley	IT- 505.11	90x45x90	* Two shelves with four sides bended down ward
	Theater Trolley	IT- 507.10	73x52x90	* Two shelves with four sides bended down ward
	Service Trolley	IT- 507.11	73x52x90	* Two shelves with four sides bended upward. * One side pull handle.
	Service Trolley 3 Shelf	IT- 507.12	73X52X90	* Two shelves with four sides bended upward. * One side pull handle.*
40	Instrument Trolley	IT- 508. 11	56.5x44x93	* Two flat shelves with four sides bended downward. * two side pull handle * Without rails
	Instrument Trolley with Rails	IT- 508-11R	56.5x44x93	* Two flat shelves with four sides bended downward. * two side pull handle*With rails
	Instrument Trolley with Drawer	IT- 508.11D	56.5x44x93	* Two flat shelves with four sides bended downward. * two side pull handle * With rails
	Instrument Trolley with 3 Shelf	IT- 508-12R	56.5x44x93	* Three flat shelves with four sides bended downward and four side rails on each shelves
41	Instrument Trolley	IT-509.10	56.5x44x93	* Two flat shelves with four sides bended downward and three side rails on each shelves
	Instrument Trolley	IT-509.11	56.5x44x93	* Two flat shelves with four sides bended downward and three side rails on each shelves * With drawers
	Monitor Trolley	MT-4W		* Top shelves with four side inclined upwards * With heavy duty stand. * On heavy duty wheels two with brakes.
42	Instrument Trolley	IT- 510.10	51x50x90	* Two flat shelves with and three side rails on each shelves
	Instrument Trolley 1 Drawer	IT- 510.11	51x50x90	* Two flat shelveswith four sides bended downward and three side rails on each shelves * With single drawers
	Instrument Trolley 2 Drawer	IT- 510.21	51x50x90	* Two flat shelveswith four sides bended downward and three side rails on each shelves * With double drawers
	Instrument Trolley 5 Drawer	IT- 510-5D	80x50x90	* Two flat shelves with four sides bended downward and three side rails on each shelves * With Five drawers
43	Cast Cutting Trolley	CC-2S-BB	80x45x90	* With single bowl sink and SS bin * Two shelves and drawers
	Single Wash Basin Trolley	WB-SB-6	42x40x85	* Two stainless steel removable basin with 6 liter capacity
	Double Wash Basin Trolley	WB- DB- 6	78x41x85	* Single stainless steel removable basin with 6 liter capacity
	Nursing Trolley	NT- 506-2D		* With single bowl sink and SS waste bin * Two shelves and drawers
44	Procedure Trolley 3 Drawer	PT- 3D	62x50x86	* Two shelves with three side raised edges and Three, sliding drawers with system
	Critical Care Cart	CCC- B6	60x50x100	* Top plain work surface. * One sliding drawer. * Six nos of removable type wire baskets. * Heavy duty 4" swivel castors, two with brakes.

P No.	Items	Code	Measurement (CM)	Discription
44	Procedure Trolley 5 Drawer	PT- 5D	62x50x90	* Two shelves with three side raised edges and five sliding drawers * Central locking system
45	Cylinder Transport Trolley	CTT- 12		* Frame made of Aluminium/stainless steel tubes & plates * Hold 12cubic feet cylinder
	Cylinder Transport Trolley	CTT 24		* Frame made of Aluminium/stainless steel tubes & plates * Hold 24 cubic feet cylinder
	Cylinder Transport Trolley	CTT- D48		* Frame made of Aluminium/stainless steel tubes & plates * Hold 48 cubic feet cylinder
	Cylinder Transport Trolley- Double	CTT- 24D		* Frame made of Aluminium/stainless steel tubes & plates * Hold 24 cubic feet cylinder-Two cylinders
46	Emergency Cart	EMC. 1	70x50x90	* Emergency Cart-Stainless Steel * Five nos of sliding drawers
	Anesthesia Cart	ANC.1	70x50x90	* Anesthesia Cart-Stainless Steel * Six nos of sliding drawers, different sizes and Tilt in trash box
	Emergency Cart with 4 Drawers & 3 Tilt Bin	EC- 4D-3T	70x50x90	* Four nos of sliding drawers, different size and Centralized locking system. * Three tilt bin on left side"
47	Medication Cart	MC-15-2D	70x50x95	* Twenty five nos of crates on both sides (Total 50) * One sliding drawer at bottom with Central locking system"
	Isolation Cart	IC- 4DP	70x50x95	Four nos of sliding drawers, different sizes with Centralized locking system.
48	CSSD Sterile Transport Trolley	CSD- 120	100x70x120	* Two nos of inner solid shelves. * Two hinged door with lock and keys.
	CSSD Sterile Transport Trolley	CSD- 90	100x70x90	* Two nos of inner solid shelves. * Two hinged door with lock and keys.
	Powder Coated Transport Trolley	CSD- 90P	100x70x90	* Two nos of inner solid shelves. * Two hinged door with lock and keys.
49	Wax Bath	WB	55x40x30	* Double walled body, well insulated stainless steel wax tank * Thermostat to control temperature from 00 to 85 deg.Cen.
	Spica Cast Table	SPCT		* With acrylic support on head portion and removable Stainless steel support and teflon rod for leg portion. * Head portion is height adjustable and base portion is length adjustable.
	Infant Blood Drawing Station	IBDS		* 50mm hight imitation leather upholstery * Powder coated/Stainless Steel frame
50	File Transport Trolley	FTT- 50	75x55x90	* Well type file storage tray 35 cm depth with file holding rods inside and two hinged lids. * With one bottom shelf.
	Suturing Trolley	SUT- TB	80x545x90	* With two compartment * One side with Four sliding drawer and Other sith one removable wash basin and suture waste bin * Two flat shelves with four sides inclined upward, Without side rails.
51	Clean Linen Cart with Hinged Door	CLC- 3H	110x50x128	* Hinged doors with lock and keys. * Two nos of height adjustable inner shelves.

P No.	Items	Code	Measurement (CM)	Discription
51	Paper Dispensing Cart on wheel	PD- 120T	130x47x93	* Four nos of tubes with the catchers to hold the papers. * One side pull handle. * One wire shelf at the bottom.
	Clean Linen Cart with Sliding Door	CLC- 3S	110x50x128	* Sliding doors with lock and keys. * Two nos of height adjustable inner shelves.
52	Loading Trolley with Wire Bottom	LT- WB	80x60x160	* Removable wire shelf and Wire base. * On heavy duty swivel castors two with brakes.
	Loading Trolley with Solid Bottom	LT- SB	80x60x160	* Removable wire shelf and Solid base. * On heavy duty swivel castors two with brakes.
	Mobile Video Cart	MVC- 5S	70x64x170	* Five nos of shelves with one sliding drawer * One control switch in the front. * Four nos of three pin plug points with individual switch.
53	S/S Dustbin 60	DBS- 60	43x31.5x81	* With wheels * Free standing, Pedal operated, Soundless closing top lid. * Less weight removable container
	S/S Dustbin 20	DBS- 20	23x28.5x57	* Free standing, Pedal operated, Soundless closing top lid. * Less weight removable container
	S/S Dustbin Round	DBRD- 10	20, H28.5	Stainless steel round dust bin with 10ltr capacity
	Powder Coated Dustbin 60	PC- DB- 60	43x31.5x81	* With wheels * Free standing, Pedal operated, Soundless closing top lid. * Less weight removable container
	Powder Coated Dustbin 20	PC-DB- 20	23x28.5x57	* Free standing, Pedal operated, Soundless closing top lid. * Less weight removable container
54	Baby Bassinet Cabinet	BBC		Stainless Steel/Powder Coated Cabinet and acrylic Baby Basinet
	Baby Bassinet Trolley	BBT		Stainless Steel/Powder Coated Frame and acrylic Baby Basinet
	Over Bed Table	OV- BT	81x40x75/106	Stainless Steel/ Powder Coated Frame with wooden top
	Sharp Container Trolley	SC- T		Capacity: As per customer requirement
55	Dirty Linen Trolley Double bag with Lid	DL- SB2L	86x47x90	Dirty Linen Trolley Double bag with Lid-Stainless Steel frame and lid with double nylon bag
	Dirty Linen Trolley Single Bag Square with Lid	DL- SB1L	50x56x93	Dirty Linen Trolley Double bag with Lid-Stainless Steel frame and lid with single nylon bag
	Linen Skip with Lid	LS- WL	40x47x90	Linen Skip with Lid-Stainless Steel frame and lid with single nylon bag
56	Dirty Linen Trolley Single Bag square	DL- SB	50x56x93	Dirty Linen Trolley Double bag without Lid-Stainless Steel frame with single nylon bag
	Dirty Linen Trolley Double Bag square	DL-SB2L	86x47x90	Dirty Linen Trolley Double bag with Lid-Stainless Steel frame and lid with double nylon bag
	Dirty Linen Trolley Double Bag Round	DL- RD	50x2-H: 93	Dirty Linen Trolley Double (Round)-Stainless Steel frame and lid with double nylon bag
	Dirty Linen Trolley Single Bag Round	DL- RS	50-H: 93	Dirty Linen Trolley Single (Round)-Stainless Steel frame and lid with double nylon bag
57	Laundry Trolley	CT- 100	120x80x70	Laundry Trolley-stainless steel with castors
	Wire Mesh Trolley	CT- W100	100x80x80	Wire mesh Trolley-stainless steel with castors

P No.	Items	Code	Measurement (CM)	Discription
57	Box Sacks Linen Cart	BSLT- 120	110x60x120	Box sacks Linen Cart-stainless steel with castors
58	Pathology Grossing Work Station	GWSU-SS-150	150x75x190	Stainless Sttel Patology grossing Station
59	Lead Apron Hanger Two Tier on wheel	AP.OWS- 100	120x60x140	* Stainless Steel/Powder Coated * With 10 apron Hanger
	Lead Apron Hanger on Wheel	AP.OW- 100	69x72.5x140	* Stainless Steel/Powder Coated * With 10 apron Hanger
	Lead Apron Hanger wall Mounted	AP.WM- 110	30x60x30	* Stainless Steel/Powder Coated * With 10 apron Hanger
60	Pharmacy Shelving with Lin Bins	2B- WT.LBD		*Powder coated *Slanted/staright shelves
	Pharmacy Shelving with Sliding Drawer	2B-180		* Powder coated * Slanted/staright shelves and drawers
	Pharmacy Shelving	ST- WM.D		*Powder coated *Slanted/staright shelves
61	Heavy Duty Rack Shelving	HD-RS		*With two shelves and heavy duty bross supporting bars *Powder coated
	Pot Rack	PR-100	120x35x15	Pot Rack- Stainless Steel
	Single Tier Wall Shelf with Hook	SW-HK-120	100x25x15	Single Tier Wall Shelf with Hook-Stainless Steel
	Wall Shelf	SW-120	120x53x183	Wall Shelf-Stainless Steel
62	Wire Shelving on wheel	WS-120T	120x53x200	* Wire Shelf-Stainless steel * F*our tire wire Shelves
	Wire Shelving	WS-120T	120x53x200	* Wire Shelf-Stainless steel * Four tire wire Shelves
	Wire Shelving	WS-120H	120x53x200	* Wire Shelf-Stainless steel with wheels * Four tire wire Shelves
	Wire Shelving	WST-5		* Wire Shelf-Stainless steel with Solid bottom * Four tire wire Shelves
63	Lin Bin Holder Panel- Wall Mounted	LBW		* Lin Bin Holder Panel- Wall Mounted- Powder coated
	Lin Bin with Holder panel	LBH		* Lin Bin with Holder panel-Powder coated
	3Lin Bin Holder on Wheel	LBH-WL		* Lin Bin Holder on Wheel-Powder coated
64	Compact Cabinet Storage System	CCSS		* Compact Cabinet Storage System-Powder coated
	Boltfree Shelving System	BFSS		* Boltfree Shelving System-Powder Coated- Four shelves
65	Examination Couch with Cabinet	EC-BC-601	185x62x75	* 50mm hight imitation leather upholstery * Powder coated/Stainless Steel bottom Cabinet
	Examination Couch	EC-601	195x62x75	* 50mm hight imitation leather upholstery * Stainless Steel frame * Knock down fabrication
	Massage Couch	EC-604	200x80x75	* 50mm hight imitation leather upholstery * Stainless Steel frame
66	Examination Couch with Side Rails	EC-602	185x62x75	* 50mm hight imitation leather upholstery * Stainless Steel frame
	Ward Bed	PB-2	210x90x55	* Ward Bed-Powder coated * Without mattress

P No.	Items	Code	Measurement (CM)	Discription
67	Gynecology Examination Couch Electrical	GEC-E	185x62x75	Gynecology Examination Couch -Electrical
	Gynecology Examination Couch	EC- 605	181x61x86	* 50mm hight imitation leather upholstery * Stainless Steel frame
68	Electrical 2 Section Couch	EEC- 610	185x62x75	Electrical 2 Section Couch-Powder coated frame structure
	Hydraulic Examination Couch -3 section	ECH-61H	185x62x75	Hydraulic Examination Couch -3 section - Powder coated Frame
69	Hydraulic Examination Couch with Side Rails	ECH-613	186x62x50-100	* Powder coated Frame * 50mm hight imitation leather upholstery
	Hydraulic Examination Couch with Padded Side Rails	ECH-614	186x62x50-100	* Powder coated Frame * 50mm hight imitation leather upholstery * Padded Side rails
70	Electrical Examination Couch- 4 section	EEC-615		Electrical 4 Section Couch-Powder coated frame structure
	Electrical Examination Couch- 3 section	EEC-614		Electrical 3 Section Couch-Powder coated frame structure
71	Baby Coat	BA- CT	132x62x75	* Powder coated Frame and Two side rails * With mattress
	Bed Side Cabinet	BSL- 50		* ABS construction with One drawer and cabinet
72	Infusion Stand with Nylon Base	IV- PB4		* With Nylon Base and 4 hooks * Height adjustable telescopic stem
	Infusion Stand with Stainless Steel Base	IV-SS4		* With Stainless Steel Base and 4 hooks * Height adjustable telescopic stem
	Heavy Duty Infusion Stand	IV- HD		* With Stainless Steel Base and 4 hooks * Height adjustable telescopic stem
73	Syringe Pump Stand	IV- SP 4		* With Stainless Steel spider leg Base * Height adjustable telescopic stem and Electric Connection
	Mobile Suction Trolley	IV-MST		* With Stainless Steel base * Suction bottle hanging hooks
	Monitor Stand with Wire Basket-MRI Compactable	MS-WB		Monitor Stand with Wire Basket-MRI Compactable 5 Wheels with base * One wire basketand Monitor holder
74	Stainless Steel Gloves Dispenser	GD-SS	38.9x9.5x24.5	Gloves Dispenser- Stainless Steel * 3 gloves dispensing holder
	PPE Holder Powder Coated	PPEH-P	38.9x9.5x24.5	PPE Holder Powder Coated * With gloves holder tissue holder and mask and apron holder
	Tissue Dispenser	TD-P	13x9.5x24.5	Tissue Dispenser- Stainless Steel * 3 Tissue dispensing holder
	Acrylic Gloves Dispenser	GD-AL	13x9.5x24.5	Gloves Dispenser-Acrylic * 3 gloves dispensing holder
	Acrylic PPE Holder	PPE-AL	40.5X15X61	PPE Holder Acrylic * With gloves holder tissue holder and mask and apron holder
75	Examination Light Halogen	EL- HA		Examination Light Halogen
	Magnifying Lamp on Wheel	ML-OW		Magnifying Lamp on Wheel
	Examination Light Wall Mounted	EL- WM		Examination Light Wall Mounted

P No.	Items	Code	Measurement (CM)	Discription
76	Single Step Footstool	FS-SR	40x24x22	* Stainless steel frame * Anti slip rubber mat
	Double Step Footstool	FS- DR	40x52x33	* Stainless steel frame * Anti slip rubber mat
	Single Step Foot Stool- Stainless Steel	FS- SS	40x24x22	* Stainless steel frame * Anti slip chequered plate mat
	Double Step Foot Stool- Stainless Steel	FS- DS	40x52x33	* Stainless steel frame * Anti slip chequered plate mat
77	Single Step Foot Stool with One Side Hand Rail	FS- SR1	40x24x22+70	* Stainless steel frame With side rail * Anti slip rubber mat
	Double Step Foot Stool with Two Side Rail	FS- 2R	40x52x33+70	* Stainless steel frame with two side rail * Anti slip rubber mat
	Triple Step Foot Stool	FS- TR	80x40x56	* Stainless steel frame * Anti slip rubber mat
	Revolving Stool SS Seat	RS- SS		* With stainless steel Seat and Height adjustable stem * Chrome plated/nylon base with/without foot ring
78	Revolving Stool - Low Backrest	RS- LB		* With imitation leather upholstery Seat and Low Backres * Height adjustable stem Chrome plated/nylon base With/without foot ring
	Revolving Stool	RS- PBF		* With imitation leather upholstery Seat without back rest * Height adjustable stem Chrome plated/nylon base With/without foot ring
	Revolving Stool with Side Rest	RS-SDR		* With imitation leather upholstery Seat with side rest * Height adjustable stem Chrome plated/nylon base With/without foot ring
	Revolving Stool with High Backrest	RS-HB		* With imitation leather upholstery Seat and high Backres * Height adjustable stem Chrome plated/nylon base With/without foot ring
79	Blood Donor Chair	BD-M	60X75X125	* Stainless steel frame with imitation leather upholstery * Height adjustable
80	X-ray flm viewer Double Standard	XV-HA 1		X-ray flm viewer Double Standard
	X-ray flm viewer Single Standard	XV- HA 2		X-ray flm viewer Single Standard
	X-ray flm viewer Double Slim	XV- TB 2		X-ray flm viewer Double Slim
	X-ray flm viewer Single Slim	XV- TB 1		X-ray flm viewer Single Slim
81	3 Fold Ward Screen	WS. 3 F	W:184, H:175	* Stainless steel /powder coated frame * Fire retardent PVC cloth
	4 Fold Ward Screen	WS. 4F	W:245, H:175	* Stainless steel /powder coated frame * Fire retardent PVC cloth
	Urine Bag Holder	UBH	L:14, H:31	Urine Bag Holder-Stainless Steel
	Brush Dispenser	BD-80	92x16x80	Brush Dispenser-Stainless Steel
82	Bed pan/Urinal Holder	BP- HW		Bed pan/Urinal Holder-Stainless Steel
	Kick Bowl	KB- 6	40,H:19	* Five leg stainless steel tube base with 5 vertical supports * Bowl with 6liters capacity. * Mobile on 5 wheels

P No.	Items	Code	Measurement (CM)	Discription
82	Kick Bucket	KBK- 20	40,H:36	* Five leg stainless steel tube base with 5 vertical supports * Bucket with lid * Mobile on 5 wheels
	Shoe Rack on Wheel	SR- OW	120x35x140	Shoe Rack -Stainless steel with Wheels
	Shoe Rack Standard	SR- FS	80x25x30	Shoe Rack Standard
83	Wire Basket	WBKC	40x30x10	* With Mounting bracket * Frame made of 8mm rods and Basket made of 6mm rods
	Wire Basket	WBK	40x30x10	* With Mounting bracket * Frame made of 8mm rods and Basket made of 6mm rods
	Perforated Sheet Basket	PSB	40x30x10	* With Mounting bracket * Made with stainless steel sheets
	Wire Basket with Wall Mounting Bracket	WB-WMB	40x30x10	* With Mounting bracket * Frame made of 8mm rods and Basket made of 6 mm rods
84	Soap Dispenser- Normal	SCD- N		Stainless steel Soap Dispenser
	Soap Dispenser Elbow Operated	SD- EL		Stainless steel elbow operated Soap Dispenser
	Soap Dispenser Table Mounted	SD-TM		Stainless steel table mounted Soap Dispenser
	Patient Chart Holder - Plastic	CHP- A4	23x34	Plastic Patient chart holder
	Patient Chart Holder- Stainless Steel	CHS-A4	23x34	Stainless steel Patient chart holder
85	Walking Forearm Crutches- Closed Type	WFC- C		
	Walking Crutches- Open Type	WFC- O		
	Walking Stick Foldable	WS- F		
	Adult Crutches - Underarm	CU-A		
	Child Crutches- Underarm	CU-C		
	Quad Walking Stick Straight Handle	QWS-S		
	Quad Working Stick Scurved Handle	QWS- C		
86	Walking Frame Foldable with Wheel	WFF-W		
	Walking Frame Foldable Standard	WFF-S		
	Walking Frame with Seat on Two Wheel	WF- S- 2W		
	Walking Frame with Seat on Four Wheel	WF- S- 4W		
87	Wheel Chair with Padded Seat Non Foldable.	WC- P-NF		
	Wheel Chair Foldable	WC- F		
88	Catheter Basket -Long	CB- L		Made of stainless steel Perforated Sheet.
	Catheter Basket - Wiremesh	CB- W		Made of stainless steel rods
	Catheter Basket - Divided	CB- DV		Made of stainless steel Perforated Sheet with center seperation

P No.	Items	Code	Measurement (CM)	Discription
89	Cadaver washing Trolley	CW.T	200x65x80	Stainless steel cadaver washing trolley * On heavy duty wheels two with breaks
	Cadaver closed Transport Trolley	CCT	200x65x80+40	Stainless steel cadaver trolley with cover * On heavy duty wheels two with breaks
	Autopsy Table	CAT-SS	200x65x90	Stainless steel Autopsy table with sink and shower * With tools hanger
90	Pass through Chamber Sliding Door	PA-CH-S	80x60x130	* Mechanically operated hatches, stainless steel framed ,verically sliding clear acrylic glass doors. * One door can be open at a time
	Pass through Chamber Hinged Door	PA-CH-H	60x60x60	* Mechanically operated hatches, stainless steel framed clear acrylic glass hinged doors. * One door can be open at a time
	Pass through Window Hinged Door	PAW-H		* Pass Through window with stainless steel framed clear acrylic glass hinged doors. * One door can be open at a time