

claricuff™

Blood Pressure Cuffs

Claricuff™ is an unrivaled level of exceptional product quality. The Neonatal Cuffs are especially made with safe materials and manufactured carefully for the health and well-being of Newborns and Babies.

Features

- Claricuff™ Neonatal Blood Pressure Cuffs are made for single patient use (disposable style)
- Claricuff™ Neonatal Blood Pressure Cuffs are made to AHA/ AAMI Standards
- Claricuff™ Neonatal Blood Pressure Cuffs have different colors to aid in sizing
- Claricuff™ Neonatal Blood Pressure Cuffs are available in all standard sizes
- Claricuff™ Neonatal Blood Pressure Cuffs are made with highest quality materials
- Claricuff™ Neonatal Blood Pressure Cuffs do not contain natural latex rubber, DEHP or PVC



Claricuff™ Neonatal

Description	Singe Tube	Double Tube	Cuff Size
Claricuff™ BP Cuff TPU Neonatal #1 with ML fitting, 10/BG	CNN1ST-10	CNN1DT-10	size 1
Claricuff™ BP Cuff TPU Neonatal #2 with ML fitting, 10/BG	CNN2ST-10	CNN2DT-10	size 2
Claricuff™ BP Cuff TPU Neonatal #3 with ML fitting, 10/BG	CNN3ST-10	CNN3DT-10	size 3
Claricuff™ BP Cuff TPU Neonatal #4 with ML fitting, 10/BG	CNN4ST-10	CNN4DT-10	size 4
Claricuff™ BP Cuff TPU Neonatal #5 with ML fitting, 10/BG	CNN5ST-10	CNN5DT-10	size 5

Cuffs come standard with the ML fitting. It is offered with single tube or double tube.



Tapmedic LLC
9424 Baymeadows Rd, Suite 250
Jacksonville, FL 32256 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com

claricuff™

Blood Pressure Cuffs

Claricuff™ is an unrivaled level of exceptional product quality. The Claricuff™ Blood Pressure Cuffs are especially made with safe materials and manufactured carefully for the health and well-being of the patient.

Features

- Claricuff™ Blood Pressure Cuffs are made for single patient use (disposable style) and multiple patient use (reusable style)
- Claricuff™ Blood Pressure Cuffs are made to AHA/ AAMI Standards
- Claricuff™ Blood Pressure Cuffs have different colors to aid in sizing
- Claricuff™ Blood Pressure Cuffs are available in all standard sizes to cover all patients
- Claricuff™ Blood Pressure Cuffs are made with highest quality materials
- Claricuff™ Blood Pressure Cuffs do not contain natural latex rubber, DEHP or PVC (reusable only)



Reusable & Disposable Cuffs

Fittings



Screw (DM)



Male Luer (ML)



Twist Lock (MQ)



Female Luer (FL)




Bayonet (HP)



Male Twist Lock (MQML)



Metal Bayonet (MHP)



Male & Female (MQMF)

Description	Singe Tube	Double Tube	Cuff Size
Claricuff™ BP Cuff Reusable Infant w/ fitting, 5/BG	RCNST0814xx-05	RCNDT0814xx-05	size 6
Claricuff™ BP Cuff Reusable Child w/ fitting, 5/BG	RCNST1320xx-05	RCNDT1320xx-05	size 8
Claricuff™ BP Cuff Reusable Small Adult w/ fitting, 5/BG	RCNST1826xx-05	RCNDT1826xx-05	size 10
Claricuff™ BP Cuff Reusable Adult w/ fitting, 5/BG	RCNST2635xx-05	RCNDT2635xx-05	size 12
Claricuff™ BP Cuff Reusable Adult Long w/ fitting, 5/BG	RCNST2938xx-05	RCNDT2938xx-05	size 13
Claricuff™ BP Cuff Reusable Large Adult w/ fitting, 5/BG	RCNST3242xx-05	RCNDT3242xx-05	size 14
Claricuff™ BP Cuff Reusable Large Adult Long w/ fitting, 5/BG	RCNST3544xx-05	RCNDT3544xx-05	size 15
Claricuff™ BP Cuff Reusable Thigh w/ fitting, 5/BG	RCNST4250xx-05	RCNDT4250xx-05	size 16
Description	Singe Tube	Double Tube	Cuff Size
Claricuff™ BP Cuff Disposable Infant w/ fitting, 10/BG	CSST0814xx-10	CSDT0814xx-10	size 6
Claricuff™ BP Cuff Disposable Child w/ fitting, 10/BG	CSST1320xx-10	CSDT1320xx-10	size 8
Claricuff™ BP Cuff Disposable Small Adult w/ fitting, 10/BG	CSST1826xx-10	CSDT1826xx-10	size 10
Claricuff™ BP Cuff Disposable Adult w/ fitting, 10/BG	CSST2635xx-10	CSDT2635xx-10	size 12
Claricuff™ BP Cuff Disposable Adult Long w/ fitting, 10/BG	CSST2938xx-10	CSDT2938xx-10	size 13
Claricuff™ BP Cuff Disposable Large Adult w/ fitting, 10/BG	CSST3242xx-10	CSDT3242xx-10	size 14
Claricuff™ BP Cuff Disposable Large Adult Long w/ fitting, 10/BG	CSST3544xx-10	CSDT3544xx-10	size 15
Claricuff™ BP Cuff Disposable Thigh w/ fitting, 10/BG	CSST4250xx-10	CSDT4250xx-10	size 16

To order cuffs for a specific monitor, add the appropriate suffix to the product number.



Claricuff™ BP Cuff
Disposable



Claricuff™ BP Cuff
Reusable



Tapmedic LLC
9424 Baymeadows Rd, Suite 250
Jacksonville, FL 32256 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com



clariflo™
Pressure Infuser

Clariflo™ Flex

Designed for rapid IV infusion and pressure monitoring procedures.



www.tapmedic.com | info@tapmedic.com

ADVANTAGES

Quality you can count on

- No natural rubber latex
- CE Mark
- Durable, state of the art nylon bag
- Leak proof, premium seams for tear resistance
- Product needs minimal operator training
- 100% leak tested after manufacturing

Product features for easy use

- Integrated fluid bag hook
- Large hand bulb for quick inflation
- Security clamp
- IV pole loop with extended suspension
- No-crimp tubing

Choices: we have what you want

- 2 lines of products: Clariflo™ and Clariflo™ Flex
- 3 bag sizes: 500ml, 1000ml, 3000ml
- Disposable / reusable versions
- Piston gauge or aneroid gauge
- Stopcock or thumbwheel valve



Reinforced **see-through** TPU sleeve allows proper label observation.

OPTIONS

Pressure Gauges

The Clariflo™ is available in 2 versions of pressure gauges:



PISTON GAUGE

- Pressure range 0-300mmHg
- 360° gauge visibility
- Bleed valve prevents over-inflation



ANEROID GAUGE

- Pressure range 0-300mmHg
- Large, easy to read display
- Exact pressure control

Pressure Control Valves

The Clariflo™ offers a choice of 2 pressure control valves:



3-WAY STOPCOCK

- Precise control of pressure
- Conveniently operated by one hand



THUMBWHEEL VALVE

- Instant control of pressure
- Open or close mode

Clariflo™ Flex

The newest line is the Flex. Clariflo™ Flex provides the separate cuff and inflator. The concept is - use and keep the inflator. The cuff can be exchanged. It is semi- reusable and a reduction in bioburden.

4 styles available:

- Semi- disposable with piston gauge
- Semi- disposable with aneroid gauge
- Semi- reusable with piston gauge
- Semi- reusable with aneroid gauge

➤ 3 bag sizes available



	Part Number	Description	Box Qty
disposable	500PGSFlex	Clariflo™ Flex Classic - 500ml bag with piston gauge and stopcock	2x10/box
	1000PGSFlex	Clariflo™ Flex Classic - 1000ml bag with piston gauge and stopcock	2x10/box
	3000PGSFlex	Clariflo™ Flex Classic - 3000ml bag with piston gauge and stopcock	2x5/box
	500AGSFlex	Clariflo™ Flex Classic - 500ml bag with aneroid gauge and stopcock	2x10/box
	1000AGSFlex	Clariflo™ Flex Classic - 1000ml bag with aneroid gauge and stopcock	2x10/box
	3000AGSFlex	Clariflo™ Flex Classic - 3000ml bag with aneroid gauge and stopcock	2x5/box
reusable	RU500PGSFlex	Clariflo™ Flex Ultimate - 500ml bag with piston gauge and stopcock	2x10/box
	RU1000PGSFlex	Clariflo™ Flex Ultimate - 1000ml bag with piston gauge and stopcock	2x10/box
	RU3000PGSFlex	Clariflo™ Flex Ultimate - 3000ml bag with piston gauge and stopcock	2x5/box
	RU500AGSFlex	Clariflo™ Flex Ultimate - 500ml bag with aneroid gauge and stopcock	2x10/box
	RU1000AGSFlex	Clariflo™ Flex Ultimate - 1000ml bag with aneroid gauge and stopcock	2x10/box
	RU3000AGSFlex	Clariflo™ Flex Ultimate - 3000ml bag with aneroid gauge and stopcock	2x5/box



Clariflo™ Flex Ultimate

Designed for semi- reuse, easy cleaning and made with a durable, reinforced blue nylon bag.



Clariflo™ Flex Classic

The durable white nylon bag is made for single patient use and semi- disposable.

Clariflo™ Flex Components



bag only
(reusable)



aneroid gauge with stopcock



piston gauge with stopcock



bag only
(disposable)



aneroid gauge with stopcock



piston gauge with stopcock

	Part Number	Description	Box Qty
disposable	500Flex	Clariflo™ Flex Classic - 500ml bag	2x10/box
	1000Flex	Clariflo™ Flex Classic- 1000ml bag	2x10/box
	3000Flex	Clariflo™ Flex Classic - 3000ml bag	2x5/box
	PGSFlex	Clariflo™ Flex Classic - piston gauge and stopcock	4x5/box
	AGSFlex	Clariflo™ Flex Classic - aneroid gauge and stopcock	4x5/box
reusable	RU500Flex	Clariflo™ Flex Ultimate - 500ml bag	2x10/box
	RU1000Flex	Clariflo™ Flex Ultimate - 1000ml bag	2x10/box
	RU3000Flex	Clariflo™ Flex Ultimate - 3000ml bag	2x5/box
	RUPGSFlex	Clariflo™ Flex Ultimate - piston gauge and stopcock	4x5/box
	RUAGSFlex	Clariflo™ Flex Ultimate - aneroid gauge and stopcock	4x5/box



Tapmedic is a global provider of high-quality medical equipment and healthcare devices, with a broad range of products used for rapid infusion, NIBP diagnostics, blood pressure cuffs and monitoring, and patient care.

Tapmedic LLC
9424 Baymeadows Rd, Suite 250
Jacksonville, FL 32256 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com

Tapmedic reserves the right to make changes to this sheet and the product at any time without notice. Clariflo™ is a registered trademark of Tapmedic.



clariflo™
Pressure Infuser

Clariflo™ Classic/ Ultimate/ Flex

Designed for rapid IV infusion and
pressure monitoring procedures.

Tapmedic is a global provider of high-quality medical equipment and healthcare devices, with a broad range of products used for rapid infusion, NIBP diagnostics, blood pressure cuffs and monitoring, and patient care.



Tapmedic LLC
9424 Baymeadows Rd, Suite 250
Jacksonville, FL 32256 USA
T +1 904-351-0244 F +1 904-351-0245
www.tapmedic.com | info@tapmedic.com



ADVANTAGES

Quality you can count on

- No natural rubber latex
- CE Mark
- Durable, state of the art nylon bag
- Leak proof, premium seams for tear resistance
- Product needs minimal operator training
- 100% leak tested after manufacturing

Product features for easy use

- Integrated fluid bag hook
- Large hand bulb for quick inflation
- Security clamp
- IV pole loop with extended suspension
- No-crimp tubing

Choices: we have what you want

- 2 lines of products: Clariflo™ and Clariflo™ Flex
- 3 bag sizes: 500ml, 1000ml, 3000ml
- Disposable / reusable versions
- Piston gauge or aneroid gauge
- Stopcock or thumbwheel valve



Reinforced **see-through** TPU sleeve allows proper label observation.

OPTIONS

Pressure Gauges

The Clariflo™ is available in 2 versions of pressure gauges:



PISTON GAUGE

- Pressure range 0-300mmHg
- 360° gauge visibility
- Bleed valve prevents over-inflation



ANEROID GAUGE

- Pressure range 0-300mmHg
- Large, easy to read display
- Exact pressure control

Pressure Control Valves

The Clariflo™ offers a choice of 2 pressure control valves:



3-WAY STOPCOCK

- Precise control of pressure
- Conveniently operated by one hand



THUMBWHEEL VALVE


- Instant control of pressure
- Open or close mode

Clariflo™ Classic

The Clariflo™ Classic features a durable white nylon bag. It is made for single patient use and completely disposable, eliminating the need to clean or repair.

disposable	Part Number	Description	Box Qty
	500PGS	Clariflo™ Classic - 500ml bag with piston gauge and stopcock	2x10/box
	1000PGS	Clariflo™ Classic - 1000ml bag with piston gauge and stopcock	2x10/box
	3000PGS	Clariflo™ Classic - 3000ml bag with piston gauge and stopcock	2x5/box
	500PGT	Clariflo™ Classic - 500ml bag with piston gauge and thumbwheel valve	2x10/box
	1000PGT	Clariflo™ Classic - 1000ml bag with piston gauge and thumbwheel valve	2x10/box
	3000PGT	Clariflo™ Classic - 3000ml bag with piston gauge and thumbwheel valve	2x5/box
	500AGS	Clariflo™ Classic - 500ml bag with aneroid gauge and stopcock	2x10/box
	1000AGS	Clariflo™ Classic - 1000ml bag with aneroid gauge and stopcock	2x10/box
	3000AGS	Clariflo™ Classic - 3000ml bag with aneroid gauge and stopcock	2x5/box
	500AGT	Clariflo™ Classic - 500ml bag with aneroid gauge and thumbwheel valve	2x10/box
	1000AGT	Clariflo™ Classic - 1000ml bag with aneroid gauge and thumbwheel valve	2x10/box
	3000AGT	Clariflo™ Classic - 3000ml bag with aneroid gauge and thumbwheel valve	2x5/box

3 bag sizes available



Clariflo™ Ultimate

The Clariflo™ Ultimate is designed for reuse, easy cleaning and made with a durable, reinforced blue nylon bag.

reusable	Part Number	Description	Box Qty
	RU500PGS	Clariflo™ Ultimate - 500ml bag with piston gauge and stopcock	2x10/box
	RU1000PGS	Clariflo™ Ultimate - 1000ml bag with piston gauge and stopcock	2x10/box
	RU3000PGS	Clariflo™ Ultimate - 3000ml bag with piston gauge and stopcock	2x5/box
	RU500PGT	Clariflo™ Ultimate - 500ml bag with piston gauge and thumbwheel valve	2x10/box
	RU1000PGT	Clariflo™ Ultimate - 1000ml bag with piston gauge and thumbwheel valve	2x10/box
	RU3000PGT	Clariflo™ Ultimate - 3000ml bag with piston gauge and thumbwheel valve	2x5/box
	RU500AGS	Clariflo™ Ultimate - 500ml bag with aneroid gauge and stopcock	2x10/box
	RU1000AGS	Clariflo™ Ultimate - 1000ml bag with aneroid gauge and stopcock	2x10/box
	RU3000AGS	Clariflo™ Ultimate - 3000ml bag with aneroid gauge and stopcock	2x5/box
	RU500AGT	Clariflo™ Ultimate - 500ml bag with aneroid gauge and thumbwheel valve	2x10/box
	RU1000AGT	Clariflo™ Ultimate - 1000ml bag with aneroid gauge and thumbwheel valve	2x10/box
	RU3000AGT	Clariflo™ Ultimate - 3000ml bag with aneroid gauge and thumbwheel valve	2x5/box

3 bag sizes available



Clariflo™ Flex

The newest line is the Flex. Clariflo™ Flex provides the separate cuff and inflator. The concept is - use and keep the inflator. The cuff can be exchanged. It is semi-reusable and a reduction in bioburden.

4 styles available:

- Semi-disposable with piston gauge
- Semi-disposable with aneroid gauge
- Semi-reusable with piston gauge
- Semi-reusable with aneroid gauge

➤ 3 bag sizes available



	Part Number	Description	Box Qty
disposable	500PGSFlex	Clariflo™ Flex Classic - 500ml bag with piston gauge and stopcock	2x10/box
	1000PGSFlex	Clariflo™ Flex Classic - 1000ml bag with piston gauge and stopcock	2x10/box
	3000PGSFlex	Clariflo™ Flex Classic - 3000ml bag with piston gauge and stopcock	2x5/box
	500AGSFlex	Clariflo™ Flex Classic - 500ml bag with aneroid gauge and stopcock	2x10/box
	1000AGSFlex	Clariflo™ Flex Classic - 1000ml bag with aneroid gauge and stopcock	2x10/box
	3000AGSFlex	Clariflo™ Flex Classic - 3000ml bag with aneroid gauge and stopcock	2x5/box
reusable	RU500PGSFlex	Clariflo™ Flex Ultimate - 500ml bag with piston gauge and stopcock	2x10/box
	RU1000PGSFlex	Clariflo™ Flex Ultimate - 1000ml bag with piston gauge and stopcock	2x10/box
	RU3000PGSFlex	Clariflo™ Flex Ultimate - 3000ml bag with piston gauge and stopcock	2x5/box
	RU500AGSFlex	Clariflo™ Flex Ultimate - 500ml bag with aneroid gauge and stopcock	2x10/box
	RU1000AGSFlex	Clariflo™ Flex Ultimate - 1000ml bag with aneroid gauge and stopcock	2x10/box
	RU3000AGSFlex	Clariflo™ Flex Ultimate - 3000ml bag with aneroid gauge and stopcock	2x5/box

Clariflo™ Flex Components

	Part Number	Description	Box Qty
disposable	500Flex	Clariflo™ Flex Classic - 500ml bag	2x10/box
	1000Flex	Clariflo™ Flex Classic - 1000ml bag	2x10/box
	3000Flex	Clariflo™ Flex Classic - 3000ml bag	2x5/box
	PGSFlex	Clariflo™ Flex Classic - piston gauge and stopcock	4x5/box
	AGSFlex	Clariflo™ Flex Classic - aneroid gauge and stopcock	4x5/box
reusable	RU500Flex	Clariflo™ Flex Ultimate - 500ml bag	2x10/box
	RU1000Flex	Clariflo™ Flex Ultimate - 1000ml bag	2x10/box
	RU3000Flex	Clariflo™ Flex Ultimate - 3000ml bag	2x5/box
	RUPGSFlex	Clariflo™ Flex Ultimate - piston gauge and stopcock	4x5/box
	RUAGSFlex	Clariflo™ Flex Ultimate - aneroid gauge and stopcock	4x5/box

Tapmedic reserves the right to make changes to this sheet and the product at any time without notice. Clariflo™ is a registered trademark of Tapmedic.