




Professional
GAS THERAPY SYSTEMS
SUCTION & OXYGEN THERAPY



Professional
GAS THERAPY SYSTEMS
SUCTION & OXYGEN THERAPY

Manufactured & Distributed by



MORRIS PRECISION TECHNOLOGY CO., LTD.

No. 7-1, Kung Yeh Rd., Kuan Tien Ind. Dist., Tainan City, Taiwan 72042

Tel: 886-6-6991119 | Fax: 886-6-6991239

Email: iosales@ioxygen.com

Website: www.ioxygen.com

MI-2024-JAN-1 © 2024 MORRIS PRECISION TECHNOLOGY CO., LTD.





Reliable
Precision
High Quality

MORRIS PRECISION TECHNOLOGY CO., LTD., a Taiwanese company set up in 1977, has been designing and manufacturing a wide range of medical oxygen therapy and suction devices since 1997.

Thanks to its strong team of engineers and technicians supported by a well-trained and responsive team of sales representatives, MORRIS has built up an enviable and worldwide reputation as a medical professional. Our markets spread all over the world, providing its ability to meet the needs of the various markets in which it operates.



The company is **ISO 9001**, **ISO 13485**, and good manufacturing practices(**QMS**) certified. The products are **CE** marked and individually numbered for perfect traceability.

We are pleased to invite you to share MORRIS' recent developments through this new catalogue allowing you to select from our wide range of products.

INDEX

_ SR SERIES VACUUM REGULATORS	1
_ SAFETY TRAP, SUCTION LINER & BOTTLES	5
_ FL SERIES FLOWMETER REGULATORS	7
_ LM SERIES FLOWGAUGE REGULATORS	17
_ PS SERIES PRESET REGULATORS	23
_ PR SERIES PRESSURE REGULATORS	31
_ MF & MFB SERIES EMS REGULATORS	35
_ TF SERIES TUBE STYLE FLOWMETERS	37
_ MFC SERIES CLICK STYLE FLOWMETERS	45
_ IOF SERIES QUICK CONNECTS & ADAPTORS	47
_ 103 & SR-AG SERIES GAUGES	53
_ HUMIDIFIERS & ACCESSORIES	55

SR603C-760A-F2-VT9003



SR SERIES


Vacuum Regulators

The purpose of IOxygen SR Series Vacuum Regulator is to measure and adjust the vacuum level. It is used to remove substances from the patient's body during surgical and other procedures. The vacuum regulator should be connected to a wall connection, using a direct probe or through a rail mounting system.

FEATURES:

VACUUM REGULATORS

1



SR603C-760A-F2-XX




SR602C-160A-F2-XX



SR603C-300A-F2-XX

HOW TO ORDER:




2

SPECIFICATIONS:

Model Number	Mode	Operating Ranges	Max Flow Rates (LPM)	Wall Connection	Patient Connection
SR603C-760(A)-F2-00	3 Mode Continuous (REG/OFF/FULL)	0 ~ -760mmHg (Anti-) clockwise	0~80	1/8 NPT Female	None
SR603C-300(A)-BP-DM	3 Mode Continuous (REG/OFF/FULL)	0 ~ -300mmHg (Anti-) clockwise	0~40	British Probe	DISS Male
SR602C-160(A)-OM-VT9003D	2 Mode Continuous (REG/OFF)	0 ~ -160mmHg (Anti-) clockwise	0~40	Ohmeda Male	VT9003D


SR888 SERIES
Vacuum Regulators

FEATURES:



SR888-760A-FM-VT9003D

HOW TO ORDER:




SPECIFICATIONS:

Model Number	Operating Ranges	Max Flow Rates (LPM)	Wall Connection	Patient Connection
SR888-760(A)-F2-DM	0 ~ -760mmHg (Anti-) clockwise	0~80	1/8 NPT Female	DISS Male
SR888-300(A)-FM-TN	0 ~ -300mmHg (Anti-) clockwise	0~80	French Male	Tubing Nipple
SR888-160(A)-JM-VT9003D	0 ~ -160mmHg (Anti-) clockwise	0~40	Japanese Male	VT9003D

Safety Trap

VACUUM SAFETY TRAP BOTTLE

The VT9000 series safety trap protects suction/vacuum regulators and central vacuum systems from accidental contamination by liquid overflow from suction collection canisters. This unique design assures positive Shut-off at low or line vacuum conditions. Various configurations are available for attachment to vacuum regulator or central vacuum system.



Suction Liner & Bottles

+ Purpose:

The bottle collects the waste liquid from the patient who is performing the surgical operation. The bottle lid has the connector for connection tube, suction tube from patient and discharging the waste liquid. It can provide more safety to the operation.



+ FEATURES:

- + Connect with vacuum equipment.
- + Single or multiple set-ups high capacity, multiple systems can save the time for replacing the bag and ensure the surgical operation won't be interrupted.

Bottle:



Liner:



Model Number	Specification	Model Number	Specification
BVT-200A-R	2000ml PC With green holder	BSD-200A	2000ml EVA With valve
BVT-150A-R	1500ml PC With green holder	BSD-150A	1500ml EVA With valve



106RM - FL SERIES

Single stage / Flowmeter regulators

A regulator is a control valve that reduces the input pressure of medical gas to a desired value at its output. The regulators are used for medical gas and can be an integral device with an output pressure setting, a restrictor and a sensor all in the one body, or consist of a separate.

FEATURES:

2" diameter inlet pressure gauge.

Maximum inlet pressure 3000psi(ANSI)B-40.1/210 Bar (CE).

Sintered filter inlet for additional safety and longer life.

Polycarbonate inner and outer tubes provide durability & safety.

Outlet flow control valve 90° or 180° are available.

Neoprene diaphragm with internal safety valve protects against over-pressurization.

Chrome-plated brass body with all brass high pressure chamber.

106RM-15FL

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:

106RM- 15 FL- XXX-SET

Flow rate

Including 1, 3.5, 8, 10, 15, 50 (LPM)

Means flowmeter regulator

Inlet connection

CGA540, CGA540HT, CGA870, CGA992, CGA992HT, CGA346, CGA950, CGA320, CGA320HT and CGA940.

SET

"SET" means the whole set including the combination of accessories.

SPECIFICATIONS:


Model	Gas	Flow Range (LPM)	Inlet Connection
106RM-3.5FL-540	OXYGEN	1/2~3.5	CGA-540 Nut and Nipple
106RM-3.5FL-540HT	OXYGEN	1/2~3.5	CGA-540 Hand tight Nut and Nipple
106RM-3.5FL-870	OXYGEN	1/2~3.5	CGA-870 Yoke (Pin index)
106RM-3.5FL-992	OXYGEN	1/2~3.5	CGA-992 Nut and Nipple
106RM-3.5FL-992HT	OXYGEN	1/2~3.5	CGA-992 Hand tight Nut and Nipple
106RM-10FL-540	OXYGEN	0.2~10	CGA-540 Nut and Nipple
106RM-10FL-540HT	OXYGEN	0.2~10	CGA-540 Hand tight Nut and Nipple
106RM-10FL-870	OXYGEN	0.2~10	CGA-870 Yoke (Pin index)
106RM-10FL-992	OXYGEN	0.2~10	CGA-992 Nut and Nipple
106RM-10FL-992HT	OXYGEN	0.2~10	CGA-992 Hand tight Nut and Nipple
106RM-15FL-540	OXYGEN	1~15	CGA-540 Nut and Nipple
106RM-15FL-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
106RM-15FL-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
106RM-15FL-992	OXYGEN	1~15	CGA-992 Nut and Nipple
106RM-15FL-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple

101RM - FL SERIES

Single stage / Flowmeter regulators

FLOWMETER REGULATORS


9



FEATURES:

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:




SPECIFICATIONS:

Model	Gas	Flow Range (LPM)	Inlet Connection
101RM-3.5FL-540	OXYGEN	1/2~3.5	CGA-540 Nut and Nipple
101RM-3.5FL-540HT	OXYGEN	1/2~3.5	CGA-540 Hand tight Nut and Nipple
101RM-3.5FL-870	OXYGEN	1/2~3.5	CGA-870 Yoke (Pin index)
101RM-3.5FL-992	OXYGEN	1/2~3.5	CGA-992 Nut and Nipple
101RM-3.5FL-992HT	OXYGEN	1/2~3.5	CGA-992 Hand tight Nut and Nipple
101RM-10FL-540	OXYGEN	0.2~10	CGA-540 Nut and Nipple
101RM-10FL-540HT	OXYGEN	0.2~10	CGA-540 Hand tight Nut and Nipple
101RM-10FL-870	OXYGEN	0.2~10	CGA-870 Yoke (Pin index)
101RM-10FL-992	OXYGEN	0.2~10	CGA-992 Nut and Nipple
101RM-10FL-992HT	OXYGEN	0.2~10	CGA-992 Hand tight Nut and Nipple
101RM-15FL-540	OXYGEN	1~15	CGA-540 Nut and Nipple
101RM-15FL-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
101RM-15FL-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
101RM-15FL-992	OXYGEN	1~15	CGA-992 Nut and Nipple
101RM-15FL-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple

10

FR120-FL SERIES


Single stage / Flowmeter regulators



FEATURES:

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:



SPECIFICATIONS:


Model	Gas	Flow Range (LPM)	Inlet Connection
FR120-3.5FL-540	OXYGEN	1/2~3.5	CGA-540 Nut and Nipple
FR120-3.5FL-540HT	OXYGEN	1/2~3.5	CGA-540 Hand tight Nut and Nipple
FR120-3.5FL-870	OXYGEN	1/2~3.5	CGA-870 Yoke (Pin index)
FR120-3.5FL-992	OXYGEN	1/2~3.5	CGA-992 Nut and Nipple
FR120-3.5FL-992HT	OXYGEN	1/2~3.5	CGA-992 Hand tight Nut and Nipple
FR120-10FL-540	OXYGEN	0.2~10	CGA-540 Nut and Nipple
FR120-10FL-540HT	OXYGEN	0.2~10	CGA-540 Hand tight Nut and Nipple
FR120-10FL-870	OXYGEN	0.2~10	CGA-870 Yoke (Pin index)
FR120-10FL-992	OXYGEN	0.2~10	CGA-992 Nut and Nipple
FR120-10FL-992HT	OXYGEN	0.2~10	CGA-992 Hand tight Nut and Nipple
FR120-15FL-540	OXYGEN	1~15	CGA-540 Nut and Nipple
FR120-15FL-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
FR120-15FL-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
FR120-15FL-992	OXYGEN	1~15	CGA-992 Nut and Nipple
FR120-15FL-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple

HOW TO ORDER:

YR88 - FL SERIES

Single stage / Flowmeter regulators

YR88- 15 FL- XXX-SET




SPECIFICATIONS:

Model	Gas	Flow Range (LPM)	Inlet Connection
YR88-3.5FL-540	OXYGEN	1/2~3.5	CGA-540 Nut and Nipple
YR88-3.5FL-540HT	OXYGEN	1/2~3.5	CGA-540 Hand tight Nut and Nipple
YR88-3.5FL-870	OXYGEN	1/2~3.5	CGA-870 Yoke (Pin index)
YR88-3.5FL-992	OXYGEN	1/2~3.5	CGA-992 Nut and Nipple
YR88-3.5FL-992HT	OXYGEN	1/2~3.5	CGA-992 Hand tight Nut and Nipple
YR88-10FL-540	OXYGEN	0.2~10	CGA-540 Nut and Nipple
YR88-10FL-540HT	OXYGEN	0.2~10	CGA-540 Hand tight Nut and Nipple
YR88-10FL-870	OXYGEN	0.2~10	CGA-870 Yoke (Pin index)
YR88-10FL-992	OXYGEN	0.2~10	CGA-992 Nut and Nipple
YR88-10FL-992HT	OXYGEN	0.2~10	CGA-992 Hand tight Nut and Nipple
YR88-15FL-540	OXYGEN	1~15	CGA-540 Nut and Nipple
YR88-15FL-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
YR88-15FL-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
YR88-15FL-992	OXYGEN	1~15	CGA-992 Nut and Nipple
YR88-15FL-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple

FLOWMETER REGULATORS

13



FEATURES:


- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.

14

YR90-FL SERIES

Single stage / Flowmeter regulators


FLOWMETER REGULATORS



FEATURES:

- + The standard preset outlet pressure is 50psi/3.5bar (ANSI) or 60psi/4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:




SPECIFICATIONS:

Model	Gas	Flow Range (LPM)	Inlet Connection
YR90-3.5FL-540	OXYGEN	1/2~3.5	CGA-540 Nut and Nipple
YR90-3.5FL-540HT	OXYGEN	1/2~3.5	CGA-540 Hand tight Nut and Nipple
YR90-3.5FL-870	OXYGEN	1/2~3.5	CGA-870 Yoke (Pin index)
YR90-3.5FL-992	OXYGEN	1/2~3.5	CGA-992 Nut and Nipple
YR90-3.5FL-992HT	OXYGEN	1/2~3.5	CGA-992 Hand tight Nut and Nipple
YR90-10FL-540	OXYGEN	0.2~10	CGA-540 Nut and Nipple
YR90-10FL-540HT	OXYGEN	0.2~10	CGA-540 Hand tight Nut and Nipple
YR90-10FL-870	OXYGEN	0.2~10	CGA-870 Yoke (Pin index)
YR90-10FL-992	OXYGEN	0.2~10	CGA-992 Nut and Nipple
YR90-10FL-992HT	OXYGEN	0.2~10	CGA-992 Hand tight Nut and Nipple
YR90-15FL-540	OXYGEN	1~15	CGA-540 Nut and Nipple
YR90-15FL-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
YR90-15FL-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
YR90-15FL-992	OXYGEN	1~15	CGA-992 Nut and Nipple
YR90-15FL-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple



106RM - LM SERIES


Single stage / Flowgauge regulators



FEATURES:

+ The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:



SPECIFICATIONS:

Model	Gas	Flow Range (LPM)	Inlet Connection
106RM-5LM-540	OXYGEN	1/2~5	CGA-540 Nut and Nipple
106RM-5LM-540HT	OXYGEN	1/2~5	CGA-540 Hand tight Nut and Nipple
106RM-5LM-870	OXYGEN	1/2~5	CGA-870 Yoke (Pin index)
106RM-5LM-992	OXYGEN	1/2~5	CGA-992 Nut and Nipple
106RM-5LM-992HT	OXYGEN	1/2~5	CGA-992 Hand tight Nut and Nipple
106RM-8LM-540	OXYGEN	1~8	CGA-540 Nut and Nipple
106RM-8LM-540HT	OXYGEN	1~8	CGA-540 Hand tight Nut and Nipple
106RM-8LM-870	OXYGEN	1~8	CGA-870 Yoke (Pin index)
106RM-8LM-992	OXYGEN	1~8	CGA-992 Nut and Nipple
106RM-8LM-992HT	OXYGEN	1~8	CGA-992 Hand tight Nut and Nipple
106RM-15LM-540	OXYGEN	1~15	CGA-540 Nut and Nipple
106RM-15LM-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
106RM-15LM-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
106RM-15LM-992	OXYGEN	1~15	CGA-992 Nut and Nipple
106RM-15LM-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple

HOW TO ORDER:

101RM-15 LM- XXX-SET

Flow rate
Including 5, 8, 15 (LPM)

Means flowgauge regulator.

Inlet connection
CGA540, CGA540HT, CGA870, CGA992, CGA992HT, CGA346, CGA950, CGA320, CGA320HT and CGA940.

SET
"SET" means the whole set including the combination of accessories

101RM - LM SERIES

Single stage / Flowgauge regulators

FLOWGAUGE REGULATORS

19

2" diameter outlet flow gauge.

Maximum inlet pressure 3000psi (ANSI)B-40.1/210Bar (CE).

Sintered filter inlet for additional safety and longer life.

Reliable external safety relief valve.

Neoprene diaphragm type.

101RM-15LM

+ Chrome-plated brass body with all brass high pressure chamber.
+ Chrome forged brass body is also available.

FEATURES:

+ The international standard inlet and outlet thread sizes are available.


SPECIFICATIONS:

Model	Gas	Flow Range (LPM)	Inlet Connection
101RM-5LM-540	OXYGEN	1/2~5	CGA-540 Nut and Nipple
101RM-5LM-540HT	OXYGEN	1/2~5	CGA-540 Hand tight Nut and Nipple
101RM-5LM-870	OXYGEN	1/2~5	CGA-870 Yoke (Pin index)
101RM-5LM-992	OXYGEN	1/2~5	CGA-992 Nut and Nipple
101RM-5LM-992HT	OXYGEN	1/2~5	CGA-992 Hand tight Nut and Nipple
101RM-8LM-540	OXYGEN	1~8	CGA-540 Nut and Nipple
101RM-8LM-540HT	OXYGEN	1~8	CGA-540 Hand tight Nut and Nipple
101RM-8LM-870	OXYGEN	1~8	CGA-870 Yoke (Pin index)
101RM-8LM-992	OXYGEN	1~8	CGA-992 Nut and Nipple
101RM-8LM-992HT	OXYGEN	1~8	CGA-992 Hand tight Nut and Nipple
101RM-15LM-540	OXYGEN	1~15	CGA-540 Nut and Nipple
101RM-15LM-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
101RM-15LM-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
101RM-15LM-992	OXYGEN	1~15	CGA-992 Nut and Nipple
101RM-15LM-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple

20

206RM-LM SERIES


Two stage / Flowgauge regulators



FEATURES:

- +The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:




SPECIFICATIONS:

Model	Gas	Flow Range (LPM)	Inlet Connection
206RM-5LM-540	OXYGEN	1/2~5	CGA-540 Nut and Nipple
206RM-5LM-540HT	OXYGEN	1/2~5	CGA-540 Hand tight Nut and Nipple
206RM-5LM-870	OXYGEN	1/2~5	CGA-870 Yoke (Pin index)
206RM-5LM-992	OXYGEN	1/2~5	CGA-992 Nut and Nipple
206RM-5LM-992HT	OXYGEN	1/2~5	CGA-992 Hand tight Nut and Nipple
206RM-8LM-540	OXYGEN	1~8	CGA-540 Nut and Nipple
206RM-8LM-540HT	OXYGEN	1~8	CGA-540 Hand tight Nut and Nipple
206RM-8LM-870	OXYGEN	1~8	CGA-870 Yoke (Pin index)
206RM-8LM-992	OXYGEN	1~8	CGA-992 Nut and Nipple
206RM-8LM-992HT	OXYGEN	1~8	CGA-992 Hand tight Nut and Nipple
206RM-15LM-540	OXYGEN	1~15	CGA-540 Nut and Nipple
206RM-15LM-540HT	OXYGEN	1~15	CGA-540 Hand tight Nut and Nipple
206RM-15LM-870	OXYGEN	1~15	CGA-870 Yoke (Pin index)
206RM-15LM-992	OXYGEN	1~15	CGA-992 Nut and Nipple
206RM-15LM-992HT	OXYGEN	1~15	CGA-992 Hand tight Nut and Nipple



106RM - PS SERIES


Single stage / Preset regulators



FEATURES:

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.


HOW TO ORDER:



SPECIFICATIONS:

Model	Gas	Inlet Connection
106RM-50PS-540(-V)	OXYGEN	CGA-540 Nut and Nipple
106RM-50PS-540HT(-V)	OXYGEN	CGA-540 Hand tight Nut and Nipple
106RM-50PS-870(-V)	OXYGEN	CGA-870 Yoke (Pin index)
106RM-50PS-346(-V)	MEDICAL AIR	CGA-346 Nut and Nipple
106RM-50PS-346HT(-V)	MEDICAL AIR	CGA-346 Hand tight Nut and Nipple
106RM-50PS-950(-V)	MEDICAL AIR	CGA-950 Yoke (Pin index)
106RM-50PS-326(-V)	N2O	CGA-326 Nut and Nipple
106RM-50PS-326HT(-V)	N2O	CGA-326 Hand tight Nut and Nipple
106RM-50PS-910(-V)	N2O	CGA-910 Yoke (Pin index)
106RM-50PS-320(-V)	CO2	CGA-320 Nut and Nipple
106RM-50PS-320HT(-V)	CO2	CGA-320 Hand tight Nut and Nipple
106RM-50PS-940(-V)	CO2	CGA-940 Yoke (Pin index)
106RM-50PS-580(-V)	Nitrogen	CGA-580 Nut and Nipple
106RM-50PS-580HT(-V)	Nitrogen	CGA-580 Hand tight Nut and Nipple
106RM-50PS-960(-V)	Nitrogen	CGA-960 Yoke (Pin index)


HOW TO ORDER:



101RM - PS SERIES
Single stage / Preset regulators

PRESET REGULATORS

25



FEATURES:

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.


SPECIFICATIONS:

Model	Gas	Inlet Connection
101RM-50PS-540(-V)	OXYGEN	CGA-540 Nut and Nipple
101RM-50PS-540HT(-V)	OXYGEN	CGA-540 Hand tight Nut and Nipple
101RM-50PS-870(-V)	OXYGEN	CGA-870 Yoke (Pin index)
101RM-50PS-346(-V)	MEDICAL AIR	CGA-346 Nut and Nipple
101RM-50PS-346HT(-V)	MEDICAL AIR	CGA-346 Hand tight Nut and Nipple
101RM-50PS-950(-V)	MEDICAL AIR	CGA-950 Yoke (Pin index)
101RM-50PS-326(-V)	N2O	CGA-326 Nut and Nipple
101RM-50PS-326HT(-V)	N2O	CGA-326 Hand tight Nut and Nipple
101RM-50PS-910(-V)	N2O	CGA-910 Yoke (Pin index)
101RM-50PS-320(-V)	CO2	CGA-320 Nut and Nipple
101RM-50PS-320HT(-V)	CO2	CGA-320 Hand tight Nut and Nipple
101RM-50PS-940(-V)	CO2	CGA-940 Yoke (Pin index)
101RM-50PS-580(-V)	Nitrogen	CGA-580 Nut and Nipple
101RM-50PS-580HT(-V)	Nitrogen	CGA-580 Hand tight Nut and Nipple
101RM-50PS-960(-V)	Nitrogen	CGA-960 Yoke (Pin index)

26

FR120-PS SERIES


Single stage / Preset regulators



FEATURES:

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.


HOW TO ORDER:



SPECIFICATIONS:


Model	Gas	Inlet Connection
FR120-50PS-540(-V)	OXYGEN	CGA-540 Nut and Nipple
FR120-50PS-540HT(-V)	OXYGEN	CGA-540 Hand tight Nut and Nipple
FR120-50PS-870(-V)	OXYGEN	CGA-870 Yoke (Pin index)
FR120-50PS-346(-V)	MEDICAL AIR	CGA-346 Nut and Nipple
FR120-50PS-346HT(-V)	MEDICAL AIR	CGA-346 Hand tight Nut and Nipple
FR120-50PS-950(-V)	MEDICAL AIR	CGA-950 Yoke (Pin index)
FR120-50PS-326(-V)	N2O	CGA-326 Nut and Nipple
FR120-50PS-326HT(-V)	N2O	CGA-326 Hand tight Nut and Nipple
FR120-50PS-910(-V)	N2O	CGA-910 Yoke (Pin index)
FR120-50PS-320(-V)	CO2	CGA-320 Nut and Nipple
FR120-50PS-320HT(-V)	CO2	CGA-320 Hand tight Nut and Nipple
FR120-50PS-940(-V)	CO2	CGA-940 Yoke (Pin index)
FR120-50PS-580(-V)	Nitrogen	CGA-580 Nut and Nipple
FR120-50PS-580HT(-V)	Nitrogen	CGA-580 Hand tight Nut and Nipple
FR120-50PS-960(-V)	Nitrogen	CGA-960 Yoke (Pin index)

HOW TO ORDER:



YR88 - PS SERIES
Single stage / Preset regulators

PRESET REGULATORS



FEATURES:

- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + The international standard inlet and outlet thread sizes are available.


SPECIFICATIONS:

Model	Gas	Inlet Connection
YR88-50PS-540(-V)	OXYGEN	CGA-540 Nut and Nipple
YR88-50PS-540HT(-V)	OXYGEN	CGA-540 Hand tight Nut and Nipple
YR88-50PS-870(-V)	OXYGEN	CGA-870 Yoke (Pin index)
YR88-50PS-346(-V)	MEDICAL AIR	CGA-346 Nut and Nipple
YR88-50PS-346HT(-V)	MEDICAL AIR	CGA-346 Hand tight Nut and Nipple
YR88-50PS-950(-V)	MEDICAL AIR	CGA-950 Yoke (Pin index)
YR88-50PS-326(-V)	N2O	CGA-326 Nut and Nipple
YR88-50PS-326HT(-V)	N2O	CGA-326 Hand tight Nut and Nipple
YR88-50PS-910(-V)	N2O	CGA-910 Yoke (Pin index)
YR88-50PS-320(-V)	CO2	CGA-320 Nut and Nipple
YR88-50PS-320HT(-V)	CO2	CGA-320 Hand tight Nut and Nipple
YR88-50PS-940(-V)	CO2	CGA-940 Yoke (Pin index)
YR88-50PS-580(-V)	Nitrogen	CGA-580 Nut and Nipple
YR88-50PS-580HT(-V)	Nitrogen	CGA-580 Hand tight Nut and Nipple
YR88-50PS-960(-V)	Nitrogen	CGA-960 Yoke (Pin index)



106RM-P SERIES


Single stage / Adjustable pressure regulators



FEATURES:

- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:



SPECIFICATIONS:


Model	Gas	Pressure Range (PSI)	Inlet Connection
106RM-50P-540	OXYGEN	0~50	CGA-540 Nut and Nipple
106RM-50P-540HT	OXYGEN	0~50	CGA-540 Hand tight Nut and Nipple
106RM-50P-870	OXYGEN	0~50	CGA-870 Yoke (Pin index)
106RM-50P-992	OXYGEN	0~50	CGA-992 Nut and Nipple
106RM-50P-992HT	OXYGEN	0~50	CGA-992 Hand tight Nut and Nipple
106RM-100P-540	OXYGEN	0~100	CGA-540 Nut and Nipple
106RM-100P-540HT	OXYGEN	0~100	CGA-540 Hand tight Nut and Nipple
106RM-100P-870	OXYGEN	0~100	CGA-870 Yoke (Pin index)
106RM-100P-992	OXYGEN	0~100	CGA-992 Nut and Nipple
106RM-100P-992HT	OXYGEN	0~100	CGA-992 Hand tight Nut and Nipple

206RM-P SERIES

Two stage / Adjustable pressure regulators

ADJUSTABLE PRESSURE REGULATORS


33



FEATURES:

- + The international standard inlet and outlet thread sizes are available.
- + Two-stage design for flow control accuracy purpose such as calibration and testing.

HOW TO ORDER:




SPECIFICATIONS:

Model	Gas	Pressure Range (PSI)	Inlet Connection
206RM-50P-540	OXYGEN	0~50	CGA-540 Nut and Nipple
206RM-50P-540HT	OXYGEN	0~50	CGA-540 Hand tight Nut and Nipple
206RM-50P-870	OXYGEN	0~50	CGA-870 Yoke (Pin index)
206RM-50P-992	OXYGEN	0~50	CGA-992 Nut and Nipple
206RM-50P-992HT	OXYGEN	0~50	CGA-992 Hand tight Nut and Nipple
206RM-100P-540	OXYGEN	0~100	CGA-540 Nut and Nipple
206RM-100P-540HT	OXYGEN	0~100	CGA-540 Hand tight Nut and Nipple
206RM-100P-870	OXYGEN	0~100	CGA-870 Yoke (Pin index)
206RM-100P-992	OXYGEN	0~100	CGA-992 Nut and Nipple
206RM-100P-992HT	OXYGEN	0~100	CGA-992 Hand tight Nut and Nipple

34

MF SERIES

EMS Regulators



FEATURES:


- + Maximum inlet pressure 2000psi (ANSI) / 140 bar (CE).
- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + Light-weight aluminum body.
- + 1-1/2" gauge with protector.
- + Specially drilled orifices provide superior liter flow accuracy with 12 flow selections.
- + Sintered filter for additional safety and longer service life.
- + The international standard inlet and outlet thread sizes are available.

SPECIFICATIONS:

Model	Flow Selection (LPM)	Inlet Connection	Outlet Connection
108MF-15LM-540(HT)	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	CGA-540 / CGA-540 Hand Tight	Barb outlet to connect with cannula / connector connect with humidifier.
108MF-15LM-870	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	CGA-870 (Pin Index)	Barb outlet to connect with cannula / connector connect with humidifier.
108MF-15LM-992(HT)	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15	CGA-992 / CGA-992 Hand Tight	Barb outlet to connect with cannula / connector connect with humidifier.

MFB SERIES

EMS Regulators



FEATURES:

- + Maximum inlet pressure 2000psi (ANSI) / 140 bar (CE).
- + The standard preset outlet pressure is 50psi (ANSI) or 4bar (CE), for stabilizing and accurate flowing.
- + Durable chrome-plated brass body.
- + 1-1/2" gauge with protector.
- + Specially drilled orifices provide superior liter flow accuracy with 12 flow selections.
- + Sintered filter for additional safety and longer service life.
- + The international standard inlet and outlet thread sizes are available.

SPECIFICATIONS:


Model	Flow Selection (LPM)	Inlet Connection	Outlet Connection
208MFB-25LM-540(HT)	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15, 25	CGA-540 / CGA-540 Hand Tight	Barb outlet to connect with cannula / connector connect with humidifier.
208MFB-25LM-870	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15, 25	CGA-870 (Pin Index)	Barb outlet to connect with cannula / connector connect with humidifier.
208MFB-25LM-992(HT)	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 12, 15, 25	CGA-992 / CGA-992 Hand Tight	Barb outlet to connect with cannula / connector connect with humidifier.

801 SERIES

Tube type / Medical flowmeters for regulators

TUBE FLOWMETERS

37



Polycarbonate inner and outer tubes provide greater accuracy, durability & safety.

Stainless Steel filter.

Chrome-plated brass body.

Outlet flow control valve 90° or 180° are available.


807 Flow selector
Designed to quickly switch to mask or humidifier.

With quick-connect adaptors.

FEATURES:

- + Uncompensated design accurately measures from 0~10 LPM or 0~15 LPM.
- + Calibrated at 50psi (ANSI) or 4bar (CE).
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:



38

822 SERIES

Tube type / Medical double flowmeters with quick-connect adaptor

Stainless Steel filter.

Polycarbonate inner and outer tubes provide greater accuracy, durability & safety.

Metal outlet flow control valve horizontal 180° to inlet connection.


Chrome-plated brass body.

PT-001 Power Take-off

FEATURES:

- + Double flowmeters type.
- + Uncompensated design accurately measures from 0~10 LPM or 0~15 LPM.
- + Calibrated at 50psi (ANSI) or 4bar (CE).
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:



825 SERIES

Tube type / Medical double flowmeters with quick-connect adaptor

826 SERIES

Tube type / Medical double flowmeters with quick-connect adaptor

TUBE FLOWMETERS

Stainless Steel filter.

Polycarbonate inner and outer tubes provide greater durability & safety.

Metal outlet flow control valve horizontal 180° to inlet connection.

Chrome-plated brass body construction for HIGH durability.

Aluminum body

HOW TO ORDER:

82xxp- OX- 15- XX-XX-SET

Outlook identification

825: Chrome-plated brass body
826: Aluminum body

None: means inlet and flow control valve are horizontal (180°).

Add "P" last: means including Power Take-off.

Flowmeter rate

O1: 1LPM
O3: 3.5LPM
O8: 8LPM
15: 15LPM

Gas Type

OX: Oxygen
CO: CO2
AI: Medical Air
N2: N2O
HE: Helium

Outlet connection

9/16: Male 9/16"
G 3/8: Male G 3/8
G 1/4: Male G 1/4
M12: Male M12
M16: Male M16
TN: Nipple Connection

SET

Means the whole set includes accessories.

Inlet connection

AM: Australian Male
BP: British Probe
CM: Chemetron Male
DH: DISS Handtight
DM: DISS Male
DN: Diss Nut and Nipple
F2: 1/8" FNPT
F4: 1/4" FNPT
FM: French Male

GM: Germany Male
JM: Japanese Male
M2: 1/8" MNPT
M4: 1/4" MNPT
MM: Medsta Male
XM: Oxequip Male
OM: Ohmeda Male
PM: Puritan Male
SM: Schrader Male

FEATURES:

- + Double flowmeters type.
- + Calibrated at 50psi (ANSI) or 4bar (CE).
- + The international standard inlet and outlet thread sizes are available.

MFC803 SERIES

Click type / Medical flowmeter with quick-connect adaptor



FEATURES:

- + Click type (Dial style)/single flowmeter.
- + Precision engineered accurately measures from 0~4LPM, 0~8LPM, 0~15LPM or 0~25LPM.
- + Calibrated at 50psi (ANSI) or 4bar (CE).
- + Optional power take-off units provide 4 bar (60psi) or 3.5 bar (50psi) additional outlet pressure.
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:

MFC803xP-OX- 15- XX-XX-SET

Outlook identification

None: means inlet and flow control valve are horizontal (180°).

V: means inlet and flow control valve are vertical (90°).

Add "P" last: means including Power Take-off.

Gas Type

OX: Oxygen
CO: CO2
AI: Medical Air
N2: N2O
HE: Helium

Flowmeter rate

O4: 4LPM
O8: 8LPM
15: 15LPM
25: 25LPM

Outlet connection

9/16: Male 9/16"
G 3/8: Male G 3/8
G 1/4: Male G 1/4
M12: Male M12
M16: Male M16
TN: Nipple Connection

SET

Means the whole set includes accessories.

Inlet connection

AM: Australian Male	GM: Germany Male
BP: British Probe	JM: Japanese Male
CM: Chemetron Male	M2: 1/8" MNPT
DH: DISS Handtight	M4: 1/4" MNPT
DM: DISS Male	MM: Medsta Male
DN: Diss Nut and Nipple	XM: Oxequip Male
F2: 1/8" FNPT	OM: Ohmeda Male
F4: 1/4" FNPT	PM: Puritan Male
FM: French Male	SM: Schrader Male

MFC823 SERIES

Click type / Medical double flowmeters with quick-connect adaptor



FEATURES:

- + Click type (Dial style)/double flowmeters.
- + Precision engineered accurately measures from 0~4LPM, 0~8LPM, 0~15LPM or 0~25LPM.
- + Calibrated at 50psi (ANSI) or 4bar (CE).
- + Optional power take-off units provide 4 bar (60psi) or 3.5 bar (50psi) additional outlet pressure.
- + The international standard inlet and outlet thread sizes are available.

HOW TO ORDER:

MFC823xP-OX- 15- XX-XX-SET

Outlook identification

None: means inlet and flow control valve are horizontal (180°).

V: means inlet and flow control valve are vertical (90°).

Add "P" last: means including Power Take-off.

Gas Type

OX: Oxygen
CO: CO2
AI: Medical Air
N2: N2O
HE: Helium

Flowmeter rate

O4: 4LPM
O8: 8LPM
15: 15LPM
25: 25LPM

Outlet connection

9/16: Male 9/16"
G 3/8: Male G 3/8
G 1/4: Male G 1/4
M12: Male M12
M16: Male M16
TN: Nipple Connection

SET

Means the whole set includes accessories.

Inlet connection

AM: Australian Male	GM: Germany Male
BP: British Probe	JM: Japanese Male
CM: Chemetron Male	M2: 1/8" MNPT
DH: DISS Handtight	M4: 1/4" MNPT
DM: DISS Male	MM: Medsta Male
DN: Diss Nut and Nipple	XM: Oxequip Male
F2: 1/8" FNPT	OM: Ohmeda Male
F4: 1/4" FNPT	PM: Puritan Male
FM: French Male	SM: Schrader Male



IOF SERIES

Quick connects & adaptors

HOW TO ORDER:

IOF- XX - XX - XX

Quick Connect Types

AU: Australian	GS: Germany
BS: British Probe	JA: Japanese
CH: Chemetron	MS: Medsta
DH: DISS Handtight	XQ: Oxequip
DM: DISS Male	OH: Ohmeda
DN: DISS Nut & Nipple	SH: Schrader
PB: Puritan Bennett	FR: French

Inlet connection

M2: 1/8" MNPT	TN1/4: 1/4" Hose Barb
F2: 1/8" FNPT	TN5/16: 5/16" Hose Barb
M4: 1/4" MNPT	DN: DISS Hex Nut
DH: DISS Handtight Nut & Nipple	

Gas Type

OX: Oxygen	N2: Nitrous Oxide
AI: Medical Air	CO: CO2
VA: Vacuum	WA: WAGD

Medical Color coded :

U.S.A Standard	ISO color coding
■ : Oxygen	□ : Oxygen
■ : Nitrous Oxide	■ : Nitrous Oxide
■ : Medical Air	■ : Medical Air
■ : Nitrogen	■ : Nitrogen
■ : Carbon Dioxide	■ : Carbon Dioxide
■ : Helium	■ : Helium

Ohmeda Type

Gases	Connections	Code
 Oxygen	 1/8" NPT Male	IOF-OH-OX-M2
	 1/4" or 5/16" Hose Barb	IOF-OH-OX-TN1/4 IOF-OH-OX-TN5/16
 Medical Air	 1/8" NPT Male	IOF-OH-AI-M2
	 1/4" or 5/16" Hose Barb	IOF-OH-AI-TN1/4 IOF-OH-AI-TN5/16
 Vacuum	 1/8" NPT Male	IOF-OH-VA-M2
	 1/4" or 5/16" Hose Barb	IOF-OH-VA-TN1/4 IOF-OH-VA-TN5/16
 Nitrous Oxide	 1/8" NPT Male	IOF-OH-N2-M2
	 1/4" or 5/16" Hose Barb	IOF-OH-N2-TN1/4 IOF-OH-N2-TN5/16

Chemetron Type

Gases	Connections	Code
 Oxygen	 1/8" NPT Male	IOF-CH-OX-M2
	 1/4" or 5/16" Hose Barb	IOF-CH-OX-TN1/4 IOF-CH-OX-TN5/16
 Medical Air	 1/8" NPT Male	IOF-CH-AI-M2
	 1/4" or 5/16" Hose Barb	IOF-CH-AI-TN1/4 IOF-CH-AI-TN5/16
 Vacuum	 1/8" NPT Male	IOF-CH-VA-M2
	 1/4" or 5/16" Hose Barb	IOF-CH-VA-TN1/4 IOF-CH-VA-TN5/16
 Nitrous Oxide	 1/8" NPT Male	IOF-CH-N2-M2
	 1/4" or 5/16" Hose Barb	IOF-CH-N2-TN1/4 IOF-CH-N2-TN5/16

French Type: AFNOR NFS 90-116

Gases	Connections	Code
 Oxygen	1/8" NPT Male	IOF-FR-OX-M2 IOF-FR-AI-M2 IOF-FR-VA-M2 IOF-FR-N2-M2
 Medical Air		
 Vacuum	1/4" or 5/16" Hose Barb	IOF-FR-OX-TN1/4 IOF-FR-AI-TN1/4 IOF-FR-VA-TN5/16 IOF-FR-N2-TN1/4
 Nitrous Oxide		

IOF SERIES

Quick connects & adaptors

Specification

DISS Male Type

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide WAGD Carbon Dioxide Nitrogen	1/8" NPT Male	IOF-DM-OX-M2 IOF-DM-AI-M2 IOF-DM-VA-M2 IOF-DM-N2-M2 IOF-DM-WA-M2 IOF-DM-CO-M2 IOF-DM-NI-M2

DISS Nut & Nipple Type

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide WAGD Carbon Dioxide Nitrogen	1/4" or 5/16" Hose Barb	IOF-DN-OX-TN1/4 IOF-DN-AI-TN1/4 IOF-DN-VA-TN5/16 IOF-DN-N2-TN1/4 IOF-DN-WA-TN5/16 IOF-DN-CO-TN1/4 IOF-DN-NI-TN1/4

DISS Handtight Type

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide WAGD Carbon Dioxide Nitrogen	1/8" NPT Male	IOF-DH-OX-M2 IOF-DH-AI-M2 IOF-DH-VA-M2 IOF-DH-N2-M2 IOF-DH-WA-M2 IOF-DH-CO-M2 IOF-DH-NI-M2
	Oxygen Medical Air Vacuum Nitrous Oxide WAGD Carbon Dioxide Nitrogen	1/4" or 5/16" Hose Barb	IOF-DH-OX-TN1/4 IOF-DH-AI-TN1/4 IOF-DH-VA-TN5/16 IOF-DH-N2-TN1/4 IOF-DH-WA-TN5/16 IOF-DH-CO-TN1/4 IOF-DH-NI-TN1/4

Puritan Bennett Type

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide WAGD	1/8" NPT Male	IOF-PB-OX-M2 IOF-PB-AI-M2 IOF-PB-VA-M2 IOF-PB-N2-M2 IOF-PB-WA-M2
	Oxygen Medical Air Vacuum Nitrous Oxide WAGD	1/4" or 5/16" Hose Barb	IOF-PB-OX-TN1/4 IOF-PB-AI-TN1/4 IOF-PB-VA-TN5/16 IOF-PB-N2-TN1/4 IOF-PB-WA-TN5/16

Japanese Type: JIS T7101

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide	1/8" NPT Male	IOF-JA-OX-M2 IOF-JA-AI-M2 IOF-JA-VA-M2 IOF-JA-N2-M2
	Oxygen Medical Air Vacuum Nitrous Oxide	1/4" or 5/16" Hose Barb	IOF-JA-OX-TN1/4 IOF-JA-AI-TN1/4 IOF-JA-VA-TN5/16 IOF-JA-N2-TN1/4

Germany Type: DIN 13260-2

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide Carbon Dioxide Nitrogen	1/8" NPT Male	IOF-GS-OX-M2 IOF-GS-AI-M2 IOF-GS-VA-M2 IOF-GS-N2-M2 IOF-GS-CO-M2 IOF-GS-NI-M2
	Oxygen Medical Air Vacuum Nitrous Oxide Carbon Dioxide Nitrogen	1/4" or 5/16" Hose Barb	IOF-GS-OX-TN1/4 IOF-GS-AI-TN1/4 IOF-GS-VA-TN5/16 IOF-GS-N2-TN1/4 IOF-GS-CO-TN1/4 IOF-GS-NI-TN1/4

British Type: BS 5682-1998

	Gases	Connections	Code
	Oxygen Medical Air Vacuum Nitrous Oxide	1/8" NPT Male	IOF-BS-OX-M2 IOF-BS-AI-M2 IOF-BS-VA-M2 IOF-BS-N2-M2
	Oxygen Medical Air Vacuum Nitrous Oxide	1/4" or 5/16" Hose Barb	IOF-BS-OX-TN1/4 IOF-BS-AI-TN1/4 IOF-BS-VA-TN5/16 IOF-BS-N2-TN1/4

103 SERIES Pressure Gauges



Specification

Model	psi/kpa	LPM	m/m	Thread
10311	4000psi		53	NPT 1/4" UL
10312	4000psi		53	NPT 1/4" UL
10313		15	53	NPT 1/4" RU
10314		15	53	NPT 1/4" RU

SR-AG SERIES Analog Gauges

Model	Specification
SR-AG-760(A)-ML	Analog Gauge with Lens 760 mmHg
SR-AG-300(A)-ML	Analog Gauge with Lens 300 mmHg
SR-AG-160(A)-ML	Analog Gauge with Lens 160 mmHg

*(A) means Anti-clockwise



SR-AG-760(A)-ML



SR-AG-300(A)-ML



SR-AG-160(A)-ML

Humidifiers/Disposable

HUMIDIFIERS

55



HM-005

- + Polyethylene container of 200ml capacity
- + Head in polycarbonate for durability
- + 9/16"x18N thread. Other thread sizes are available.
- + With 2PSI safety valve.
- + Green color
- + Reuse

HM-006

- + Polyethylene container of 150ml capacity
- + Head in polycarbonate for durability
- + 9/16"x18N thread. Other thread sizes are available.
- + Reuse



HM-007

- + Humidifier equipped with 2PSI safety valve, and standard DISS(9/16") thread connector for flow meter
- + Container's volume: 200ml



HM-007A

- + Humidifier equipped with 14PSI safety valve, and standard DISS(9/16") thread connector for flow meter
- + Container's volume: 250ml



HM-008

- + Humidifier equipped with 4PSI safety valve, and standard DISS(9/16") thread connector for flow meter
- + Container's volume: 500ml



Accessories




HM-029 OXYGEN Mask

- + Simple Oxygen mask for adult
- + Medium concentration with 7" tubing elastic strap, adjustable nose clip.
- + PVC Material
- + CE2460

HM-017 Disposable Cannula

- + Simple oxygen cannula tube for nose use.
- + Medium concentration with 7" tubing.
- + Material: polypropylene, polyethylene, polyvinylchloride.
- + CE2460



YK-801 YOKE

- + Made of forged brass, fully chrome-plated, 1/4"NPT inlet thread

Model	Gas
YK-801	OXYGEN
YK-802	AIR
YK-803	NITROUS OXIDE
YK-804	CARBON DIOXIDE

HM-CG 、HM-CW Flow Meter Connector

- + Diss (9/16") thread for follow meter and swivel 4-6 barb style connector for tubing.
- + Material: Nylon
- + Disposable



56



Note

