

WEAVER
and company

565 Nucla Way, Unit B
Aurora, Colorado 80011 USA

Ten20[®]

CONDUCTIVE



8oz (228g) Jars

Sold as a 3-pack
ITEM #10-20-8



4oz (114g) Jars

Sold as a 3-pack
ITEM #10-20-4



4oz (114g) Tubes

Sold as a 3-pack
ITEM #10-20-4T



Single-Use Cups

Sold as a 24-pack
ITEM #10-20-1

About Ten20 Conductive

Ten20 Conductive is an easy-to-use, adhesive, conductive paste designed to hold EEG electrodes in place, thereby eliminating disruptions of reattaching the electrodes during neurodiagnostic examinations.

Expert technologists and technicians around the world rely on the optimal balance of adhesiveness and conductivity of Ten20 Conductive Paste to achieve quality test results.

Ten20 Conductive Paste's washable, non-drying formula makes for a very easy clean up.

Target Markets

EEG Examinations ●

Sleep Diagnostics ●

Evoked Potential Procedures ●

Audiology ●

Biofeedback ●

Contact Us

(800) 525-2130

www.weaverandcompany.com
info@weaverandcompany.com

WEAVER
and company

565 Nucla Way, Unit B
Aurora, Colorado 80011 USA

Nuprep[®]

Skin Prep Gel



4oz (114g) Tubes

Sold as a 3-pack
ITEM #10-30



25g Tubes

Sold as a 6-pack
ITEM #10-61



25g Tubes

Sold as a 20-pack
ITEM #10-65



Single Use Cups

Sold as a 24-pack
ITEM #10-91

Target Markets

- EEG Examinations
- Sleep Diagnostics
- Evoked Potential Procedures
- Intraoperative Monitoring
- Audiology
- Biofeedback
- Cardiac Stress Tests
- Holter Monitoring
- Cardiac Rehabilitation Monitoring

About Nuprep Skin Prep Gel

For quality tracings, proper skin preparation at every electrode site is crucial.

Nuprep Skin Prep Gel gently abrades and moistens the top layer of skin. This lowers the skin impedance with minimal skin irritation and patient discomfort.

Nuprep virtually eliminates the undesirable effects of diaphoresis and unwanted artifact, thus providing the best opportunity to achieve accurate neurodiagnostic and cardiadiagnostic tracings.

Contact Us

(800) 525-2130

www.weaverandcompany.com
info@weaverandcompany.com